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Optilene® Mesh

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Instructions for use

Optilene® Mesh

Description
Optilene® Mesh is a mesh implant for re-inforcement of connective tissue structures. It is constructed from monofilament polypropylene, that has been knitted to a thin and elastic shape-stable mesh. After implantation the Optilene® Mesh adapts to the longitudinal and latitudinal expansions taking place in the connective tissue.

Optilene® Mesh is coloured blue with copper phthalocyanine (Phthalocyaninato (2-) copper) for a better visibility.

Optilene® Mesh does not possess any independent pharmacological properties. The polypropylene mesh is biostable and is not degraded in the body.

Indications
Re-inforcement material in connective tissue weakness:
- for prosthetic hernioplasty
- for reconstruction of the chest wall
- for re-inforcement of fascial tissue if a non-absorbable re-inforcement material is required

Optilene® Mesh is suitable for both conventional and laparoscopic application.

Application
Select Optilene® Mesh from the product range according to the size of the defect and cut it to fit the defect.

When repairing inguinal hernias, the mesh should, in order to prevent recurrence, be large enough to project beyond the pubic tubercle and pass reliably round the spermatic cord at the inner ring.

Non-absorbable suture material or clips are recommended for fixing Optilene® Mesh. They should be placed \geq 10 mm from the edge of the mesh. Atraumatic round-bodied needles make suturing still possible without greater damage to the implant.

Mode of action
The healing-in process is primarily influenced by the porosity of the implant. When Optilene® Mesh is implanted, connective

tissue cells soon begin to migrate into the pores of the implant. Here there can be a brief, mild tissue reaction, which, however, supports a rapid, problem-free healing-in process. The defined pore size brings about even, tensionfree penetration by the tissue. The Optilene® Mesh is characterized by its thin but extremely shape-stable mesh wall. It adapts optimally to the longitudinal and latitudinal motion of the patient as a result of body movement and muscle contraction.

Contra-indications
Optilene® Mesh should not be implanted:
- in contaminated and infected areas
- in children during the growth phase

Warning notes
- Direct contact between Optilene® Mesh and the viscera (intestines etc) must be avoided to prevent adhesion.
- Only use non-absorbable material (suture material, staples, tack devices) together with Optilene® Mesh.
- Do not reuse the product: Infection hazard for patients and/or users and impairment of product functionality due to reuse. Risk of injury, illness or death due to contamination and/or impaired functionality of the product.
- No special storage conditions are required.

Notes / precautionary measures
- Before implanting Optilene® Mesh the surgeon must be familiar with the surgical technique, the specific application and the properties of Optilene® Mesh in vivo.
- Optilene® Mesh may only be used if the expiring date is undamaged and, hence, the sterility is guaranteed.
- Optilene® Mesh may not be used after the expiry date given on the packaging.
- Restertilisation of Optilene® Mesh is not recommended.

Side effects
Like all non-absorbable meshes, implantation of Optilene® Mesh can be associated with a limited period of local irritation in the area of the wound; a transient foreign body reaction can occasionally take the form of an inflammatory reaction. Infections and adhesions on contact with the viscera can occur as a result of the intervention.

Sterilization
Optilene® Mesh is sterilized using ethylene oxide.
Date of information: 03/2011

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