Gebruiksaanweisung

Histoacryl® LapFix – Cannula

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Histoacryl® LapFix – Cannula

Description

Histoacryl® LapFix – Cannula is a single-use medical device. The cannula is designed for exclusive use with Histoacryl® and Injekt® 2mL Luer Lock Solo to fix meshes during laparoscopic inguinal hernia repair, especially with the TAPP or TEP technique. The cannula is supplied sterile packaged with easy peel opening.

**Indications**

Histoacryl® LapFix –Cannula is indicated to deliver Histoacryl® for the laparoscopic fixation of surgical meshes, especially in inguinal hernia surgery (TAPP and TEP).

Operation

See schematic view.

– Preparation – Open the Histoacryl® ampoule by twisting off the tip (fig. 1). Prepare the syringe by pulling back the green plunger until you create 2mL space in the barrel (fig. 2).

– Transfer of the Histoacryl® solution to the syringe – Insert the open ampoule tip into the syringe body and firmly press the Histoacryl® ampoule the required times until the content of the ampoule is transferred (fig. 3). Repeat the transfer with the second Histoacryl® ampoule.

– Priming – Remove the air from the syringe until the Histoacryl® reaches the tip of the syringe (fig. 4).

– Assembly and ready to use – Attach the cannula to the purged syringe by screwing the Luer lock. Prime the device by pushing the green syringe plunger until the first drop of Histoacryl® appears in the cannula tip (fig 5. and 6 ).

– Application – Protect the tip of the cannula with the tube and introduce it through the 5mm trocar (or correspondent) (fig. 7 ). To deliver the adhesive, hold the tip of the applicator on the mesh at the desired fixation point and push the plunger carefully to deliver the amount of glue needed. After application of the adhesive, remove the cannula from the mesh surface to avoid the stickiness of the tip.

– Extraction of the device – Extract the device ensuring that the protective tube covers the tip of the cannula in order to prevent trocar contamination with the remaining Histoacryl® from the tip.

If the cannula tip becomes blocked during usage, extract the cannula from inside the patient and release the clogging by rubbing the tip with a gauze or cutting a portion of it with scissors.

Recommendations on product use

Contact of the cannula tip with accumulation of body fluid may lead to the blockage of the tip. It is recommended to avoid contact with said fluids, surrounding tissue or the mesh while the device it is not delivering the adhesive.

Contraindications

The only intended use of Histoacryl® LapFix – Cannula is to apply Histoacryl® in the laparoscopic inguinal hernia procedures. The device must not be used for other intended uses.

Histoacryl® LapFix – Cannula must be used solely with Histoacryl®, it must not be used with other adhesives.

**Warnings and precautions**

The cannula is designed to be used with a trocar of 5 mm or bigger. The cannula will not fit inside a trocar of less than 5 mm diameter.

If the cannula tip becomes clogged, the tip can be cut off or the blockage removed by rubbing with a sterile gauze. The operation must be performed outside the patient to avoid the release of particles inside the body.

The device is not designed for being used as a substitute of other surgical instruments during surgical procedures (e.g., needles, forceps, tweezers).

The same warnings and measures of caution indicated in the instructions for use of Histoacryl® must be considered when Histoacryl® is used in combination with the Histoacryl® LapFix –Cannula for Laparoscopic Hernia Repair.

Do not use the product if the packaging is opened or damaged. The applicator must not be re-sterilized. The user should be familiar with the surgical procedures and techniques involved in laparoscopic hernia repair using meshes as implants.

The device is designed for a single procedure use. Do not reuse it. Infection hazard for patients and/or users and impairment of product functionality due to reuse. Risk of illness or death due to contamination of the product.

Do not reuse the same cannula with a new reload of Histoacryl during the same procedure as there is a risk of the tip becoming clogged and stopping the flow of Histoacryl during the application.

The product must not be used after the expiry date. Prior to use, check the expiry date displayed on the product label. The cannula must be solely used with Histoacryl® and Injekt® 2mL Luer Lock Solo. The efficacy and safety usage of the cannula has not been proven with other components.

Adverse effects

Same adverse effects indicated in the instructions for use of Histoacryl® must be considered when Histoacryl® is used in combination with the Histoacryl® LapFix – Cannula for Laparoscopic Hernia Repair.

Sterilization

Histoacryl® LapFix –Cannula for Laparoscopic Hernia Repair has been sterilized with ethylene oxide gas.

Storage

No special storage conditions are required.

**Additional notes**

The instrument is supplied sterile for use in a single procedure only.

The applicator must not be re-sterilized.

Instructions for use of Histoacryl® must be considered when using the Histoacryl® LapFix –Device for Laparoscopic Hernia Repair.

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Title: FOLLETO HISTOACRYL LAPFIX CANNULA Initiator: Susana ? Calvillo

This document is signed electronically in compliance with the B. Braun electronic signature policies and procedures by following persons:

UserName: Calvillo, Susana (calvsues)  
Title: SUT-PR-ES04-Ayudante Técnico Sistema Impresión  
Date: Tuesday, 29 November 2016, 10:11 W. Europe Daylight Time  
Meaning: Document signed as Author  
=====

UserName: García, Guillermo (garcgues)  
Title: SU-PR-ES04-Technical support manager  
Date: Tuesday, 29 November 2016, 15:08 W. Europe Daylight Time  
Meaning: Precheck of Document  
=====

UserName: Piera, Eva (piermaes)  
Title: Senior Regulatory Affairs  
Date: Wednesday, 30 November 2016, 15:51 W. Europe Daylight Time  
Meaning: Approve Document  
=====

UserName: Turón, Pau (turopaes)  
Title: R&D, Regulatory Affairs and Quality Management Director CoE CT  
Date: Thursday, 01 December 2016, 13:40 W. Europe Daylight Time  
Meaning: Final Release of the Document  
=====