# Aesculap® Acculan 4

## en Instructions for use/Technical description

Oscillating saw GA331

**USA** Note for U.S. users

This Instructions for Use is NOT intended for United States users. Please discard. The Instructions for Use for United States users can be obtained by visiting our website at [www.aesculapusa.com](http://www.aesculapusa.com). If you wish to obtain a paper copy of the Instructions for Use, you may request one by contacting your local Aesculap representative or Aesculap's customer service at 1-800-282-9000. A paper copy will be provided to you upon request at no additional cost.

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

## Oscillating saw GA331

### Legend

- 1 Product (oscillating saw)
- 2 Trigger (for oscillation frequency control)
- 3 Safety catch
- 4 Battery bay
- 5 Locking pin
- 6 Sterile funnel
- 7 Battery
- 8 Lid
- 9 Cover release element
- 10 Battery removal aid
- 11 Rapid Action saw blade
- 12 Saw head
- 13 Button for unlocking saw blade
- 14 Saw blade coupling

The depictions are schematic only.

### Symbols on product and packages

	Caution Observe important safety information such as warnings and precautions in the instructions for use.
YYYY-MM	Maintenance label Indication of the next maintenance appointment (Date: Year-Month)
	Two-dimensional machine readable code The code contains a unique serial number that can be used for electronically tracking single instruments. The serial number is based on the worldwide standard sGTIN (GS1).
	Manufacturer
	Date of manufacture
	Follow the instructions for use
	Labeling of electrical and electronic devices pursuant to directive 2002/96/EG (WEEE)
	Classification Type BF
	Manufacturer's batch designation
	Manufacturer serial number
	Manufacturer order number

	Temperature limits during transport and storage
	Air humidity limits during transport and storage
	Atmospheric pressure limits during transport and storage
	Operating mode
	According to US federal law, this product may only be sold by a physician or by the order of a physician
	CE mark according to Regulation (EU) 2017/745
	Medical device
	Oscillation frequency control

### Contents

1.	About this document .....	3
1.1	Scope .....	3
1.2	Warnings .....	3
2.	General information .....	3
2.1	Intended use .....	3
2.2	Main functions and design characteristics .....	3
2.2.1	Operating mode .....	3
2.3	Indications .....	3
2.4	Absolute contraindications .....	3
2.5	Relative contraindications .....	3
3.	Safe handling .....	4
4.	Product description .....	4
4.1	Scope of supply .....	4
4.2	Components required for operation .....	4
4.3	Operating principle .....	4
5.	Preparation .....	4
6.	Working with the device .....	4
6.1	System set-up .....	5
6.1.1	Connecting the accessories .....	5
6.1.2	Inserting the rechargeable battery .....	5
6.1.3	Intraoperative battery change .....	5
6.1.4	Removing the rechargeable battery .....	5
6.1.5	Protection against inadvertent activation .....	5
6.1.6	Attaching and detaching the saw blade .....	6
6.2	Function checks .....	6
6.3	Safe operation .....	6
6.3.1	Operating the product .....	6

7.	Validated reprocessing procedure .....	7
7.1	General safety notes .....	7
7.2	General information .....	7
7.3	Reusable products .....	7
7.4	Preparations at the place of use.....	7
7.5	Preparation before cleaning .....	7
7.6	Product-specific safety instructions for the processing procedure .....	7
7.7	Manual cleaning with wipe disinfection .....	8
7.8	Automatic cleaning/disinfection with manual pre-cleaning ..	9
7.8.1	Manual pre-cleaning with a brush.....	9
7.8.2	Mechanical alkaline cleaning and thermal disinfection.....	10
7.9	Inspection, maintenance and checks.....	10
7.10	Packaging .....	10
7.11	Steam sterilization.....	10
7.12	Storage.....	10
8.	Maintenance .....	11
9.	Troubleshooting list.....	11
10.	Technical Service .....	12
11.	Accessories/Spare parts.....	12
12.	Technical data .....	13
12.1	Classification according to Regulation (EU) 2017/745.....	13
12.2	Performance data, information about standards .....	13
12.3	Operating mode .....	13
12.4	Ambient conditions.....	13
13.	Disposal .....	13

## 1. About this document

### 1.1 Scope

These instructions for use apply to the products listed below:

Art. no.	Name
GA331	Oscillating saw

► For instructions for use of specific items and information on material compatibility, see B. Braun eIFU at [eifu.bbraun.com](http://eifu.bbraun.com)

### 1.2 Warnings

Warnings indicate risks to patients, users, and/or the product that could develop while using the product. Warnings are marked as follows:

#### DANGER

Indicates a potential danger. If it is not prevented, death or severe injuries may result.

#### WARNING

Indicates a potential danger. Minor or moderate injuries may result if this is not prevented.

#### CAUTION

Indicates a potential damage to property. The product may become damaged if it is not prevented.

## 2. General information

### 2.1 Intended use

#### Task/Function

When combined with the appropriate tool, the Oscillating saw GA331, is used to work on hard tissue, cartilage and related materials, as well as bone replacement materials.

#### Application Environment

The product fulfills the requirements for type BF pursuant to IEC/DIN EN 60601-1 and is used in operating rooms in sterile environments of explosion risk areas (such as areas with pure oxygen or anesthesia gases).

### 2.2 Main functions and design characteristics

Oscillation frequency

min. 0 min<sup>-1</sup> to max. 13 000 min<sup>-1</sup>

#### 2.2.1 Operating mode

Operation with non-periodic load and speed changes (type S9 pursuant to IEC EN 60034-1)

- 30 second application, 30 second pause
- 7 repetitions
- 30 min cooling time
- Max. Temperature 48 °C

Electrical systems generally heat up during continual operation. It is advised to give the system a break after use to cool down, as listed in the table on operating mode.

Heating depends on the tool used and the load. After a certain number of repetitions, the system should cool down. This procedure prevents the system overheating as well as possible injury to the patient or user.

The user is responsible for the use and adherence to the pause sequence described.

### 2.3 Indications

Type and area of application depend on the selected tool.

### 2.4 Absolute contraindications

The product is not approved for use in the central nervous system or central circulatory system.

### 2.5 Relative contraindications

The safe and effective use of the product depends greatly on factors that only operators themselves can control. The information provided therefore only represents general conditions.

The clinically successful use of the product depends on the knowledge and experience of the surgeon. They must decide which structures can be sensibly treated and must take the safety information and warnings in the instruction for use into account.

### 3. Safe handling

#### ⚠ WARNING

Risk of injury and material damage if this product is not used as intended!

- Use the product only for its intended purpose.

#### ⚠ WARNING

Risk of injury and damage to property due to improper handling of the product!

- Follow the instructions for use of all products used.

- General risk factors associated with surgical procedures are not described in this documentation.
- It is the operating surgeon's responsibility to ensure that the surgical procedure is performed correctly.
- The operating surgeon must have a thorough understanding of both the hands-on and conceptual aspects of the established operating techniques.
- Remove the transport packaging and clean the new product, either manually or mechanically, prior to its initial sterilization.
- Check the functionality and proper condition before using the product.
- Observe "Guidance on Electromagnetic Compatibility (EMC) for Acculan Components" TA022450, see B. Braun eIFU at [eifu.bbraun.com](http://eifu.bbraun.com)
- To prevent damage caused by improper setup or operation, and in order not to compromise warranty and manufacturer liability:
  - Use the product only according to these instructions for use.
  - Follow the safety and maintenance instructions.
  - Only combine Aesculap products with each other.
- Ensure that the product and its accessories are operated and used only by persons with the requisite training, knowledge, or experience.
- Keep the instructions for use accessible for the user.
- Comply with valid standards.
- Ensure that the electrical installation of the room is consistent with the requirements of IEC/DIN EN.
- Do not operate the product in explosion-hazard areas.
- Sterilize product before use.
- When using ECCOS holder systems, observe relevant instructions for use TA009721, see B. Braun eIFU at [eifu.bbraun.com](http://eifu.bbraun.com)

#### Note

*The user must report all serious incidents occurring in connection with the product to the manufacturer and to the competent authority of the country in which the user is registered.*

### 4. Product description

#### 4.1 Scope of supply

Art. no.	Name
GA331	Oscillating saw
GA678	Sterile funnel
TA014543	Instructions for use for oscillating saw GA331 (flyer)

#### 4.2 Components required for operation

- NiMH battery long GA676 (charged)
- Sterile funnel GA678
- Cap GA675
- Saw blade (according to indication)

#### 4.3 Operating principle

The product 1 contains an electric motor, which is powered by a replaceable battery 7.

The charged non-sterile battery 7 is inserted into the product 1 via the sterile funnel 6 and sealed with the lid 9 so that it is sterile.

The oscillation frequency is controlled electronically and can be continuously regulated with the trigger for oscillation frequency control 2.

The product 1 has a coupling, which allows for the attachment of various saw blades 11. These saw blades 11 lock automatically when attached to the product.

The saw blade 11 can be removed by activating a push button 13.

### 5. Preparation

Aesculap assumes no liability if the following rules are not followed:

- Do not use products from open or damaged sterile packaging.
- Prior to use, inspect the product and its accessories for any visible damage.
- Use the products and their accessories only if they are in perfect technical condition.

### 6. Working with the device

#### ⚠ WARNING

Risk of infection and contamination!

Product is delivered unsterilized!

- Sterilize the product before use pursuant to the operating instructions.

#### ⚠ WARNING

Risk of injury and material damage due to accidental activation of the product!

- Products which are not being actively used must be secured against accidental activation (position OFF).

#### ⚠ WARNING

Risk of injury and material damage due to inappropriate use of tools!

- Always follow the safety advice and information given in the instructions for use.
- When coupling/uncoupling, handle tools with cutting edges with care.

**⚠ WARNING**

**Damage to the product if dropped!**

- Use the products only if they are in perfect technical condition, see Function check

**⚠ WARNING**

**Risk of burns to skin and tissue through blunt tools or if product has not been maintained properly!**

- Use tools only if they are in perfect condition.
- Replace blunt tools.
- Maintain the product properly, see maintenance guide.

## 6.1 System set-up

### 6.1.1 Connecting the accessories

Combinations of accessories that are not mentioned in the present instructions for use may only be employed if they are specifically intended for the respective application, and if they do not compromise the performance and safety characteristics of the products.

All configurations must fulfill the fundamental standard IEC/DIN EN 60601-1. The person connecting the devices with each other is responsible for the configuration and must ensure that the fundamental standard IEC/DIN EN 60601-1 or relevant national standards.

- Follow the instructions for use of individual accessories.
- In case of questions, contact your B. Braun/Aesculap partner or Aesculap Technical Service, address see Technical Service.

### 6.1.2 Inserting the rechargeable battery

**⚠ CAUTION**

**The product does not function when using the wrong battery!**

- Do not use battery GA666 (with red base) or GA346.
- Product 1 with the battery shaft 4 must be turned upwards and the sterile funnel 6 (sterile) attached, see Fig. A.
- Battery 7 (non-sterile) is inserted into the battery shaft 4 by a second (non-sterile) person.

**Note**

*After the battery has been inserted, multiple signals will sound, which indicate that the product is ready to use.*

- After insertion of the battery, have the sterile funnel 6 (non-sterile) removed by a second person.
- Attach the lid 8 (sterile) in such a manner that it clicks into place with both release catches 9.

**Note**

*The sterility of the product is only guaranteed with a correctly placed lid.*

### 6.1.3 Intraoperative battery change

**⚠ WARNING**

**Risk of injury from hot battery!**

**Battery may be hot after use in the machine.**

- Remove the battery with the battery removal device and let it cool.

The battery removal aid is to safely replace the battery intraoperatively in compliance with sterile conditions.

- Turn the product 1 with the battery shaft 4 upwards.
- Press both safety catches 9 on the lid 8 at the same time and remove the lid 8.
- Attach the sterile battery removal aid 10, see Fig. B.
- Gently shake the product 1 with the attached battery removal aid 10 with battery compartment 4 facing downwards.
- The battery 7 slides easily into the battery removal aid 10.
- Pass on the battery removal aid 10 including dead battery 7 to a non-sterile person.
- Insert a charged battery 7, see Inserting the rechargeable battery.

### 6.1.4 Removing the rechargeable battery

**⚠ WARNING**

**Risk of injury from hot battery!**

**Battery may be hot after use in the machine.**

- Let the battery cool in the machine and then remove.
- or –
- Remove the battery with the battery removal device and let it cool.

**⚠ CAUTION**

**Damage to the battery caused by knocking it against hard objects!**

- Only remove battery by tapping the product on the palm of the hand.

**⚠ CAUTION**

**Damage to, or destruction of the batteries due to processing!**

- Do not sterilize batteries.

Rechargeable battery must be removed after each surgical procedure and prior to processing of the device.

**Note**

*The battery removal aid can be used to remove the battery more easily, see Fig. B.*

- Turn the product 1 with the battery shaft 4 upwards.
- Press both safety catches 9 on the lid 8 at the same time and remove the lid 8.
- Tightly grasp the battery shaft 4 at the lower end.
- Tap the battery shaft 4 on the palm of the hand until the battery 7 slides out of the battery shaft 4 and can be removed.

### 6.1.5 Protection against inadvertent activation

To prevent inadvertent activation of the product during tool change/alignment of the saw head 12, the trigger 2 can be locked.

To lock trigger 2:

- Twist the safety catch 3 to position OFF.

The trigger 2 is blocked and the product 1 cannot be operated.

To unlock trigger 2:

- Twist the safety catch 3 to position ON.

The trigger 2 is unlocked and the product 1 can be operated.

## 6.1.6 Attaching and detaching the saw blade

### ⚠ WARNING

Risk of injury when attaching/removing sawblades in the ON position through inadvertent activation of the product!

- ▶ Only attach/remove saw blades in the OFF position.

### Coupling the saw blade

- ▶ Secure product 1 from inadvertent activation with the safety catch 3, see Protection against inadvertent activation.
- ▶ Insert the saw blade 11 with its coupling end in the slot of saw blade coupling 14 as far as it will go. Ensure that the locking pin engages in the saw blade fenestration, and that the lateral stops of the saw blade rest on the coupling, see Fig. C.
- ▶ If necessary press button for saw blade release/fixation 13.

### Uncoupling the saw blade

- ▶ Secure product 1 against inadvertent activation with the safety catch 3. see Protection against inadvertent activation
- ▶ Fully press button for saw blade release/fixation 13.
- ▶ Pull out saw blade 11 from the saw blade coupling 14, see Fig. C.

### Aligning the saw head

Saw head 12 can be turned to eight different positions, in steps of 45°.

- ▶ Secure product 1 from inadvertent activation with the safety catch 3, see Protection against inadvertent activation.
- ▶ Slightly pull out saw head 12 and turn it to the required position (in steps of 45°).

Saw head 12 noticeably engages at each step, see Fig. D.

## 6.2 Function checks

The function checks must be carried out prior to each use and after each intraoperative battery change.

- ▶ Make certain there is a battery in the battery compartment.
- ▶ Ensure that the lid has completely clicked into place.
- ▶ Pull on the saw blade to check that it is securely coupled.
- ▶ Check that the cutting edges of the saw blade do not exhibit any mechanical damage.
- ▶ Check secure seating of the saw head: turn without pulling out.
- ▶ Release product for use (ON position).
- ▶ Briefly operate product at maximum oscillation frequency.
- ▶ Do not use the product if it is damaged or defective.
- ▶ Set aside the product if it is damaged.

## 6.3 Safe operation

### ⚠ WARNING

Coagulation of patient tissue or risk of burns for patients and user through hot product!

- ▶ Cool the tool during operation.
- ▶ Keep product/tool out of the reach of patients.
- ▶ Let the product/tool cool down.
- ▶ Use a cloth to protect against burns when changing the tool.

### ⚠ WARNING

Risk of infection from aerosol formation!

Risk of injuries caused by particles being released from the tool!

- ▶ Use suitable protection (such as waterproof protective clothing, face mask, safety gases, suction).

### ⚠ WARNING

Risk of injury and/or malfunction!

- ▶ Always carry out a function check prior to using the product.

### ⚠ WARNING

Risk of injury when using the product beyond the field of view!

- ▶ Apply the product only under visual control.

### ⚠ WARNING

Risk of injury and damage to the tool/system!

The tool may get caught in drapes (such as textiles).

- ▶ Do not let the tool come into contact with drapes (such as textiles) during operation.

### ⚠ WARNING

Risk of injury if the saw blade is released unintentionally!

- ▶ Do not press the button to unlock the saw blade during operation.
- ▶ Check that the tool is securely in position each time it is changed.

### Note

*A magnetic sensor system is used to operate the drive motor of the product. To prevent the motor being activated inadvertently, the product should not be exposed to any magnetic fields (such as magnetic instrument pads).*

### 6.3.1 Operating the product

#### Note

*The slight whistling noise when the product is starting up is due to the construction.*

- ▶ Activate trigger 2.

The oscillation frequency of the product 1 is continuously variable.

## 7. Validated reprocessing procedure

### 7.1 General safety notes

#### Note

*Adhere to national statutory regulations, national and international standards and directives, and local, clinical hygiene instructions for reprocessing.*

#### Note

*For patients with Creutzfeldt-Jakob disease (CJD), suspected CJD, or possible variants of CJD, observe the relevant national regulations concerning the reprocessing of products.*

#### Note

*Mechanical reprocessing should be favored over manual cleaning as it gives better and more reliable results.*

#### Note

*It should be noted that successful reprocessing of this medical device can only be guaranteed following prior validation of the reprocessing method. The operator/reprocessing technician is responsible for this.*

#### Note

*If there is no final sterilization, then a virucidal disinfectant must be used.*

#### Note

*See also B.Braun eifu at [eifu.bbraun.com](http://eifu.bbraun.com) for current information on processing and compatibility of material*

*The validated steam sterilization procedure was carried out in the Aesculap sterile container system.*

### 7.2 General information

Dried or affixed surgical residues can make cleaning more difficult or ineffective and lead to corrosion. Therefore the time interval between application and processing should not exceed 6 h; also, neither fixating pre-cleaning temperatures >45 °C nor fixating disinfecting agents (active ingredient: aldehydes/alcohols) should be used.

Excessive neutralizing agents or basic cleaners may result in a chemical attack and/or fading and the laser marking becoming unreadable either visually or by machine.

On stainless steel, residues containing chlorine or chloride (such as surgical residues, drugs, saline solutions in water for cleaning, disinfection and sterilization) may lead to corrosion (pitting corrosion, tensile corrosion) and thus to the destruction of the product. These must be removed by rinsing thoroughly with demineralized water and then drying.

Perform additional drying, if necessary.

Only process chemicals that have been tested and approved (e.g. VAH or FDA approval or CE mark) and which are compatible with the product's materials according to the chemical manufacturers' recommendations may be used for processing the product. All the chemical manufacturer's application specifications must be strictly observed. Failure to do so can result in the following problems:

- Material damage (such as corrosion, cracks, breaks, premature aging or swelling).
- Do not use metal cleaning brushes or other abrasives that would damage the product surface and could cause corrosion
- Further detailed advice on hygienically safe and material-/value-preserving reprocessing can be found at [www.a-k-i.org](http://www.a-k-i.org) link to "AKI-Brochures", "Red brochure".

### 7.3 Reusable products

The service life of the product is limited by damage, normal wear, type and duration of the application, handling, storage and transportation of the product.

Careful visual and functional testing prior to next use is the best way to identify a malfunctioning product.

### 7.4 Preparations at the place of use

- Remove all attached components from the product (tool an accessories).
- Remove any visible surgical residues as much as possible with a damp, lint-free cloth.
- Place the dry product in a sealed waste container and forward it on for cleaning and disinfection within 6 hours.

### 7.5 Preparation before cleaning

- Before the machine is cleaned/disinfected for the first time: Assemble ECCOS holder in a suitable screen basket (e.g. GB243800).
- Insert products correctly into ECCOS brackets, see Fig E.

### 7.6 Product-specific safety instructions for the processing procedure

#### ⚠ CAUTION

**Damage to the product due to inappropriate cleaning/disinfecting agents and/or excessive temperatures!**

- Use cleaning and disinfecting agents according to the manufacturer's instructions which
  - are approved for plastic material and high-grade steel,
  - and which do not affect softeners (e.g., in silicone).
- Do not use cleaning agents that contain acetone.
- Observe information concerning concentration, temperature and exposure time.
- Do not exceed the maximum temperature of 60 °C in chemical cleaning and/or disinfection.
- Do not exceed maximum temperature during thermal disinfection using 96 °C DI water.
- Dry the product for at least 10 minutes at a maximum of 120 °C.

#### ⚠ CAUTION

**Damage or destruction of the batteries due to processing!**

- Protect battery from moisture.

#### Note

*The drying time listed is only indicative. It has to be checked taking into account the specific conditions (e.g. load) and adjusted if necessary.*

## 7.7 Manual cleaning with wipe disinfection

Phase	Step	T [°C/°F]	t [min]	Conc. [%]	Water qual- ity	Chemicals
I	Pre-cleaning	RT (cold)	≥2	-	D-W	Until visually clean
II	Cleaning with enzyme solution	RT (cold)	≥2	0.8	D-W	pH-neutral*
III	Intermediate rinse	RT	≥5	-	D-W	-
IV	Drying	RT	-	-	-	-
V	Wipe disinfection	-	>1	-	-	Meliseptol HBV wipes 50 % Propan-1-ol
VI	Final rinse	RT (cold)	0.5	-	FD-W	-
VII	Drying	RT	-	-	-	-

D-W: Drinking water

FD-W: Fully de-ionized water (de-mineralized, microbiological, with minimum quality of drinking water)

RT: Room temperature

\* Suitable enzyme solution: Helizyme, Cidezyme (the latter for validation purposes)

- Do not clean the product in a ultrasonic bath and do not immerse the product in any fluids. Let any fluid incursions drain out immediately, otherwise there is a danger of corrosion and loss of function.

### Phase I

- Mobilize non-rigid components during cleaning.
- Clean the product under running water, using a suitable cleaning brush until all visible residues have been removed from the surfaces.
- Brush difficult to access surfaces with a suitable plastic cleaning brush for at least 1 min.

### Note

For details on difficult to access surfaces, see Acculan pre-cleaning and care information TA016000 (available in Aesculap Extranet at <https://extranet.bbraun.com>).

### Phase II

- Follow the operating instructions of the enzyme cleaner with regard to correct concentration, dilution, temperature and water quality.
- Spray products with a pH neutral enzyme solution, let soak in for at least 2 minutes and then wipe off.
- Contamination should be removed with a lint-free cloth or soft brush moistened with enzyme cleaner.
- Flexible components should be rinsed for 20 seconds with the water pistol (cold water, at least 2.5 bar).
- After manual cleaning, check visible surfaces and areas of flexible components for residues.

### Phase III

- Mobilize non-rigid components during cleaning.
- Rinse product under running tap water for at least 5 minutes.
- If necessary, repeat the cleaning process (phase 1 to 3).

### Phase IV

- Dry the product in the drying phase with suitable equipment (such as lint-free cloths, pressurized air).

### Phase V

- Wipe all surfaces of the product with a single-use disinfecting wipe.

### Phase VI

- Rinse disinfected surfaces after the prescribed reaction time for at least 1 minute under running demineralized water.
- Drain any remaining water fully.

### Stage VII

- Dry the product in the drying phase with suitable equipment (such as lint-free cloths, pressurized air).

## 7.8 Automatic cleaning/disinfection with manual pre-cleaning

### Note

The cleaning and disinfection device must fundamentally have a tested efficacy (such as FDA approval or CE label pursuant to DIN EN ISO 15883).

### Note

The cleaning and disinfection machine used for processing must be serviced and checked at regular intervals.

### 7.8.1 Manual pre-cleaning with a brush

Phase	Step	T [°C/°F]	t [min]	Conc. [%]	Water quality	Chemistry/Note
I	Rinsing	RT (cold)	-	-	D-W	until visually clean
II	Brushes	RT (cold)	-	-	D-W	until visually clean

D-W: Drinking water

RT: Room temperature

- ▶ Do not clean the product in a ultrasonic bath and do not immerse the product in any fluids. Let any fluid incursions drain out immediately, otherwise there is a danger of corrosion and loss of function.

### Phase I

- ▶ Mobilize non-rigid components during cleaning.
- ▶ Thoroughly clean the product under running water.

### Phase II

- ▶ Mobilize non-rigid components during cleaning.
- ▶ Brush difficult to access surfaces with a suitable plastic cleaning brush for at least 1 min.
- ▶ After manual preparation, check visible surfaces for residue and repeat the pre-cleaning process as needed.

### Note

For details on difficult to access surfaces, see Acculan pre-cleaning and care information TA016000 (available in Aesculap Extranet at <https://extranet.bbraun.com>).

## 7.8.2 Mechanical alkaline cleaning and thermal disinfection

Machine type: single-chamber cleaning/disinfection device without ultrasound

Phase	Step	T [°C/°F]	t [min]	Water quality	Chemical
I	Pre-rinse	<25/77	3	D-W	-
II	Cleaning	55/131	10	FD-W	<ul style="list-style-type: none"> <li>■ Concentrate, alkaline:           <ul style="list-style-type: none"> <li>- pH ~ 13</li> <li>- &lt;5 % anionic surfactant</li> </ul> </li> <li>■ 0.5 % working solution           <ul style="list-style-type: none"> <li>- pH ~ 11*</li> </ul> </li> </ul>
III	Intermediate rinse	>10/50	1	FD-W	-
IV	Thermal disinfection	90/194	5	FD-W	-
V	Drying	max. 120/248	min. 10min	-	-

DW: Drinking water

FD-W: Fully desalinated water (demineralized, low microbiological contamination: drinking water quality at least)

\*Recommended: BBraun Helimatic Cleaner alkaline

- ▶ Insert the product in its correct position into the ECCOS holder.
- ▶ Set the saw head **12** between two locking positions for cleaning.
- ▶ After automatic cleaning/disinfection, check visible surfaces for residues and repeat the cleaning/disinfection process as needed.

## 7.9 Inspection, maintenance and checks

- ▶ Allow the product to cool down to room temperature.

### Note

Aesculap additionally recommends occasionally spraying moving parts (such as knobs, coupling, lid caps) with STERILIT Power Systems oil spray.

- ▶ Inspect the product after every cleaning and disinfecting cycle for: cleanliness, damage, function, abnormal operation noise, excessive heat or strong vibration.
- ▶ Inspect the saw blades for broken, damaged or blunt cutting edges.
- ▶ Set aside the product if it is damaged.

## 7.10 Packaging

- ▶ Observe the instructions for use of the respective packaging and holders used (e.g. Instructions for use TA009721 for Aesculap ECCOS holder system).
- ▶ Insert the products in the correct position into the ECCOS holders, see Fig. E.
- ▶ Pack trays appropriately for the sterilization process (e.g. in Aesculap sterile containers).
- ▶ Ensure that the packaging will prevent a recontamination of the product.

## 7.11 Steam sterilization

### ⚠ CAUTION

Damage or destruction of the batteries due to processing!

- ▶ Do not sterilize the rechargeable battery.

### Note

Remove all attached components from the product (tools, accessories) before sterilization.

- ▶ Check to ensure that the sterilizing agent will come into contact with all external and internal surfaces (e.g. by opening any valves and faucets).
- ▶ Use a validated sterilization method:
  - Steam sterilization in fractionated vacuum process
  - Steam sterilizer in accordance with DIN EN 285 and validated in accordance with DIN EN ISO 17665
  - Sterilization in fractionated vacuum process at 134°C, holding time 5 min

When sterilizing multiple products in one steam sterilizer:

- ▶ Ensure that the maximum permissible load of the steam sterilizer according to the manufacturer's specifications is not exceeded.

## 7.12 Storage

- ▶ Store sterile products in germ-proof packaging, protected from dust, in a dry, dark, temperature-controlled area.

## 8. Maintenance

To ensure reliable operation, the product must be maintained in accordance with the maintenance labeling or at least once a year.

YYYY-MM

For technical service, please contact your national B. Braun/Aesculap agency, see Technical Service.

## 9. Troubleshooting list

- Have defective products repaired by Aesculap Technical Service, see Technical Service.

Issue	Cause	Detection	Remedy
Product not running	No battery	No battery in the battery compartment	Insert battery.
	Battery not charged	No beep when inserting the battery	Charge battery in charger.
	Battery defective	No signal tone when inserting the battery	Have the battery repaired by the manufacturer.
	Product in secured OFF position	Safety catch is in OFF position	Turn the safety catch to the ON position.
	Product defective	Product not running	Have the product repaired by the manufacturer.
Product is getting too hot	Excessive use	Product is heating up	Observe instructions for use (nominal operating mode).
	Reprocessing/maintenance carried out incorrectly	Product heating	Follow operating instructions (preparation, care). Preventative: lubricate product before every sterilization.
	Fall damage, product defective	Product heating	Have the product repaired by the manufacturer.
	Blunt tool	Tool and product are heating up	Change tool.
Insufficient power	Product defective	Inadequate product power	Follow operating instructions (preparation, care). Preventive measure: Lubricate product before every sterilization. Have the product repaired by the manufacturer.
		Large amount of heat generated after a short time	Follow operating instructions (operating mode). Have the manufacturer repair the product
	Blunt tool	Tool blades worn out	Replace tool.
Loud running noise	Defective product drive/ball bearings	Loud, perceptible noise during operation	Observe instructions for use (processing, care). Preventive measure: Lubricate product before every sterilization. Have the manufacturer repair the product
Lid cannot be attached/removed	Lid not compatible	Lid does not click into place	Use the compatible lid for GA331.
	Lid deformed/defective	Lid is difficult or cannot be mounted/removed	Have the manufacturer repair the lid.
	Connection to product deformed/defective	Lid is hard or impossible to fit/remove	Have the product repaired by the manufacturer.
Battery cannot be inserted/removed	Battery not compatible	Battery cannot be completely inserted into the battery compartment	Use the compatible battery for GA331.
	Battery deformed/defective	Battery is difficult to insert/remove or cannot be inserted/removed	Have the manufacturer repair the battery.
	Battery compartment on product deformed/defective	Battery is difficult or cannot be mounted/removed	Have the product repaired by the manufacturer

Issue	Cause	Detection	Remedy
Sterile funnel or battery removal aid cannot be plugged in	Sterile funnel or battery removal device not compatible	Sterile funnel or battery removal aid not plugged into battery compartment	Use suitable sterile funnel or battery removal aid for GA331.
	Sterile funnel or battery removal device deformed/defective	Sterile funnel or battery removal device difficult or not attachable	Replace sterile funnel or battery removal aid.
	Battery compartment on product deformed/defective	Difficult or impossible to plug in the sterile funnel or battery removal aid	Have the product repaired by the manufacturer.
Speed control knob cannot be actuated	Product in locked OFF position	Safety catch is in OFF position	Turn the safety catch into the ON position.
	Speed control knob jammed/defective	Speed control knob not working	Have the manufacturer repair the product.
Saw blade not moving	Loud running noise from motor	Broken gear unit	Have repaired by manufacturer
Insufficient cutting performance of the saw blade	Saw blade teeth worn	Blunt saw blade	Replace saw blade
	Battery too weak	Insufficient motor power/speed	Charge battery
	Battery worn/defective		Replace battery
	Saw blade overheated	Chips not removed adequately	When sawing, move back and forth more and rinse sufficiently
The saw blade cannot be coupled	Coupling components or saw blades deformed	Deformation	Have repaired by manufacturer
	-	Incompatible saw blade	Insert a compatible Aesculap saw blade, see Accessories/Spare parts
Saw blade cannot be uncoupled	Locking pin not fully releasing the saw blade	Push button for unlocking saw blade not fully pushed	Fully press the saw blade release button, see Attaching and detaching the saw blade

## 10. Technical Service

### ⚠ DANGER

Danger to life of patients and users in case of malfunctions and/or failure of protective measures!

- Do not perform any servicing or maintenance work under any circumstances while the product is being used on a patient.

### ⚠ CAUTION

Modifications of medical devices may result in a loss of potential guarantee/warranty claims and forfeiture of applicable licenses.

- Do not modify the product.
- Contact national B. Braun/Aesculap representative for service and repair.

### Service addresses

Aesculap Technischer Service

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78532 Tuttlingen / Germany

Phone: +49 (7461) 95 -1601

Fax: +49 (7461) 14 -939

E-Mail: ats@aesclap.de

Other service addresses can be obtained from the address indicated above.

## 11. Accessories/Spare parts

Art. no.	Name
GA675	Lid
GA676	NiMH battery long
GA678	Sterile funnel
GA679	Battery removal aid
GB243R GB244R	ECCOS screen basket with holder for Acculan 4
GB600	STERILIT Power Systems oil spray
GB495R	ECCOS fixation for machine
GB487R	ECCOS holder for lid
GB488R	ECCOS holder for sterile funnel
GB489R	ECCOS holder for battery removal aid
TA014542	Instructions for use for oscillating saw GA331 (A4 for ring binder)
TA014543	Instructions for use for oscillating saw GA331 (flyer)

## 12. Technical data

### 12.1 Classification according to Regulation (EU) 2017/745

Art. no.	Designation	Class
GA331	Oscillating saw	IIa

### 12.2 Performance data, information about standards

Max. power	approx. 250 W
Max. oscillation frequency	13 000 min <sup>-1</sup>
Weight (ready for operation)	1,6 kg ±10 %
Dimensions (L x W x H, ready for operation)	211 mm x 177 mm x 55 mm ±5 %
Application part	Type BF
EMC	IEC/EN 60601-1-2
Conforming to standard	IEC/EN 60601-1

After 500 preparatory cycles, the product was tested by the manufacturer and passed.

### 12.3 Operating mode

Operation with non-periodic load and speed changes (type S9 pursuant to IEC EN 60034-1)

- 30 second application, 30 second pause
- 7 repetitions
- 30 min cooling time
- Max. Temperature 48 °C

### 12.4 Ambient conditions

	Operation	Transport and storage
Temperature	10 °C to 27 °C	-10 °C to 50 °C
Relative humidity	30 % to 75 %	10 % to 90 %
Atmospheric pressure	700 hPa to 1 060 hPa	500 hPa to 1 060 hPa

## 13. Disposal

### ⚠ WARNING

Risk of infection from contaminated products!

- Observe national regulations when disposing of or recycling the product, its components, and their packagings.

### Note

The operator must process the product before disposal, see *Validated reprocessing procedure*.

The recycling pass can be downloaded from the Extranet as a PDF document under the respective article number. (The recycling pass includes disassembling instructions for the product, as well as information for proper disposal of components harmful to the environment.)

Products carrying this symbol are subject to separate collection of electrical and electronic devices. Within the European Union, disposal is taken care of by the manufacturer as a free-of-charge service.

- For questions regarding the disposal of the product, please contact your local B. Braun/Aesculap representative, see Technical Service.

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Aesculap AG | Am Aesculap-Platz | 78532 Tuttlingen | Germany