Histoacryl®

D	Gebrauchsanweisung	- 3	I	Istruzioni per l'uso	- 46
GB	Instructions for use	- 7	GR	Οδηγίες χρήσης	- 50
CN	使用说明 组织胶水 蓝色	- 11	CZ	Návod k použití	- 54
F	Mode d'emploi	- 14	PL	Instrukcja użytkowania	- 58
NL	Gebruiksaanwijzing	- 18	TR	Kullanım Kılavuzu	- 62
E	Instrucciones de uso	- 22	SK	Návod na použitie	- 66
P	Instruções de uso	- 26	BG	Инструкции за употреба	- 70
DK	Bruksanvisning	- 30	HU	Használati utasítás	- 74
S	Bruksanvisning	- 34	RO	Instrucţiuni de utilizare	- 78
FIN	Käyttöohje	- 38	RU	Инструкция для применения	- 82
N	Bruksanvisning	- 42	CN	使用说明 组织胶水 无色	- 88

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- **Maluf-Filho F., Sakai P., Ishioka S., Matuguma S.E. (2001)**
Endoscopic sclerotherapy versus cyanoacrylate endoscopic injection for the first episode of variceal bleeding: a prospective, controlled, and randomized study in Child-Pugh class C patients. *Endoscopy*. 2001; 33(5): 421-7.
- **Binmoeller K.F., Soehendra N. (1995)**
Superglue: the answer to variceal bleeding and fundal varices? *Endoscopy*. 1995; 27: 392-6.
- **Gotlib J.-P. (1990)**
Endoscopic obturation of esophageal and gastric varices with a cyanoacrylate tissue adhesive. *Can J Gastroenterol*. 1990; 4(9): 637-8.
- **Seewald S., Tiing Leong A., Imazu H., Naga M., Omar S., Groth S., Seitz U., Yan Zhong, Thonke F., Soehendra N. (2008)**
A standardized injection technique and regimen ensures success and safety of N-butyl-2-cyanoacrylate injection for the treatment of gastric fundal varices. *Gastrointest Endosc*. 2008, 68(3): 447-54.

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Histoacryl®

Description

The tissue adhesive Histoacryl® consists of enbucrilate. For easy assessment of the thickness of the layer which has been applied Histoacryl® is coloured blue with the dye D&C violet No. 2. A translucent version, which does not contain dye, is also available.

Ingredients:

Enbucrilate 1 ampoule of 0.5 ml contains 0.5 g

Indications

- Closure of minimum-tension wounds from clean surgical incisions and simple, thoroughly cleansed, trauma-induced lacerations
- Sclerotherapy of large oesophageal or fundal varices
- Fixation of hernia meshes, especially in inguinal hernia surgery

Mode of application and dosage for closure of smooth and fresh skin wounds

Open the aluminium pouch and take out the ampoule containing the adhesive, in sterile conditions. Hold the ampoule by its cannula, upside down, and shake it vigorously to ensure that no adhesive remains within the cannula.

Then open the ampoule by twisting off the ribbed tip of the cannula. It is advisable to hold the yellow part of the ampoule between two fingers, holding it vertically with the thin end pointing upwards, while opening the ampoule. This will prevent any leakage of Histoacryl® from the ampoule on opening.

Appose tissue edges with forceps and hold in apposition while applying Histoacryl®, for approximately 30 seconds, to allow Histoacryl® to cure and to prevent seepage between wound edges.

Unless otherwise prescribed, as little Histoacryl® tissue adhesive as possible should be applied; the amount applied is already sufficient if slight colouring is visible. A sparse application of thin layers or spots is required for undisturbed wound healing. Heavy application may cause thermal damage to tissues, and delayed healing may result. If too much adhesive is accidentally applied it can be removed within the first few seconds using a dry swab.

Mode of action for closure of smooth and fresh skin wounds

In the presence of tissue moisture Histoacryl® immediately polymerises into a solid substance which attaches firmly to the tissue.

Mode of application for sclerotherapy of large oesophageal or fundal varices

Open the aluminium pouch and take out the ampoule containing the adhesive, in sterile conditions. Hold the ampoule by its cannula, upside down, and shake it vigorously to ensure that no adhesive remains within the cannula. Then open the ampoule by twisting off the ribbed tip of the cannula. For convenient transfer of the adhesive into a syringe the opening of the cannula can be enlarged by cutting the cannula with scissors. It is advisable to hold the yellow part of the ampoule between two fingers, holding it vertically, with the thin end pointing upwards, while opening the ampoule. This will prevent any leakage of Histoacryl® from the ampoule on opening.

Sclerotherapy with Histoacryl® should only be performed by a physician who is familiar with its application for this indication. The method of application depends on the user's individual experience in endoscopic varix therapy with cyanoacrylate adhesives. It is described in detail in the relevant literature, for instance in the following articles:

- **Maluf-Filho F., Sakai P., Ishioka S., Matuguma S. E. (2001)**

Endoscopic sclerosis versus cyanoacrylate endoscopic injection for the first episode of variceal bleeding: a prospective, controlled, and randomized study in Child-Pugh class C patients. *Endoscopy*. 2001; Vol 33 No 5: 421-427.

- **Binmoeller K. F., Soehendra, N. (1995)**

Superglue: the answer to variceal bleeding and fundal varices?

Endoscopy. 1995; Vol. 27: 392-396.

- **Gotlib J.-P. (1990)**

Endoscopic obturation of esophageal and gastric varices with a cyanoacrylic tissue adhesive. *Can J Gastroenterol* 1990; Vol 4 No 9: 637-638.

- **Seewald S., Tiing Leong A., Imazu H., Naga M., Omar S., Groth S., Seitz U., Yan Zhong, Thonke F., Soehendra N. (2008)**

A standardized injection technique and regimen ensures success and safety of N-butyl-2-cyanoacrylate injection for the treatment of gastric fundal varices.

Gastrointestinal endoscopy, 2008, Vol. 68, no 3: 447-454.

The application of Histoacryl® should be limited mainly to large oesophageal or stomach fundus varices. In cases of small-calibre varices, administration must be carried out strictly intravascularly, since paravariceal injections can cause consider-

able tissue ulceration.

In sclerotherapy, the success of the treatment depends on immediate haemostasis and on preventing recidivist haemorrhaging caused by the removal of the varices. The limiting factor in this is that for some weeks there is a high risk of recidivist haemorrhaging, until the varices are completely thrombosed. The idea behind using cyanoacrylates is first to induce a rapid eradication of the varices.

Mode of action: sclerotherapy of large oesophageal or fundal varices

Varically injected Histoacryl® polymerises intravascularly to form a plastic cylinder and effects immediate obturation or thrombosing of the vessel. However, there is also an inflammatory reaction, which then leads to fibrosis. The mechanical obliteration and fibrosis induced in this way act synergically in cases of oesophageal or stomach fundus varices and reduce the risk of recidivist haemorrhage. Following injection into submucosal varices, the Histoacryl® cylinders are expelled from the oesophagus/stomach wall into the lumen within 1-4 months.

Mode of application for the fixation of meshes in the body

Open the aluminium pouch and take out the ampoule containing the adhesive, in sterile conditions. Hold the ampoule by its cannula, upside down, and shake it vigorously to ensure that no adhesive remains within the cannula. Then open the ampoule by twisting off the ribbed tip of the cannula.

It is advisable to hold the yellow part of the ampoule between two fingers, holding it vertically, with the thin end pointing upwards, while opening the ampoule. This will prevent any leakage of Histoacryl® from the ampoule on opening.

Mesh fixation with Histoacryl® should only be performed by a physician who is familiar with its application for this indication. The method of application depends on the user's individual experience in endoscopic or open mesh fixation with cyanoacrylate adhesives. It is described in detail in the relevant literature, for instance in the following article:

- **Hebeling C., Schlumpf R. (2003)**
Sutureless Lichtenstein: first results of a prospective randomised trial
Hernia 2003; Vol 7 No 2: 80-84.

Mode of action: fixation of meshes in the body

The adhesive should be applied in drops, without allowing direct contact of the application tip with the tissue, so that the adhesive does not harden

on the tip and remains applicable for the duration of the surgery. When applied on the placed and adapted mesh the adhesive flows through the mesh and hardens upon contact with the underlying tissue. The mesh is fixed by the hardened adhesive. The adhesive hardens within 5-7 seconds. After 3 minutes it has hardened completely and its surface is no longer adherent.

In Lichtenstein repair the mesh has to be slit. If it is necessary to make a strip of mesh as narrow as 2 cm or less, the mesh has to be fixed with one or two stitches in this area because fixation with adhesive is not secure in this case.

Contraindications

Histoacryl® is not indicated for holding together wound surfaces on internal organs or on the brain surface, or in the central or peripheral nervous system, as such action could cause tissue damage and scarring, with the resultant problems. Application to the intima and media of blood vessels must also be avoided, since this would involve the risk of thrombosis and vessel wall damage. In sclerosation of oesophageal and fundal varices, these effects are desired. Histoacryl® is not to be applied to areas with dense natural hair, or within the conjunctival sac of the eye. Histoacryl® is not to be applied to the surface of the eye. Histoacryl® is not to be applied to wounds showing evidence of infection or gangrene, or to decubitus ulcers. Histoacryl® is not to be used on patients with known preoperative systemic infections, uncontrolled diabetes, or diseases or conditions that are known to interfere with the wound healing process. Histoacryl® is not indicated on patients with a known hypersensitivity to cyanoacrylate, formaldehyde, or the dye D&C Violet No. 2.

Warning note

During closure of smooth and fresh skin wounds Histoacryl® must not be introduced into the wound, since this would interfere with wound healing.

Following varix obliteration with Histoacryl®, haemorrhaging from necroses/ulcerations, bacteraemia, fever and chronically scarred oesophageal stenoses can occur in rare cases. Furthermore, it cannot be entirely ruled out that Histoacryl® embolisations may, for example, get washed into the lungs, potentially resulting in pulmonary embolism. The use of a dissolution medium can lead to early polymerisation. However, this must not be attributed to the quality of Histoacryl®; rather, it depends on the mixing ratio, water content and other factors in connection with the dissolution medium.

Only product that is liquid in consistency must be used. The condition of the adhesive can be judged with the plastic ampoule unopened.

In mesh fixation only the amount of adhesive necessary for secure fixation should be used. Heavy application may cause thermal damage to tissues, and delayed healing may result. Coverage of complete areas of the reinforcement material is not necessary. If the application tip is used for pushing the mesh, the tip must be removed immediately after application of a Histoacryl® drop on the mesh. Otherwise the application tip may become adhered to the mesh or the tissue.

When using cyanoacrylates as fixation material, the peritoneum must be closed in order to avoid contact of the adhesive with the bowel or any other intra-abdominal organ. Cyanoacrylates, like other fixation materials, may induce the formation of adhesions.

Precautionary measures on application

During closure of smooth and fresh skin wounds it should be ensured that the edges of the wound have been thoroughly cleaned, debrided and properly apposed before applying Histoacryl®. Local anaesthesia might be necessary. Histoacryl® generates a small amount of heat during polymerisation and should not be applied to tissues that may be affected by such heat. Histoacryl® may be used in conjunction with, but not in place of, dermal sutures. Assure haemostasis, close the dermis as needed, and assure that surface edges are easily apposable before applying Histoacryl®. For wounds at risk of tension, before applying Histoacryl® to the skin surface relieve potential stress along the wound line by approximating the wound edges with subcuticular sutures. Use of Histoacryl® may result in localised sensitisation or irritation reactions. The tissue areas to be joined can be coated directly by cautiously pressing the wall of the plastic ampoule. Particularly fine control of compression can be achieved by gripping the plastic ampoule in the jaws of a pair of blunt forceps. The areas to be joined should be as dry as possible. The skin surfaces should be held together for 30 seconds after application of Histoacryl®. It is not possible to make corrections once the adhesive has hardened. Wounds that are under tension should be secured with holding sutures; wounds longer than 3 cm should also be secured by sutures.

The user must be familiar with the particular properties of a liquid tissue adhesive when applying Histoacryl® tissue adhesive.

When used for varix obliteration, liquid Histoacryl® can damage the endoscope through adhesion. Prior to applying Histoacryl®, the biopsy duct of the endoscope should be irrigated with silicone oil.

If Histoacryl® gets into the eye by accident, adhesion can be avoided by rinsing with water immediately. If adhesions or substance depositions still occur, especially in the eye, these will dissolve, without further action, after some days. Do not force the eye open.

Accidental bonding of unwanted skin may occur. Do not pull the skin apart. If necessary, adhesions in less sensitive skin areas can be carefully dissolved with acetone or by soaking in warm water until the skin can be separated.

Care must be taken to ensure that instruments, cloths, swabs or gloves do not come into direct contact with the adhesive applied, as they might adhere to the tissue.

Histoacryl® must not be introduced into the wound, since this could lead to long-term blue staining of the tissue. For facial application it is recommended to use undyed Histoacryl®. Histoacryl® has not been evaluated in patients with a history of hypertrophic scarring or keloid formation.

Ampoules which get clogged, i.e. during mesh fixation, can be reopened with a sterile scalpel or scissors. The physician should take care with regard to possible sharp edges on the reopened ampoule, and should avoid any cut-offs or flakes becoming transferred to the site.

Side effects

The use of this product leads to an exothermic reaction. During the closure of smooth and fresh skin wounds the inappropriate application of too thick a layer of adhesive can lead, upon polymerisation, to thermal damage of the tissue. Connective tissue healing can be impeded when too much tissue adhesive either in thickness or in area is applied. Cyanoacrylates can be associated with a limited period of local irritation in the application area; a transient foreign body reaction can occasionally take the form of an inflammatory reaction.

Sterilisation

The sterility of Histoacryl® tissue adhesive is assured by membrane filtration and aseptic filling procedures.

Additional notes

Instruments that are contaminated with Histoacryl® can be cleaned with dimethylformamide or acetone.

Histoacryl® is intended for single use and should be applied immediately after opening the ampoule. Opened unused ampoules should be discarded. Do not re-use: Infection hazard for patients and/or users and impairment of product functionality due to re-use. Risk of injury, illness or death due to contamination and/or impaired functionality of the product.

Histoacryl® should be stored at ambient temperature below +22 °C. The ampoule containing the adhesive should only be removed from the aluminium pouch immediately prior to application.

In order to discard any remaining Histoacryl® without the risk of contaminating the surroundings, it is recommended that, prior to discarding the ampoule, the ribbed tip which was removed on opening the ampoule be replaced over the narrow end of the ampoule.

Histoacryl® must not be re-sterilised.

Symbols used on labelling

DIM

Volume

REF

Cat. No.

LOT

Batch Number

Use until Year + Month

STERILE A

Sterile unless package is opened or damaged.

Method of sterilisation: aseptic filling

Upper limit of temperature

Do not reuse!

See Instructions for Use!

CE mark and identification number of notified body. Product conforms to the essential requirements of the Medical Device Directive 93/42/EEC.

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- **Maluf-Filho F., Sakai P., Ishioka S., Matuguma S. E. (2001)**
Endoscopic sclerosis versus cyanoacrylate endoscopic injection for the first episode of variceal bleeding: a prospective, controlled, and randomized study in Child-Pugh class C patients.
Endoscopy. 2001; Vol 33 No 5: 421-427.
- **Binmoeller K. F., Soehendra, N. (1995)**
Superglue: the answer to variceal bleeding and fundal varices?
Endoscopy. 1995; Vol. 27: 392-396.
- **Gotlib J.-P. (1990)**
Endoscopic obturation of esophageal and gastric varices with a cyanoacrylic tissue adhesive.
Can J Gastroenterol 1990; Vol 4 No 9: 637-638.
- **Seewald S., Tiing Leong A., Imazu H., Naga M., Omar S., Groth S., Seitz U., Yan Zhong, Thonke F., Soehendra N. (2008)**
A standardized injection technique and regimen ensures success and safety of N-butyl-2-cyanoacrylate injection for the treatment of gastric fundal varices.
Gastrointestinal endoscopy, 2008, Vol. 68, no 3: 447-454.

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• Helbling C., Schlumpf R. (2003)

Sutureless Lichtenstein: first results of a prospective randomised trial
Hernia. 2003; Vol 7 No 2: 80-84.

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- **Maluf-Filho F., Sakai P., Ishioka S., Matuguma S. E. (2001)**
Endoscopic sclerosis versus cyanoacrylate endoscopic injection for the first episode of variceal bleeding: a prospective, controlled, and randomized study in Child-Pugh class C patients. *Endoscopy*. 2001; Vol 33 No 5: 421-427.
- **Binnmoeller K. F., Soehendra, N. (1995)**
Superglue: the answer to variceal bleeding and fundal varices? *Endoscopy*. 1995; Vol. 27: 392-396.
- **Gottlieb J.-P. (1990)**
Endoscopic obturation of esophageal and gastric varices with a cyanoacrylic tissue adhesive. *Can J Gastroenterol* 1990; Vol 4 No 9:637-638
- **Seewald S., Tiing Leong A., Imazu H., Naga M., Omar S., Groth S., Seitz U., Yan Zhong, Thonke F., Soehendra N. (2008)**
A standardized injection technique and regimen ensures success and safety of N-butyl-2-cyanoacrylate injection for the treatment of gastric fundal varices. *Gastrointestinal endoscopy*, 2008, Vol. 68, no 3: 447-454.

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- **Maluf-Filho F., Sakai P., Ishioka S., Matuguma S. E. (2001)**
Endoscopic sclerosis versus cyanoacrylate endoscopic injection for the first episode of variceal bleeding: a prospective, controlled, and randomized study in Child-Pugh class C patients. *Endoscopy*. 2001; Vol 33 No 5: 421-427.
- **Binmoeller K. F., Soehendra, N. (1995)**
Superglue: the answer to variceal bleeding and fundal varices? *Endoscopy*. 1995; Vol. 27: 392-396.
- **Gotlib J.-P. (1990)**
Endoscopic obturation of esophageal and gastric varices with a cyanoacrylic tissue adhesive. *Can J Gastroenterol* 1990; Vol 4 No 9: 637-638.
- **Seewald S., Tiing Leong A., Imazu H., Naga M., Omar S., Groth S., Seitz U., Yan Zhong, Thonke F., Soehendra N. (2008)**
A standardized injection technique and regimen ensures success and safety of N-butyl-2-cyanoacrylate injection for the treatment of gastric fundal varices. *Gastrointestinal endoscopy*, 2008, Vol. 68, no 3: 447-454.

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• Helbing E, Schlumpf R. (2003)

Suturless Lichtenstein: first results of a prospective randomised trial
Hernia 2003; Vol 7 No 2: 80-84.

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- **Maluf-Filho F., Sakai P., Ishioka S., Matuguma S. E. (2001)**

Endoscopic sclerosis versus cyanoacrylate endoscopic injection for the first episode of variceal bleeding: a prospective, controlled, and randomized study in Child-Pugh class C patients.

Endoscopy. 2001; Vol 33 No 5: 421-427.

- **Binnmoeller K. F., Soehendra N. (1995)**
Superglue: the answer to variceal bleeding and fundal varices?

Endoscopy. 1995; Vol. 27: 392-396.

- **Gottlieb J.-P. (1990)**
Endoscopic obturation of esophageal and gastric varices with a cyanoacrylic tissue adhesive. Can J Gastroenterol 1990; Vol 4 No 9: 637-638.

- **Seewald S., Tilg Leong A., Imazu H., Naga M., Omar S., Groth S., Seitz U., Yan Zhong, Thonke F., Soehendra N. (2008)**

A standardized injection technique and regimen

ensures success and safety of N-butyl-2-cyanoacrylate injection for the treatment of gastric fundal varices.
Gastrointestinal endoscopy, 2008, Vol. 68, no 3: 447-454.

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- **Maluf-Filho F., Sakai P., Ishioka S., Matuguma S. E. (2001)**
Endoscopic sclerosis versus cyanoacrylate endoscopic injection for the first episode of variceal bleeding: a prospective, controlled, and randomized study in Child-Pugh class C patients.
Endoscopy. 2001; Vol 33 No 5: 421-427.
- **Binmoeller K. F., Soehendra, N. (1995)**
Superglue: the answer to variceal bleeding and fundal varices?
Endoscopy. 1995; Vol. 27: 392-396.
- **Gottlieb J.-P. (1990)**
Endoscopic obturation of esophageal and gastric varices with a cyanoacrylic tissue adhesive.
Can J Gastroenterol 1990; Vol 4 No 9: 637-638.
- **Seewald S., Tiing Leong A., Imazu H., Naga M., Omar S., Groth S., Seitz U., Yan Zhong, Thonke F., Soehendra N. (2008)**
A standardized injection technique and regimen ensures success and safety of N-butyl-2-cyanoacrylate injection for the treatment of gastric fundal varices.
Gastrointestinal endoscopy, 2008, Vol. 68, no 3: 447-454.

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- **Maluf-Filho F., Sakai P., Ishioka S., Matuguma S. E. (2001)**
Endoscopic sclerosis versus cyanoacrylate endoscopic injection for the first episode of variceal bleeding: a prospective, controlled, and randomized study in Child-Pugh class C patients.
Endoscopy. 2001; Vol 33 No 5: 421-427.
- **Binmoeller K. F., Soehendra, N. (1995)**
Superglue: the answer to variceal bleeding and fundal varices?
Endoscopy. 1995; Vol. 27: 392-396.
- **Gotlib J.-P. (1990)**
Endoscopic obturation of esophageal and gastric varices with a cyanoacrylic tissue adhesive.
Can J Gastroenterol 1990; Vol 4 No 9: 637-638.
- **Seewald S., Tiing Leong A., Imazu H., Naga M., Omar S., Groth S., Seitz U., Yan Zhong, Thonke F., Soehendra N. (2008)**
A standardized injection technique and regimen ensures success and safety of N-butyl-2-cyanoacrylate injection for the

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- **Maluf-Filho F., Sakai P., Ishioka S., Matuguma S. E. (2001)**
Endoscopic sclerosis versus cyanoacrylate endoscopic injection for the first episode of variceal bleeding: a prospective, controlled, and randomized study in Child-Pugh class C patients.
Endoscopy. 2001; Vol 33 No 5: 421-427.
- **Binmoeiller K. F., Soehendra, N. (1995)**
Superglue: the answer to variceal bleeding and fundal varices?
Endoscopy. 1995; Vol. 27: 392-396.
- **Gotlib J.-P. (1990)**
Endoscopic obturation of esophageal and gastric varices with a cyanoacrylic tissue adhesive.
Can J Gastroenterol 1990; Vol 4 No 9: 637-638.
- **Seewald S., Tiing Leong A., Imazu H., Naga M., Omar S., Groth S., Seitz U., Yan Zhong, Thonke F., Soehendra N. (2008)**
A standardized injection technique and regimen ensures success and safety of N-butyl-2-cyanoacrylate injection for the treatment of gastric fundal varices.
Gastrointestinal endoscopy, 2008, Vol. 68, no 3: 447-454.

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- **Maluf-Filho F., Sakai P., Ishioka S., Matuguma S. E. (2001)**
Endoscopic sclerosis versus cyanoacrylate endoscopic injection for the first episode of variceal bleeding: a prospective, controlled, and randomized study in Child-Pugh class C patients.
Endoscopy. 2001; Vol 33 No 5: 421-427.
- **Binmoeller K. F., Soehendra, N. (1995)**
Superglue: the answer to variceal bleeding and fundal varices?
Endoscopy. 1995; Vol. 27: 392-396.
- **Gotlib J.-P. (1990)**
Endoscopic obturation of esophageal and gastric varices with a cyanoacrylic tissue adhesive.
Can J Gastroenterol 1990; Vol 4 No 9: 637-638.

- **Seewald S., Tiing Leong A., Imazu H., Naga M., Omar S., Groth S., Seitz U., Yan Zhong, Thonke F., Soehendra N. (2008)**

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- **Maluf-Filho F., Sakai P., Ishioka S., Matuguma S. E. (2001)**
Endoscopic sclerosis versus cyanoacrylate endoscopic injection for the first episode of variceal bleeding: a prospective, controlled, and randomized study in Child-Pugh class C patients.
Endoscopy. 2001; Vol 33 No 5: 421-427.
- **Binmoeller K. F., Soehendra, N. (1995)**
Superglue: the answer to variceal bleeding and fundal varices?
Gastroenterology. 1995; Vol. 27: 392-396.
- **Gotlib J.-P. (1990)**
Endoscopic obturation of esophageal and gastric varices with a cyanoacrylic tissue adhesive.
Can J Gastroenterol 1990; Vol 4 No 9: 637-638.
- **Seewald S., Tiing Leong A., Imazu H., Naga M., Omar S., Groth S., Seitz U., Yan Zhong, Thonke F., Soehendra N. (2008)**
A standardized injection technique and regimen ensures success and safety of N-butyl-2-cyanoacrylate injection for the

treatment of gastric fundal varices.
Gastrointestinal endoscopy, 2008, Vol. 68,
no 3: 447-454.

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- **Maluf-Filho F., Sakai P., Ishioka S., Matuguma S. E. (2001)**
Endoscopic sclerosis versus cyanoacrylate endoscopic injection for the first episode of variceal bleeding: a prospective, controlled, and randomized study in Child-Pugh class C patients.
Endoscopy. 2001; Vol 33 No 5: 421-427.
- **Birmoeller K. F., Soehendra N. (1995)**
Superglue: the answer to variceal bleeding and fundal varices?
Endoscopy. 1995; Vol. 27: 392-396.
- **Gorb J.-P. (1990)**
Endoscopic obturation of esophageal and gastric varices with a cyanoacrylic tissue adhesive.
Can J Gastroenterol 1990; Vol 4 No 9: 637-638.
- **Seewald S., Tiing Leong A., Imazu H., Naga M., Omar S., Groth S., Seitz U., Yan Zhong, Thonke F., Soehendra N. (2008)**
A standardized injection technique and regimen ensures success and safety of N-butyl-2-cyanoacrylate injection for the treatment of gastric fundal varices.
Gastrointestinal endoscopy, 2008, Vol. 68, no 3, 447-454.

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B. Braun Surgical, S.A.