en	Instructions for use/Technical description Hip endoprosthesis cup
USA	Note for U.S. users This Instructions for Use is NOT intended for United States users. Please discard. The Instructions for Use for United States users can be obtained by visiting our website at www.aesculapimplantsystemsifus.com . If you wish to obtain a paper copy of the Instructions for Use, you may request one by contacting your local Aesculap representative or Aesculap's customer service at 1-866-229-3002. A paper copy will be provided to you upon request at no additional cost.
de	Gebrauchsanweisung/Technische Beschreibung Hüftendoprothesenpfanne
fr	Mode d'emploi/Description technique Cotyle pour endoprothèse de la hanche
es	Instrucciones de manejo/Descripción técnica Acetábulo de endoprótesis de cadera
it	Istruzioni per l'uso/Descrizione tecnica Cotile per endoprotesi coxo-femorale
pt	Instruções de utilização/Descrição técnica Acetábulo de endoprótese coxofemoral
nl	Gebruiksaanwijzing/Technische beschrijving Heupprothesekom
da	Brugsanvisning/Teknisk beskrivelse Hofteendoproteseskål
sv	Bruksanvisning/Teknisk beskrivning Höftendoprotesledpanna
fi	Käyttöohje/Tekninen kuvaus Lonkan endoproteesin kuppi
lv	Lietošanas instrukcijas/tehniskais apraksts Gūžas endoprotēzes acetabulārais komponents
lt	Naudojimo instrukcija/techninis aprašas Klubo sąnario endoprotezo gūžduobė
ru	Инструкция по применению/Техническое описание Чашка эндопротеза тазобедренного сустава
cs	Návod k použití/Technický popis Jamka kyčelní endoprotézy
pl	Instrukcja użytkowania/Opis techniczny Panewka protezy stawu biodrowego
sk	Návod na použitie/Technický opis Kĺbová jamka bedrovej endoprotézy
hu	Használati útmutató/Műszaki leírás Csípő endoprotézis vápa
sl	Navodila za uporabo/Tehnični opis Ponev kolčne proteze
hr	Upute za uporabu/Tehnički opis Čašica endoproteze kuka
ro	Manual de utilizare/Descriere tehnică Cupă de endoproteză de șold
bg	Упътване за употреба/Техническо описание Чашка за тазобедрена еднопротеза
tr	Kullanım Kılavuzu/Teknik açıklama Kalça endoprotezi asetabulumu
el	Οδηγίες χρήσης/Τεχνική περιγραφή Κοτύλη ενδοπρόθεσης ισχίου
zh	使用说明 / 技术说明 Plasmafit® 髋臼杯

Aesculap AG | Am Aesculap-Platz | 78532 Tuttlingen | Germany
Phone +49 (0) 7461 95-0 | Fax +49 (0) 7461 95-26 00 | www.bbraun.com

AESCULAP® – a B. Braun brand

TA013800 2021-05 Change No. AE0060729

AESCULAP® Plasmafit®

Hip endoprosthesis cup

1. About this document

Note

General risk factors associated with surgical procedures are not described in these instructions for use.

1.1 Scope

These instructions for use apply to Plasmafit® hip endoprosthesis cups.

Note

For the instructions for use and further information on B. Braun / AESCULAP products please refer to the B. Braun eIFU website at eifu.bbraun.com.

1.2 Safety messages

Safety messages make clear the dangers to patient, user and/or product that could arise during the use of the product. Safety messages are labeled as follows:

⚠ WARNING

Indicates a possible threat of danger. If not avoided, minor or moderate injury may result.

2. Clinical use

Note

The short summary of safety and clinical performance of the product is available in the European Database for Medical Products (EUDAMED).

2.1 Product description

Additional information on Aesculap implant systems is available from B. Braun/Aesculap or from your local B. Braun/Aesculap agency.

2.1.1 Materials

The materials used in the implant are listed on the packaging:

- Outer cup:
 - ISOTAN®_F titanium forged alloy Ti6Al4V in accordance with ISO 5832-3
 - PLASMAPORE® surface coating pure titanium in accordance with ISO 5832-2 or
 - PLASMAPORE®_μ-CaP surface coating pure titanium according to ISO 5832-2 with additional calcium phosphate coating
- Central sealing cap:
 - ISOTAN®_F titanium forged alloy Ti6Al4V according to ISO 5832-3
- Modular polyethylene insert:
 - ultra-high-molecular low-pressure polyethylene according to ISO 5834-2
 - VITELENE® UHMWPE-XE (highly crosslinked vitamin-E stabilized)
 - Pure titanium x-ray markers in accordance with ISO 5832-2 for insert with shoulder
 - X-ray markers made of pure titanium in accordance with ISO 5832-2 or tantalum in accordance with ISO 13782 for asymmetric insert (size dependent)
- Modular ceramic insert:
 - BIOLOX®_{delta} ceramic Al₂O₃/ZrO₂ in accordance with ISO 6474-2

ISOTAN®, PLASMAPORE®, Plasmafit® and Vitelene® are registered trademarks of Aesculap AG, 78532 Tuttlingen / Germany.

BIOLOX® is a registered trademark of CeramTec GmbH, 73207 Plochingen / Germany.

2.2 Areas of use and limitations of use

2.2.1 Intended use

The implant is used:

- As a component of a human hip endoprosthesis: Hip endoprosthesis cup, consisting of outer cup Plasmafit® Poly or Plasmafit® Plus, possibly central screw plug, possibly anchoring screws and modular Plasmafit® inserts (standard, asymmetrical or with shoulder)
- In combination with Aesculap hip endoprosthesis components
- In combination with implant components explicitly approved by Aesculap
- For implantation without bone cement

Note

The options of patient-specific care depend on the available implant components. Implant dimensions and any possible combinations in individual cases can be found in the operating technique instructions for the individual systems.

2.2.2 Indications

Use in cases of acute diseases of the joint that cannot be treated using other therapies:

- Primary and secondary arthrosis
- Rheumatic arthritis
- Joint fractures
- Femoral head necrosis

Note

The manufacturer is not responsible for any use of the product against the specified indications and/or the described applications.

2.2.3 Absolute contraindications

Do not use:

- In patients for whom reconstructive surgery to treat the joint condition is an option
- In cases of acute or chronic infections near the joint
- In cases of secondary diseases influencing function of the joint implant
- In cases of severe osteoporosis or osteomalacia with cementless implants
- In cases of poor bone quality and osseous malformations or diseases in the area of the implant fixation, which may primarily or subsequently affect the stability of the anchorage of the joint replacement
- In case of known hypersensitivity to the implant materials
- In all of the areas of application not listed in the item 'indications'

2.2.4 Relative contraindications

The states listed below, may, individually or combined, result in delayed healing or compromise the success of the operation and have a negative influence on the lifetime of the implants:

- Medical or surgical conditions (e.g. comorbidities, systemic infections) that could prevent the operation from being successful.
 - Anticipated overload on the joint implant
 - Pronounced osseous deformities
 - (prescription) drug abuse or alcoholism
 - Anticipated lack of patient cooperation
- In the presence of relative contraindications, the user decides individually regarding the use of the product.

2.3 Risks, adverse effects and interactions

Within the scope of the legal obligation to provide information, reference is made to the typical risks, interactions and side effects listed below.

Possible risks, side effects and interactions of the application currently known to the manufacturer are:

- Dislocation, loosening, wear, corrosion, disconnection and breakage of the implant components
- Joint dislocation and postoperative changes in leg length
- Primary and secondary infections
- Periprosthetic fractures
- Thromboses, embolisms
- Tissue reaction to implant materials
- Injury to the surrounding tissue, including nerve and vascular damage
- Hematomas and wound healing disorders
- Periarticular calcification
- Reduced joint mobility and flexibility
- Limited load-bearing ability of the joint and joint pain

2.4 Safety information

2.4.1 Clinical user

General safety information

To prevent damage caused by improper setup or operation, and to not compromise the manufacturer warranty and liability:

- ▶ Use the product only according to these instructions for use.
- ▶ Always follow the safety advice and information given in the instructions for use.
- ▶ Ensure that the product and its accessories are operated and used only by persons with the requisite training, knowledge and experience.
- ▶ Store any new or unused products in a dry, clean, and safe place.
- ▶ Keep instructions for use accessible to the user.

Note

The user is obligated to report all severe events in connection with the product to the manufacturer and the responsible authorities of the state in which the user is located.

Notes on surgical procedures

It is the user's responsibility to ensure that the surgical procedure is performed correctly.

Appropriate clinical training as well as a theoretical and practical proficiency of all the required operating techniques, including the use of this product, are prerequisites for the successful use of this product.

Aesculap is not responsible for complications caused by:

- Incorrect indication or implant selection
- Incorrect surgical technique
- Incorrect combination of implant components
- Combination with components of other manufacturers not approved by Aesculap
- Exceeding the limitations of the treatment method or non-observance of essential medical precautions

The user is required to obtain information from the manufacturer if there is an unclear preoperative situation regarding the use of the product.

2.4.2 Product

Product-specific safety information

The operation technique instructions describe the surgical technique and the possible combinations of implant components in individual cases.

Furthermore, all necessary implant components have to be available and all the instruments including special Aesculap implant system instruments have to be present and in working order.

- ▶ Do not use instruments belonging to another Aesculap hip system or made by another manufacturer.

The surgeon shall be responsible for compiling the implant components and their implantation with or without bone cement.

The implant components used, including their article numbers, the name of the implant, the batch number and the serial number, if applicable, must all be documented in the patient records.

- ▶ Observe the instructions for use of the individual implant components.
- ▶ Note further restrictions on the combined implant components in the relevant instructions for use.
- ▶ Only combine modular implant components as specified for the intended use.
- ▶ Observe material, friction coupling diameter and taper specifications.
- ▶ Avoid damage to the implant, in particular in the neck or taper area, caused by using instruments (e.g. surgical HF high frequency devices) close to the implant.
- ▶ Do not under any circumstances damage the implant surfaces.
- ▶ Do not use damaged or surgically removed implant components.
- ▶ Postoperatively, individual patient information, in addition to mobility and muscle training, is of particular importance.

Revising an endoprosthetic implant is a difficult procedure with individual prerequisites that have to be assessed accordingly by the surgeon.

Loosening components, bone or implant fractures and other serious complications cannot be ruled out if the force-transmitting bone structures are damaged.

- ▶ To identify such sources of error as early as possible: Periodically check the condition of the artificial joint by taking appropriate measures.

⚠ WARNING

The implant has not been evaluated for safety and compatibility in the MR environment. It has not been tested for heating, migration or image artifact in the MR environment. Scanning a patient who has this implant may result in patient injury.

Safety instructions for ceramic implant components

If a ceramic implant breaks in vivo, all implant components, including the cup articulation surface, the prosthesis taper (in the event of a failure of the ceramic head) and the modular anchoring surfaces of the cup implant (in the event of a failure of the ceramic insert), will be damaged if loads are continued to be applied to the joint.

If a fracture event is noticed or suspected, the patient must be immediately examined and treated at a specialist endoprosthetic clinic.

- ▶ Inform the patient preoperatively that the fracture of a ceramic implant component is possible in rare cases and that the sudden fracture of a ceramic implant is usually noticeable.
- ▶ Inform the patient preoperatively that if a fracture of a ceramic implant component is suspected, loads may no longer be applied to the affected hip joint to prevent the endoprosthetic components from being destroyed.
- ▶ Inform the patient preoperatively that a noticeable grinding or audible noise of the hip endoprosthesis may be an indication of surface wear of the joint, which must be examined immediately by a specialist.

2.4.3 Sterility

The product is EO sterilized or sterilized using irradiation and sterile packaged.

- ▶ Store implant components in their original packaging. Remove them from their original protective packaging only just prior to implantation.
- ▶ Do not use products from open or damaged sterile packaging.
- ▶ Do not use the product after its use-by date.

The reprocessing of the product affects its functionality. Risk of injury, illness or death due to soiling and/or impaired functionality of the product.

- ▶ Do not reprocess the product.

2.4.4 Storage and transport

- ▶ Store sterile packaged implant components protected from dust in a dry, dark and temperature-controlled room and do not expose them to extreme temperatures or ionizing radiation.
- ▶ Transport and handle implants with care to prevent the packaging and the products from damage.

2.5 Patient education

Within the framework of the patient education, the relevant circumstances needed for consent must be explained to the patient in accordance with their level of understanding, pre-existing knowledge and need for information. This includes:

- Diagnosis, procedure and risk clarification
- Operative procedure
- Advantages and disadvantages of the procedure
- All alternative procedures that can be considered

The patient must be properly informed about the procedure and in particular about the following information:

- The functionality of the prosthetic joint is essentially inferior to that of the natural one.
- The prosthetic joint can bring about only limited improvement in the patient's condition in relation to the pre-operative condition.
- The prosthetic joint can come loose due to excessive load, wear and tear, or infection.
- The lifetime of the prosthetic joint depends on the patient's body weight and the load to which the joint is subjected.
- The prosthetic joint must not be overstrained by extreme loads, hard physical labor or sporting activities.
- Revision surgery may become necessary if the implant loosens.
- In the event that revision surgery is performed, it may not be possible under certain circumstances to restore joint mobility and flexibility.
- The revision of a hip endoprosthesis stem is a complex joint replacement operation.
- Generally, the revision joint replacement will be inferior in its performance to the primary joint replacement.
- The patient must undergo medical follow-up examinations of the prosthetic joint at regular intervals.
- Damage to the weight-bearing bone cement and/or bone structures can cause loosening of the components, fracture of the bone or implant and other serious complications. To ensure the earliest possible detection of such catalysts of implant dysfunction, the prosthetic joint must be checked periodically, using appropriate techniques.

2.6 Application

2.6.1 Documentation

The user shall devise an operation plan that specifies and accurately documents the following:

- Selection of the implant components and their dimensions
- Positioning of the implant components in the bone
- Location of intraoperative landmarks

Important information concerning the implanted product and the operation can be noted down on the patient ID. The patient ID can be ordered separately from the manufacturer.

Each package contains additional labels showing the designation, article and lot number and – if applicable – the individual serial number of the product.

- Use these labels for documentation in the patient's file (for the hospital) and the patient ID (for the patient).

The surgical procedure has been explained to the patient, and the patient's consent has been documented.

2.6.2 Implantation

Note
For detailed information of implantation steps see implant system specific operating technique instructions.

The preparation of the implant bed and the implantation require the steps listed below:

- use spherical acetabulum reamers to prepare the acetabulum to the required implant size.
The nominal diameter of the Plasmafit® outer cup shall match the size of the acetabular reamer last used.
- Do not remove any bone structures that are required to anchor the implant.
- The implant bed has to be suitable for cement-free implantation.
- If necessary, use trial cups to check the implant bed.
- Observe the cup orientation:
The arrow on the Plasmafit® exterior cup has to point caudally to the acetabular incision.
- If the primary stability is insufficient, use anchoring screws or switch to the next larger implant or use a cemented cup.

- ⚠ **WARNING**
Bone fractures in the implant bed will affect the anchoring of the implants!
 - **Avoid bone fractures by using cautious surgical technique.**
 - **Treat bone fractures with appropriate intra- and postoperative measures.**
 - **Observe the correct application of the implant components.**

- ⚠ **WARNING**
Loose anchoring screws or a loose central sealing cap can damage implant components!
 - **If anchoring screws or a central sealing cap are used, they must be firmly tightened.**

- Make sure that the implant position permits the full joint mobility and that dislocations are precluded.
Joint dislocation can also occur with modular inserts (asymmetric type or type with rim) if the cup is incorrectly positioned or if there is muscular instability. When inserts with rims are used, the range of motion is limited.
- Make sure the insert matches the nominal diameter of the prosthesis head.
- Ensure that the insert matches the cup in relation to the mating dimension, labeled by the letters A–M.

- ⚠ **WARNING**
The implant components may break/fail!
 - **Only use polyethylene inserts for Plasmafit® Poly cups.**
 - Thoroughly clean and dry before inserting the modular insert:
 - outer taper of the insert,
 - inner taper of the exterior cup.
 - Carefully remove all foreign bodies, such as tissue parts, bone or cement particles, especially from the taper clamping area.
 - Perform a trial reduction and check joint mobility, joint stability and leg length.
 - If a BIOLOX® ceramic insert is used:
 - Observe the instructions for use for BIOLOX® ceramic inserts (TA009464).
 - Prior to inserting the insert, make sure the taper in the exterior cup is not damaged in any way.
 - Use impactor with plastic head for insertion
 - Only use BIOLOX® insert with BIOLOX® hip endoprosthesis heads having a suitable diameter.
 - Only use BIOLOX® ball heads. (Do not use balls heads e.g. metal, zirconium oxide or aluminum oxide of other manufacturers).
 - Insert the cup inserts into the exterior cup only once.
 - If the cup insert has to be removed, replace it with a new implant.
 - Before closing the wound, check the positioning and orientation of the implant components, if necessary by means of fluoroscopy.
 - Remove all exposed foreign bodies such as bone cement and bone residues prior to closing the wound, because bone cement particles and bone fragments that get caught between the articulation surfaces can cause abnormal wear of the artificial joint.
- If the ceramic prosthesis head is broken:
 - Remove all ceramic fragments carefully and completely.
 - If possible, clean the site using pulsed irrigation/suction systems (jet lavage)
 - The polyethylene or ceramic bearing partner must be replaced, even if the cup is firmly seated.

Note
Observe the relevant instructions for use when replacing the components.

- ⚠ **WARNING**
Ceramic splinters can injure the patient/surgical team!
 - **Take care to remove any splinters using suction or other means.**
 - **Do not perform syringe irrigation.**

- ⚠ **WARNING**
Risk of wear and tear due to ceramic fragments inside the joint following revision surgery that involved the failure of ceramic components!
 - **Do not implant a metal head (risk of third-body wear).**
 - **Replace the prosthesis stem or treat the joint using a modular BIOLOX® Option system hip endoprosthesis head.**
 - **Replace the cup insert or the PE cup even if the cup is firmly seated.**
 - **The ceramic-ceramic bearing option is preferable to that of ceramic-PE.**

3. Disposal

- ⚠ **WARNING**
Risk of infection due to contaminated products!
 - **Adhere to national regulations when disposing of or recycling the product, its components and its packaging.**

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