
# AESCULAP®

(DE) **Gebrauchsanweisung**  
PG150  
Flow 50

(ES) **Instrucciones de manejo**  
PG150  
Flow 50

(EN) **Instruction for use**  
PG150  
Flow 50

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(FR) **Mode d'emploi**  
PG150  
Flow 50

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Copyright © AESCULAP AG

**EN**

These instructions for use contain information that is subject to copyright. All rights reserved. These instructions for use may not be photocopied, duplicated on microfilm, or otherwise copied or distributed, completely or in part, without the approval of AESCULAP AG.

We reserve the right to technical changes without prior notification due to the continuous further development of our products. Function or design may partially differ from the description in the instructions for use. Please contact us for additional information about this or any of our other products.

Some of the parts and equipment referred to in these instructions for use are associated with registered trademarks but are not identified as such. It should therefore not be assumed that the absence of the trademark symbol indicates that any given designation is not subject to trademark protection.

Users of this product should not hesitate to point out to us any errors or issues concerning these instructions for use.

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## Table of contents

<b>1</b>	<b>Important Usage Notes.....</b>	<b>3</b>
<b>2</b>	<b>Safety Information.....</b>	<b>4</b>
<b>3</b>	<b>General Information .....</b>	<b>5</b>
3.1	Device Description .....	5
3.2	Intended Use and Contraindications .....	5
3.2.1	Intended Use.....	5
3.2.2	Contraindications .....	5
3.3	Warnings .....	6
3.3.1	General CO <sub>2</sub> Insufflation Warnings .....	6
3.3.2	Warnings Standard and Bariatric Indications .....	9
3.3.3	Warnings Pediatric Indication .....	11
3.3.4	Warnings Vessel Harvesting Mode (Cardio Indication).....	13
3.3.5	Warnings Fieldflooding Mode (Cardio Indication).....	14
3.4	General Precautions .....	14
<b>4</b>	<b>Initial Device Set-up.....</b>	<b>16</b>
4.1	Scope of Delivery .....	16
4.2	Device Setup and Connection.....	17
4.3	Front of Device .....	18
4.4	Rear of the Device .....	19
4.5	Gas Connection.....	19
4.5.1	Connecting a Gas Bottle.....	20
4.5.2	Connecting to Central Gas Supply.....	21
4.6	Gas Supply.....	22
4.7	Switching Device On.....	22
4.8	Switching Device Off .....	23
4.9	Demo Mode.....	23
<b>5</b>	<b>Operating the Device - General.....</b>	<b>24</b>
5.1	Monitor with Touch Screen .....	24
5.2	Insufflation Tube Sets .....	24
5.3	Connecting an Insufflation Tube Set .....	25
5.4	Gas Heating.....	25
5.5	Using the Smoke Evacuation Function .....	27
5.5.1	Replacing the Smoke Evacuation Filter Cassette .....	27
5.5.2	Replacing the Disposable Smoke Evacuation Tube Set .....	28
5.6	Displaying>Selecting Insufflation Indication.....	28
5.7	Setting the Nominal Pressure - All Indications Except Fieldflooding .....	29
5.8	Setting the Nominal Gas Flow Rate - All Indications .....	29
5.8.1	Veress Insufflation .....	30
5.8.2	Highflow Insufflation .....	30
5.8.3	Profiles.....	30
5.9	Function Test before Device Use .....	31
5.10	Gas Consumption Display.....	32
5.11	Starting/Stopping Insufflation .....	32
<b>6</b>	<b>Indications of the Flow 50 Insufflator .....</b>	<b>33</b>
6.1	Standard Indication .....	33
6.2	Bariatric Indication .....	35
6.3	Pediatric Indication .....	38
6.4	Cardio Indication .....	40
6.4.1	Fieldflooding .....	40
6.4.2	Vessel Harvesting .....	40
<b>7</b>	<b>Menu .....</b>	<b>43</b>
7.1	Profiles.....	44
7.1.1	Creating Profiles .....	44
7.1.2	Change Profile .....	44
7.1.3	Delete Profile .....	44
7.2	Indication Setting .....	45
7.2.1	Insufflation .....	45
7.2.2	Smoke Evacuation .....	45
7.3	System .....	45
7.4	Service.....	46
7.5	Protocol .....	46
7.6	User Menu (Overview) .....	46
<b>8</b>	<b>Safety Functions.....</b>	<b>50</b>
8.1	Automatic Venting System .....	50
8.2	Occlusion .....	50
8.3	Contamination .....	51
8.4	Leakage Detection .....	51
8.5	Gas Heating - Overtemperature .....	51

---

<b>9</b>	<b>Care and Maintenance.....</b>	<b>52</b>
9.1	Cleaning the Device .....	52
9.2	Maintenance Intervals.....	52
9.3	Maintenance Carried out by Authorized Service Technician .....	52
9.4	Replacing the Fuse .....	53
9.5	Care and Maintenance of Reusable Tube Sets (NOT APPROVED FOR SALE IN THE UNITED STATES) .....	54
9.5.1	Cleaning and Disinfecting.....	55
9.5.2	Sterilization of Reusable Tube Sets (NOT APPROVED FOR SALE IN THE U.S.).....	57
<b>10</b>	<b>Annual Inspection .....</b>	<b>59</b>
10.1	Electrical Safety Test.....	59
10.2	Basic Function Test (in Standard Indication) .....	59
10.3	Pressure Sensor Test (in Standard Indication and Veress Mode) .....	60
10.4	Pressure Monitoring Test (in Standard Indication and Veress Mode).....	60
10.5	Venting Valve Test .....	61
10.6	Device Maximum Pressure Test .....	61
10.7	Gas Flow Rate Test.....	62
<b>11</b>	<b>Electromagnetic Compatibility.....</b>	<b>63</b>
11.1	Impact of Mobile and Portable HF Communication Devices.....	63
11.2	Electrical Connections.....	63
11.3	Guidance and Manufacturer's Declaration - Electromagnetic Immunity .....	64
11.4	Guidance and Manufacturer's Declaration - Electromagnetic Emissions .....	65
11.5	Recommended Separation Distances between Portable and Mobile RF Communications Equipment and the Flow 50.....	66
<b>12</b>	<b>Error and Warning Messages.....</b>	<b>67</b>
<b>13</b>	<b>Technical Data .....</b>	<b>69</b>
<b>14</b>	<b>Accessories .....</b>	<b>70</b>
<b>15</b>	<b>Glossary.....</b>	<b>71</b>
<b>16</b>	<b>Appendix .....</b>	<b>72</b>
16.1	Test Log .....	72
	<b>Index .....</b>	<b>73</b>

## 1 Important Usage Notes

Read the instructions for use carefully and become familiar with the operation and function of the device and the accessories before use during surgical procedures. Non-observance of the instructions listed in this manual can lead

- to life-threatening injuries of the patient,
- to severe injuries of the surgical team, nursing staff or service personnel, or
- to damage or malfunction of device and/or accessories.

The manufacturer reserves the right to modify the appearance, graphics, and technical data of the product through continued development of its products.

**Subject to technical changes**

Paragraphs marked with the words WARNING, CAUTION, and NOTE carry special meanings. Read these paragraphs with special attention.

**Please note**

### **WARNING!**

**The safety and/or health of the patient, user, or a third party are at risk. Comply with this warning to avoid injury to the patient, user, or third parties.**

### **CAUTION!**

**These paragraphs include information provided to the operator concerning the intended and proper use of the device or accessories.**

### **NOTE!**

**These paragraphs contain information to clarify the instructions or provide additional useful information.**

**Federal Law (only for U.S. market)**

## 2 Safety Information

**CAUTION: Federal Law restricts this device to sale by or on the order of a physician.**

**Exclusion of liability**

The manufacturer is not liable for direct or consequential damages, and the warranty becomes null and void if:

- the device and/or the accessories are improperly used, prepared, or maintained;
- the instructions and rules in the instructions for use are not adhered to;
- unauthorized persons perform repairs, adjustments, or alterations on the device or accessories;
- unauthorized persons open the device;
- the prescribed inspection and maintenance schedule is not adhered to.

The handing over of technical documents does not constitute authorization to make repairs or alterations to the device or accessories.

**WARNING: Modifying the device Flow 50 is not permitted.**

**Authorized service technician**

Only an authorized service technician may perform repairs, adjustments, or alterations on the device or accessories and use the service menu. Any violation will void the manufacturer's warranty. Authorized service technicians are only trained and certified by the manufacturer.

**Care and maintenance**

The service and maintenance of the device and its accessories has to be carried out as per instructions to ensure the safe operation of the device. For the protection of the patient and the operating team, check that the device is complete and functional before each use.

Maintenance of the device may not be performed during the operation.

**Contamination**

Before shipping, decontaminate device and accessories in order to protect the service personnel. Follow the instructions listed in this manual. If this is not possible,

- the product must be clearly marked with a contamination warning and
- is to be double-sealed in safety foil.

The manufacturer has the right to reject contaminated products for repair.

**Disposal**

The product must be prepared by the operator for disposal. When disposing or recycling the product or its components, make sure to comply with national rules and regulations at all times! A product labeled with this symbol must be disposed separately according to local rules about the waste collection of electrical and electronic devices. The disposal is carried out free by the manufacturer if within the European Union. Please contact your local B. Braun/Aesculap representative or distributor for additional information about the disposal of your product; see chapter 9 Care and Maintenance.

### 3 General Information

#### 3.1 Device Description

The device Flow 50 is an insufflator and serves to create a cavity through the inflation of CO<sub>2</sub> with a preset flow rate and pressure. The gas will be transported to the patient via a tube connected with a trocar. The insufflator limits the pressure to maximum 30 mmHg (mm mercury column) and the flow rate to maximum 50 l/min.

**The device measures the actual pressure in the cavity and compares it with the set nominal pressure. The function of this device is to maintain the nominal pressure. Any overpressure within the cavity is lowered to the preset nominal pressure by the automatic venting system.**

The user and patient comfort can be increased by additional accessories:

1. Tube set with integrated gas heating to heat the insufflation gas.
2. Tube set with evacuation function to remove surgical smoke from the cavity.

These functions are available when using the appropriate tubing sets; see chapter 5.2 Insufflation Tube Sets.

The insufflator is a pneumatic device that works through the port on the back, with central gas as well as bottled gas. It is operated via buttons and indicators or displays on the front of the device.

#### 3.2 Intended Use and Contraindications

##### 3.2.1 Intended Use

The device Flow 50 is a CO<sub>2</sub> insufflator intended for use during diagnostic and/or therapeutic endoscopic procedures to distend a cavity by filling it with CO<sub>2</sub> gas. The indications **Standard, Bariatric and Pediatric** of the device are indicated to fill and distend a peritoneal cavity with gas during a laparoscopic procedure. The **Pediatric** indication is indicated for pediatric laparoscopic procedures. **Vessel Harvesting** is indicated for use during endoscopic vessel harvesting procedures to create a cavity along the saphenous vein or radial artery. The **Fieldflooding** is indicated for use in open and endoscopically assisted cardiac surgery to displace ambient air from the surgical site.

##### 3.2.2 Contraindications

The device should not be used to fill a cavity with CO<sub>2</sub> if an endoscopy is contraindicated. Please consult the manual of your endoscope for absolute and relative contraindications. The device is not suitable for hysteroscopic insufflation procedures, i.e., it may not be used to distend the uterus.

### 3.3 Warnings

#### 3.3.1 General CO<sub>2</sub> Insufflation Warnings

**WARNING!**

**Condensation / Water penetration**

**Protect device from moisture. Do not use if moisture has penetrated the device.**

**WARNING!**

**Check all factory settings**

**Factory settings are not mandatory settings for the physician. The physician is responsible for all settings affecting the surgical procedure.**

**WARNING!**

**Technique and procedures**

**Only the physician can evaluate the clinical factors involved with each patient and determine if the use of this device is indicated. The physician must determine the specific technique and procedure that will accomplish the desired clinical effect.**

**WARNING!**

**Not explosion-proof**

**The device is not explosion-proof. Do not operate the device in the vicinity of explosive anesthetic gases and not in the vicinity of oxygen-enriched environments.**

**WARNING!**

**Risk of electrical shock**

**To prevent electrical shock, do not open this device. Never open this device yourself. Refer servicing to qualified service personnel.**

**WARNING!**

**Risk of electrical shock**

**To avoid the risk of electrical shock, only use this device when connected to a properly grounded power supply network.**

**WARNING!**

**Replace the fuse**

**Replace the fuse only with a fuse provided by the manufacturer (see chapter 14 Accessories).**

**WARNING!**

**Professional qualification**

**This manual does not include descriptions or instructions for surgical procedures/techniques. It is also not suitable for training physicians in the use of surgical techniques. Medical accessories and devices may be used only by physicians and medical assistants under the direction of a physician with the appropriate technical qualification.**

**WARNING!****Automatic device self-test**

The device self-test must be performed prior to each surgery. Because the device self-test is performed during initial start up, the unit must be power cycled (off/on) prior to each surgery.

**WARNING!****Sterile mediums and accessories**

Always work exclusively with sterile substances and mediums, sterile fluids, and sterile accessories if so indicated.

**WARNING!****Replacement device and accessories**

In case the device or any of the accessories fail during surgery, a replacement device and replacement accessories should be kept within close proximity to be able to finish the operation with the replacement components.

**WARNING!****Cleaning the device**

**Do not sterilize the device.**

**WARNING!****Device defect**

If a device defect is suspected or confirmed, do not use it. Make sure the device can no longer be used until a qualified service technician conducts the appropriate tests and repairs.

**WARNING!****Positioning the patient**

Always position the patient lower than the device to prevent body fluids from leaking into the insufflation tube. Actual pressure may increase and fluid may penetrate the insufflation tube if the patient is repositioned during surgery. If this occurs, immediately disconnect the insufflation tube. When the patient is repositioned onto his or her side, internal tissue may block the insufflation channel. Always insufflate through the elevated side of the patient.

**WARNING!****Removing the insufflation tube**

Always disconnect the insufflation tube after ending surgery and before switching off the device to prevent backflow of bodily fluids. Fluid may penetrate the insufflation tube whenever you change the gas bottle and/or when you stop the gas flow during the operation. If this happens, you must immediately disconnect the insufflation tube from the trocar or from the device.

**WARNING!**

**Backflow**

**Body secretions or contaminated gas can flow counter to the insufflation direc-**

**tion and block the insufflation filter, if**

- the actual pressure is higher than the nominal pressure or**
  - the automatic venting valve is activated.**
- 

**WARNING!**

**Gas flow**

**A high gas flow can occur due to large leaks within the surgical system or instru-**

**ment. This can result in a false actual pressure reading, which in turn may endan-**

**ger the patient. In case of a disrupted gas flow, you should therefore inspect**

**device, tube, and instruments immediately. Surgical procedures should be per-**

**formed with a gas flow of 4 to 10 l/min. An even lower gas flow is recommended**

**for diagnostic purposes.**

---

**WARNING!**

**Keep filled CO<sub>2</sub> bottle on hand**

**Always keep a filled CO<sub>2</sub> bottle on hand ready for replacement. This avoids hav-**

**ing to interrupt surgery due to a lack of insufflation gas (see chapter 4.5.1 Con-**

**nnecting a Gas Bottle).**

---

**WARNING!**

**Gas supply**

**Maintain adequate gas supply at all times.**

---

**WARNING!**

**Contamination**

**Do not use device and/or accessories if signs of contamination are detected.**

**Make sure the device or/and accessories can no longer be operated until a qual-**

**ified service technician conducts the appropriate tests and repairs.**

---

**WARNING!**

**Fatigue symptoms**

**When there is a high level of CO<sub>2</sub> consumption, you should make sure to supply**

**the operating area with enough fresh air, since an increasing CO<sub>2</sub> level in the air**

**can cause the medical personnel to suffer fatigue symptoms, an inability to con-**

**centrate, unconsciousness, or even death.**

---

**WARNING!**

**The insufflation of CO<sub>2</sub> should be done carefully and while monitoring the pa-**

**tient's response. The user, particularly the anesthetist, should be informed about**

**possible cardiovascular and respiratory problems of the patient and monitor**

**these intra-operatively.**

---

**WARNING!****Medically pure CO<sub>2</sub>**

Make sure to use only medically pure CO<sub>2</sub>. Other gases (e.g., helium, N<sub>2</sub>O, argon), mixtures of gases, gases with entrapped liquids, or polluted gases must not be used with this device.

**WARNING!**

For the safety of the patient please fill the tube set with CO<sub>2</sub> gas prior to beginning the insufflation by activating the insufflation for a few seconds and then turning it off again before introducing the insufflation instrument to the cavity and beginning the surgery.

**WARNING!****Lowest flow and pressure**

Depending on age and health condition of the patient, the smallest possible flow and pressure for establishing the pneumoperitoneum or pneumorectum should be selected. It is recommended not to exceed insufflation pressures of 15 mmHg in transanal minimally invasive surgery procedures.

**WARNING!**

The venting rate of the automatic venting system is limited. Always monitor the actual pressure when using additional insufflation sources.

**WARNING!****Connecting the tube**

Always use the proper tube set for the device. The tube outlet may only be connected to instruments which are intended for intra-abdominal CO<sub>2</sub> insufflation.

**WARNING!****Electronic device control**

Do not close the valve at the trocar sleeve during surgery. The electronic control unit of the device adjusts the actual pressure as desired.

### 3.3.2 Warnings Standard and Bariatric Indications

**WARNING!****Idiosyncratic reactions**

Patients with sickle cell anemia or pulmonary insufficiency may have a higher risk of metabolic imbalance related to excessive CO<sub>2</sub> absorption (idiosyncratic reaction).

**WARNING!****CO<sub>2</sub> absorption**

CO<sub>2</sub> is absorbed during insufflation (intravasation). This means the body absorbs part of the CO<sub>2</sub> gas used for insufflation. CO<sub>2</sub> concentrations in the blood or respiratory system that are too high can lead to death of the patient in extreme cases. To lower this risk, always carefully and closely monitor the patient's vital signs during the entire insufflation process and make sure patient is breathing well. Sufficient respiration can help avoid or limit problems with CO<sub>2</sub>. High pressure or a high gas flow promotes CO<sub>2</sub> absorption. The abdomen is sufficiently distended using a pressure between 10 to 15 mmHg. Pressure values above 15 mmHg are required for only a few cases but do increase the risk of intravasation. Never exceed the maximum intra-abdominal pressure of 30 mmHg.

**WARNING!****Metabolic and cardiac reactions**

Insufflating CO<sub>2</sub> may result in metabolic acidosis. This can lead to cardiac irregularities expressed with the following symptoms:

- Reduced respiration with restricted diaphragm function
- Hypercapnia
- Reduction of venous reflux
- Reduced cardiac output
- Metabolic acidosis

**WARNING!****Hypothermia/monitoring body temperature**

The gas flow can lead to a lowering of the patient's body temperature during insufflation. Hypothermia during insufflation can cause heart and cardiovascular problems. The risk for hypothermia can be significantly reduced with the use of gas that is pre-warmed to body temperature. Therefore, you must always monitor the patient's body temperature during the entire insufflation process. Make especially sure that the following, hypothermia-fostering, surgical conditions are avoided as best as possible:

- High gas flow due to large leaks
- Long surgeries
- Use of irrigation and infusion solutions that are not preheated to body temperature.

**WARNING!****Dehydration**

Insufflation can lead to dehydration of the tissue. This can result in organ tissue damage and cardiovascular reactions of the patient. Long surgeries and large leaks increase the risk of dehydration (especially at the insertion points of the trocars or when changing instruments).

**WARNING!****Embolism/insufflation of internal organs**

Improper placement of the insufflation instrument could cause gas penetrating a vessel or an internal organ, resulting in air or CO<sub>2</sub> embolisms. To reduce the risk, use a low flow rate for the first insufflation and ensure that the insufflation instrument is correctly positioned. Check the position of the insufflation instrument immediately if the actual pressure rapidly reaches the nominal pressure value. CO<sub>2</sub> embolisms can also be caused by a high intra-abdominal pressure. Avoid high-pressure settings and close damaged blood vessels at once.

**WARNING!****Subcutaneous emphysema**

Incorrect placement of a cannula or a trocar into subcutaneous tissue may lead to emphysema. To reduce the risk, use a low gas flow rate for the first insufflation and ensure that the insufflation instrument is correctly positioned. Long surgeries, the use of many access points, duration and size of leaks at these points may also contribute to emphysema. Be sure to close leakages in trocar accesses immediately.

**WARNING!****Subcutaneous emphysema**

When puncturing the thicker abdominal wall of morbidly obese patients with the Veress cannula or the trocar, carefully monitor the correct position of the instrument in the abdomen.

**WARNING!****Additional insufflation sources/automatic venting system**

Make sure the automatic venting system is activated (see chapter 7 Menu). The use of additional insufflation sources increases the intra-abdominal pressure. Continuously monitor intra-abdominal pressure over the course of the entire insufflation if additional sources are used.

### 3.3.3 Warnings Pediatric Indication

**WARNING!****Gas flow limit**

The gas flow may not exceed 14 l/min when performing a laparoscopy on newborns or patients weighing less than 25 kg (approximately 55 US pounds).

**WARNING!****Recommended work settings**

The flow values listed for laparoscopic procedures performed on newborns, infants, and children are only suggested values. The selection of the suitable flow and pressure values is solely the responsibility of the attending physician.

**WARNING!****Pneumolabium/pneumoscrotum**

Children are at risk of a pneumolabium or pneumoscrotum.

**WARNING!**

Only specially trained and qualified personnel may use this device on children.

---

**WARNING!**

**Increased airway pressure**

**When laparoscopic procedures are performed on children, the increased intra-abdominal pressure also increases the risk for higher airway pressures. Always strictly monitor respiration and airway function when performing laparoscopic procedures on children younger than 12 years of age.**

---

---

**WARNING!**

**Compression of the vena cava**

**When insufflating the abdomen of a child with medical CO<sub>2</sub>, an increased risk of compressing the vena cava exists. This risk can be reduced by monitoring the systolic and diastolic blood pressure during the entire surgery.**

---

---

**WARNING!**

**Haemodynamic stability**

**A laparoscopy performed on children younger than 12 years of age may result in problems of the haemodynamic system due to the CO<sub>2</sub> content in the blood. It is recommended to increase the breathing rate of the patient and to work with low flow values and pressure values not exceeding 12 mmHg. The patient's circulatory system should be monitored at all times.**

---

---

**WARNING!**

**Hypothermia**

**The insufflation gas flow usually drops significantly after the target pressure has been reached and it is then only required to maintain the abdominal pressure. However, leaks within the abdomen or the instrument can lead to a constant gas flow of above 1 l/min. When operating on children younger than 12, a gas flow of more than 1 l/min poses an increased risk of hypothermia for the patient. Corresponding measures to prevent hypothermia include the use of blankets or pre-warmed gas. The patient's body temperature has to be monitored at all times during surgery.**

---

---

**WARNING!**

**Laparoscopy with CO<sub>2</sub> should not be performed on children with cardiovascular problems.**

---

---

**WARNING!**

**Because pediatric patients are especially susceptible to hypercapnia, it is recommended to establish an end-tidal CO<sub>2</sub> monitoring routine.**

---

### 3.3.4 Warnings Vessel Harvesting Mode (Cardio Indication)

#### **WARNING!**

##### **CO<sub>2</sub> absorption**

Due to the special surgical procedures - start of the heart bypass operation, and the endoscopic removal of the vessel - special care has to be taken as CO<sub>2</sub> is always absorbed through the tissue of the patient during insufflation (intravasation). This means the body absorbs part of the CO<sub>2</sub> gas used for insufflation. CO<sub>2</sub> concentrations in the blood or respiratory system that are too high can lead to death of the patient in extreme cases. To lower this risk, always carefully and closely monitor the patient's vital signs during the entire insufflation process and make sure patient is breathing well. Sufficient respiration can help avoid or limit problems with CO<sub>2</sub>. High pressure or a high gas flow promotes CO<sub>2</sub> absorption.

#### **WARNING!**

##### **Metabolic and cardiac reactions**

Insufflating CO<sub>2</sub> may result in metabolic acidosis. This can lead to cardiac irregularities expressed with the following symptoms:

- Reduced respiration with restricted diaphragm function
- Hypercapnia
- Reduction of venous reflux
- Reduced cardiac output
- Metabolic acidosis

#### **WARNING!**

##### **Dehydration**

Insufflation can lead to dehydration of the tissue. This can result in organ tissue damage and cardiovascular reactions of the patient. Long surgeries and large leaks increase the risk of dehydration (especially at the insertion points of the trocars or when changing instruments).

#### **WARNING!**

##### **Embolism/insufflation of internal organs**

Improper placement of the insufflation instrument could cause gas penetrating a vessel or an internal organ, resulting in air or CO<sub>2</sub> embolisms. To reduce the risk, use a low flow rate for the first insufflation and ensure that the insufflation instrument is correctly positioned. Check the position of the insufflation instrument immediately if the actual pressure rapidly reaches the nominal pressure value. CO<sub>2</sub> embolisms can also be caused by a high intra-abdominal pressure. Avoid high-pressure settings and close damaged blood vessels at once.

#### **WARNING!**

Only specially trained and qualified personnel may use this device for the endoscopic vessel harvesting procedure.

EN

**WARNING!**

Before using the insufflator to endoscopically harvest vessels, please check whether the instrument used is intended for CO<sub>2</sub> insufflation.

**WARNING!****Pneumoperitoneum**

When a vessel is harvested from the leg of a patient with a perforated groin, it is possible for CO<sub>2</sub> to reach the abdomen and cause a pneumoperitoneum. Make sure the abdomen does not fill with CO<sub>2</sub> during surgery.

**3.3.5 Warnings Fieldflooding Mode (Cardio Indication)****WARNING!**

A gas flow of 10 l/min should not be used longer than up to 30 seconds.

**WARNING!**

When using accessories not specified in the user manual, make sure the CO<sub>2</sub> is not directly applied to soft tissue focusses.

**3.4 General Precautions****CAUTION!****Continuous operation**

After 24 hours of continuous operation, a device self-test must be carried out.

- Switch device off and on again.

**CAUTION!****Endoscope**

The device may only be connected with endoscopes designed for and featuring the technical specification permitting such a combined use. Any utilized endoscopes must comply with the most recent versions of IEC 60601-2-18 and ISO 8600. Combining/connecting with other devices generates a medical electrical system (MES). The system configurator is responsible for compliance with the standard IEC 60601-1 / EN 60601-1 in its latest version.

**CAUTION!****Electrical Interference**

(See Chapter 11 Electromagnetic Compatibility). Care was taken during the development and testing of this device that electrical interference from this device or other devices or instruments was practically eliminated. However, if you still detect or suspect such interference, please follow these suggestions:

- Move this, the other or both devices to a different location
- Increase distance between used devices
- Consult an electro-medical expert

---

**CAUTION!****Service connection**

Access to the service menu is restricted to authorized service personnel. The connected equipment must comply with the standard EN 60950 in the currently valid version. Do not connect a device to the service connection during surgery.

EN

---

**CAUTION!****Peripheral devices**

Additional peripheral equipment connected to interfaces of the medical monitor has to meet the requirements of the following specifications in the respective current valid version: IEC 60601-2-18 / EN 60601-2-18 for endoscopic devices and IEC 60601-1 / EN 60601-1 for electrical medical devices. All configurations have to comply with IEC 60601-1 / EN 60601-1 specifications. Whoever connects additional equipment to signal output or signal input is considered the system configurator and as such is responsible for complying with requirements of the standard IEC 60601-1 / EN 60601-1.

---

**CAUTION!**

Any power cables employed by the user that are not provided by the manufacturer must meet the safety requirements of the national standards in the respective current valid version.

---

**CAUTION!****Use of other accessories, other transducers and cables**

The use of accessories, transducers, and cables other than those specified, with the exception of transducers and cables sold by the MANUFACTURER of the ME EQUIPMENT or ME SYSTEM (see Chapter 15 Glossary) as replacement parts for internal components, may result in increased EMISSIONS or decreased IMMUNITY of the ME EQUIPMENT or ME SYSTEM.

## 4 Initial Device Set-up

The operation of the device Flow 50 is reserved for persons with the relevant professional qualifications and trained to use the device.

**NOTE!**

**Locations**

**The device is only to be used in a professional facility healthcare environment.**

### 4.1 Scope of Delivery

- PG150 Insufflator Flow 50
- TA014996 Instructions for use

**Delivery inspection**

Always check all parts and optional accessories of the device immediately after receiving the shipment. The manufacturer considers only replacement claims that have been immediately submitted or reported to a sales representative or an authorized service company.

**Returning the device**

Use the original packaging material if a device has to be returned. The manufacturer is not responsible for transport damages due to insufficient or unsuitable shipping packaging.

Please make sure that all required information has been supplied:

- Name of owner
- Address of owner
- Device type and model:
- Serial number of the device (see type plate)
- Detailed description of defect

**Technical service address:**

Aesculap Technical Service  
Am Aesculap-Platz  
D-78532 Tuttlingen / Germany  
Phone: +49 7461 95-1601  
Fax: +49 7461 14-939  
E-Mail: [ats@aesclap.de](mailto:ats@aesclap.de)

Please contact the address above for additional service addresses.

Or in the US:

Aesculap Inc.  
Attn. Aesculap Technical Services  
615 Lambert Pointe Drive  
Hazelwood  
MO, 63042 USA  
Aesculap Repair Hotline  
Phone:+1 (800) 214-3392  
Fax:+1 (314) 895-4420  
Other service addresses can be obtained from the address indicated above.

#### 4.2 Device Setup and Connection

Place the device on a flat surface free of vibration located in a dry environment. The ambient temperature and humidity must meet the requirements mentioned in Chapter 13 Technical Data.

#### Setting up

##### **WARNING!**

##### **ME System (Medical Electrical System)**

**Use only parts and/or devices from ME systems (see Chapter 15 Glossary) in patient environments in compliance with the standard IEC60601-1 in the respective currently valid version.**

##### **CAUTION!**

##### **ME Device in Rack**

**If it should be necessary to operate the device close to stacked with other devices, the ME device or ME system (see chapter 15 Glossary) should be monitored to ensure it works properly as configured.**

##### **CAUTION!**

##### **ME systems**

**The medical electrical (ME) device is suitable for integration in ME equipment systems (see Chapter 15 Glossary). Operation of the ME device in vicinity of non-ME devices may result in voiding the intended use of the ME device.**

##### **CAUTION!**

**Equipment should be positioned such that power cord can be easily disconnected.**

##### **CAUTION!**

##### **Ventilation of the device**

**Avoid device overheating. Ensure free air circulation especially to the bottom and rear of the device (rear panel distance of at least 10 cm).**

##### **CAUTION!**

**Position the device in such a way that it is easy to operate and switch off.**

##### **CAUTION!**

**Place the device outside the sterile field.**

##### **CAUTION!**

##### **User position**

**To avoid a malfunction, the user must be positioned correctly**

- **within a display viewing angle of  $\pm 50^\circ$  to operate the device**
- **up to 2 m from the device front for monitoring the actual values**

EN

**Mains connection**

Ensure the connection data and power supply comply with DIN VDE or national requirements. The mains connection cable may be plugged only into a properly installed, grounded safety wall socket (shockproof socket) (see DIN VDE 0100-710). Read the device label located in rear of device (type plate) to determine the operating voltage of the device.

**CAUTION!**

**Check to make sure the available mains voltage matches the data listed on the label attached to the back of the device. Incorrect voltage can cause errors and malfunctions and may destroy the device.**

**Grounding contact**

The power connection must be equipped with a grounding contact. Use the original power cable (if included in scope of delivery) to establish a connection between the mains wall socket and the non-heating device plug located in the rear of the device.

**Only for U.S. operators**


Only use a certified (UL-listed), removable mains connection cable, type SJT, minimal 18 AWG, 3 leads. The plug connectors must comply with NEMA 5-15 or IEC 60320. Grounding will only be reliable if the equipment is connected to a corresponding hospital grade socket.

**Potential equalization**

The equipotential bonding is used as a protective measure against the failure of the protective conductor according to requirements of IEC 60601-1 in the respectively valid version. The installation must be according to the relevant local safety regulations.


**4.3 Front of Device****Fig. 4-1 Device front**

- ① ON/OFF switch
- ② Monitor with touch screen
- ③ Footswitch connection
- ④ Gas heater connection
- ⑤ Cassette for Smoke evacuation
- ⑥ Connection for insufflation tube set
- ⑦ Eject button

Familiarize yourself with the control and function elements at the front of the device.

#### 4.4 Rear of the Device

**Fig. 4-2 Device rear**

Familiarize yourself with the connection elements at the rear of the device.

#### 4.5 Gas Connection

##### WARNING!

**Medically pure CO<sub>2</sub>**

**Make sure to use only medically pure CO<sub>2</sub>. Other gases (e.g., helium, N<sub>2</sub>O, argon), mixtures of gases, gases with entrapped liquids, or polluted gases must not be used with this device.**

##### WARNING!

**Gas supply**

**Maintain adequate gas supply at all times.**


Use a high-pressure tube to connect a CO<sub>2</sub> gas cylinder to the rear gas inlet connection or connect to a centralized CO<sub>2</sub> gas supply.

The device is equipped with a universal gas connection, which can be configured with different adapters to support a central gas supply as well as the supply of bottled gas.

The universal gas connection can be operated with a supply pressure from 2.7 bar (39 psi) to 80 bar (1160 psi).

The gas connector has an exchangeable sinter filter to protect the insufflator from dirt particles. This filter has to be checked during the annual test and replaced at least every two years. During the maintenance the sinter filter has to be exchanged. The sinter filter is available as a spare part (see chapter 14 Accessories). The sinter filter may be blocked or soiled due to unclean gas or gas supply residues. Heavy soiling is especially noticeable during higher flow rates. As a result, the actual gas flow is falling short of the nominal gas flow during active insufflation or the device is showing gas supply warnings although gas supply is sufficient. In this case the sinter filter needs to be exchanged.

##### Exchangeable sinter filter


1. Remove the connection ① using 17 mm and 19 mm open-end wrenches.
2. Unscrew the sinter filter ②; if necessary, use a screwdriver.
3. Hand-tighten the new sinter filter and reattach the connection to the device.

An angled connector is available to reduce the depth needed for the installation

##### Angled gas connector

of the device. This connector includes the exchangeable sinter filter and allows the use of the different gas supply adapters described below as well.

#### 4.5.1 Connecting a Gas Bottle

##### **WARNING!**

**Keep filled CO<sub>2</sub> bottle on hand**

**Always keep a filled CO<sub>2</sub> bottle on hand ready for replacement. This avoids having to interrupt surgery due to a lack of insufflation gas.**

##### **CAUTION!**

**Always use a high-pressure hose to connect gas bottle and device.**

The gas bottle must be in a vertical position. The gas bottle pressure may not exceed 80 bar/1160.3 psi or be less than 15 bar/217.5 psi.

##### **CAUTION!**

**Gas bottles with riser pipe may release dirt and oily fluids into the device. Do not use a gas bottle with riser pipe.**

#### High pressure hoses

The following high-pressure hoses are available:

Article	Description
PG056	High pressure hose US / bottle DIN
PG057	High pressure hose US / bottle ISO
PG058	High pressure hose US / bottle PIN index

#### Installation

##### **High pressure hose with PIN connection**

- Use a 14 mm open-end wrench to attach this connection on the device or to detach it.
- The gas bottle is screwed on or off by hand.

##### **High pressure tube with DIN connection**

- Use a 14 mm open-end wrench to attach this connection on the device or to detach it.
- The gas bottle is screwed on or off with a 30 mm open-end wrench.

##### **High pressure tube with ISO connection**

- Use a 14 mm open-end wrench to attach this connection on the device or to detach it.
- The gas bottle is screwed on or off with a 32 mm open-end wrench.

## Switching valve

Two gas bottles can be connected with the insufflator via the PG068 switching valve. This ensures continuous insufflation, especially during replacement of the gas bottle.

### 4.5.2 Connecting to Central Gas Supply

The following low pressure hoses are available for connecting to the central gas supply:

#### Low pressure hoses

Article	Description
PG061	Hose for CO <sub>2</sub> central gas supply DIN:US, 3 m
PG091	Hose for CO <sub>2</sub> central gas supply DIN:US, 5 m
PG092	Hose for CO <sub>2</sub> central gas supply NF:US, 5 m
PG093	Hose for CO <sub>2</sub> central gas supply UNI:US, 5 m
PG094	Hose for CO <sub>2</sub> central gas supply AGA:US, 5 m
PG065	Hose for CO <sub>2</sub> central gas supply DISS:DISS, 3.0 m

- Connect or loosen to the device using a 14 mm open-end wrench.
- Connect or loosen to the central gas supply receptacle on the wall.


#### Installation

An adapter and a low pressure hose is required for connecting to the central gas supply with DISS connection:

- DISS Adapter for gas connection
- DISS Hose for central gas supply

1. Remove the connector (1) using 17 mm and 19 mm open-end wrenches.
2. Check if the sinter filter (3) needs a change (residues or blocked pores).
3. Exchange the sinter filter using a screwdriver.
4. Install the new connector (2) using the open-end wrenches.

#### Exchange of the gas supply adapter

For older NIST central gas supply hoses there is a NIST adapter available:

- NIST Adapter for gas connection

Set in the user menu the type of gas supply: central gas or bottle gas supply (see chapter 7 Menu).

#### 4.6 Gas Supply

The status of the gas supply is monitored by the device and indicated with symbols and acoustic signals (beeps).

##### Gas supply with gas bottle

The following gas bottle pressures are displayed:

	> 50 bar/725.8 psi
	50 bar/725.8 psi - 30 bar/435.1 psi
	30 bar/435.1 psi - 15 bar/217.5 psi
	15 bar/217.5 psi - 5 bar /72.5 psi. Acoustic signals can be heard and the message <b>Change gas tank</b> is displayed. User is advised to prepare for changing the gas bottle.
	< 5 bar /72.5 psi. Acoustic signals can be heard and the message <b>Change gas tank</b> is displayed. Replace the gas bottle.
	0 bar / 0 psi. Acoustic signals can be heard and the message <b>Check gas supply! Insufflation paused!</b> is displayed. Immediately replace the gas bottle.



##### Central gas supply

The following symbols are displayed for the house supply gas pressures:

	<b>Central gas supply pressure OK</b>
	<b>Central gas supply pressure too low</b> Beeps are emitted and the following message appears at a supply pressure of <2.7 bar: <b>Check gas supply!</b>

#### 4.7 Switching Device On

1. Connect the device with the mains power.
2. Connect the gas supply to the gas connection port and open the gas supply.
3. Press the ON/OFF switch. The power LED lights up green when the device is turned on.

4. The device performs an automatic self-test after being switched on. If the self-test is not successful, the corresponding error message will be displayed (see chapter 12 Error and Warning Messages).
5. The language selection is displayed when starting the device the first time.

- Select the desired language and confirm with the Continue button on the top right.
6. The display depicts an overview of the insufflation indications. Press the respective function button to select the desired indication (e.g. **Standard**).

#### 4.8 Switching Device Off

Disconnect device from power supply

**NOTE!**

A safe and all-pole disconnection of the product from the power supply is only ensured by disconnecting the power cord.

j

1. Switch device off: press the ON/OFF switch (1) (Fig. 4-2 Device rear).
2. Remove the power cord from the device mains socket (14) (Fig. 4-2 Device rear).

The operation of the device is safely terminated.

#### 4.9 Demo Mode

The device has a demo mode. If the device is switched on without gas, a "Demo" button appears.

All functions are accessible; however, the message **Demo Mode** appears permanently in the status line.

Changes cannot be saved.


To exit demo mode, the device must be restarted.

## 5 Operating the Device - General

### 5.1 Monitor with Touch Screen

**Fig. 5-1 Screen displays**

- (1) Back button
- (2) Indication display
- (3) Actual pressure display
- (4) Actual gas flow display
- (5) Menu button
- (6) Message display
- (7) Smoke evacuation field with button OFF and buttons level 6 l/min and 12 l/min
- (8) Volume counter
- (9) Gas consumption display with reset button
- (10) START/STOP button
- (11) Increase nominal gas flow
- (12) Nominal gas flow
- (13) Reduce nominal gas flow
- (14) Actual gas flow bar
- (15) Increase nominal pressure
- (16) Nominal pressure
- (17) Reduce nominal pressure
- (18) Actual pressure bar
- (19) Gas supply icon
- (20) Heating state icon

The touch screen depicted above shows all display and function buttons when the insufflator has been stopped.

Additional explanations for individual elements are listed in the subsequent control element descriptions.

### 5.2 Insufflation Tube Sets

Different insufflation tube set types can be connected to the insufflation tube connection at the front of the device.

PG012	Disposable tube set incl. CO <sub>2</sub> filter, sterile
PG014	Reusable tube set for insufflation
PG096SU	Disposable heating tube set for insufflation, sterile
PG097	Reusable heating tube set for insufflation
PG019	Disposable CO <sub>2</sub> filter, sterile

The tube sets have been sterilized with ethylene oxide according to procedures validated by ISO 11135-1 and ISO 10993-7.

### 5.3 Connecting an Insufflation Tube Set

#### WARNING!

##### Visual inspection of the tube set

Before the operation, perform a visual inspection of the tube set and its packaging. Damaged tube sets or tube sets from damaged packagings may not be used.

#### WARNING!

##### Air in the insufflation tube set

##### Air in the insufflation tube set can cause an embolism.

- Start the insufflation before the tube set is connected to the patient to push the air out of the insufflation tube set. Stop the insufflation when at least 1 liter of CO<sub>2</sub> has been insufflated.

- Always connect a hydrophobic filter to the reusable tube or use the disposable tube set with filter.
- Connect the insufflation tube to the insufflation tube connection and to the Veress cannula/trocars.

#### ISO 5356-1 connection


A male ISO connector is located on the device (ISO-M). Connect the female ISO connector of the filter or tube (ISO-F) firmly with the ISO-M connector.

### 5.4 Gas Heating

It is possible to insufflate gas with a body temperature (approximately 37 °C / 99°F) with the help of tubes with heating function (see 5.2 Insufflation Tube Sets).

#### CAUTION!

**Do not subject the heater tube to direct heat (e.g., heating blanket) or high room temperatures.**

**Fig. 5-2 Connecting the gas heating**

- (1) Smoke evacuation tube
- (2) Insufflation tube
- (3) Heating cable

- Connect the filter and the heating tube or a heating tube with integrated filter.
- Connect the plug of the heating tube with the gas heating connection.

#### Connecting the gas heating

As soon as the heating cable is connected, the display (6) depicts for 3 seconds the message **Gas heating OK**.

The reusable heating tube is preheated automatically to 33 °C when insufflation is turned off.

**Gas heating ON/OFF**

Start gas heating: Press the Start/Stop button. The gas is heated automatically to approximately 37 °C within approximately 10 minutes.

**Gas heating error**

The message **Gas heating defective!** is displayed for 2 seconds in case of a function defect of the gas heater (e.g., cable flawed or broken, loose plug).

EN

Check gas heating using a different tube. Should the error message be displayed again, you can continue using the device without gas heating by observing the risks for hypothermia. Check gas heating after surgery using a different tube.

If the message **Gas heating defective!** persists, call an authorized service technician to check the device.

**Gas temperature exceeds 42°C**

The message display depicts the following displays in case of excess temperature:

**Gas heating - over-temperature!****WARNING!**

**If the temperature sensor measures a gas temperature > 42° C, unplug the heating plug from the device. Hot gas in the abdomen can lead to serious injuries.**

A warning signal is audible. Gas heating and insufflation are deactivated.

1. Disconnect the insufflation tube from the trocar or Veress needle.
2. Unplug the plug for the gas heating from the device.
3. Press the Start/Stop button. The device insufflates without heating the gas.
4. Let hot gas escape until the tube is only warm to the touch.
5. Continue operation without connecting gas heating.
6. Check gas heating after surgery using a different tube. Turn device off and on again. Gas heating is reactivated.
7. Should the error message be displayed again, you can continue using the device without gas heating by observing the risks for hypothermia.
8. Call an authorized service technician to check the gas heating.

The following gas heater states are possible and displayed on the touch screen of the device as a icon and corresponding message:

Icon	Description	Information
	Gas heating connected	Gas heating function is available.
	Gas heating defective	Insufflation can be continued in this case but the heating function is not available.  Restart the device if you want to use the heating function. If the gas heating is still defective, have an authorized service technician check or repair the device.
	Over temperature (>42°C / 108°F)	Insufflation is automatically stopped and the heating function is deactivated. There are three options to continue insufflation: <ul style="list-style-type: none"> <li>• Wait until the temperature has decreased below 43°C, insufflation and heating will resume automatically.</li> <li>• Disconnect and reconnect tube set. If temperature is below 43°C insufflation with heating can be resumed (manual start).</li> <li>• Exchange tube set with a new one.</li> </ul>

## 5.5 Using the Smoke Evacuation Function

### WARNING!

**Increased concentration of CO<sub>2</sub> in the operating room**

**Smoke evacuation function may increase the CO<sub>2</sub> concentration in the operating room.**

- Ensure sufficient supply of fresh air and use the smoke extraction function only if it is necessary.

The device is equipped with a vacuum pump. The vacuum pump is used for smoke evacuation by means of a disposable tube set and a smoke evacuation filter cassette.

**Fig. 5-3 Smoke evacuation display**

- (7) Smoke evacuation levels  
6 l/min and 12 l/min

An additional feature of the device is smoke evacuation from the body cavity. For this purpose, use a disposable tube set for the smoke evacuation and a smoke evacuation filter cassette.

PG106 Smoke evacuation filter cassette (hydrophobic and activated charcoal filter for the smoke filtering)


PG107SU Disposable tube set for smoke evacuation, sterile with filter, which serves to remove moisture or liquid from the evacuated gas.

To enable this function, press the 6 l / min or 12 l / min button in the **Smoke Evacuation** field depending on the intensity of evacuation that you need. Smoke evacuation can be set manually to 6 l / min or 12 l / min.

The **Smoke Evacuation** function can be switched on with the foot control switch or switched off (see chapter 14 Accessories). Connect the foot control switch to the connection for the foot control switch on the front panel (see Fig. 4-1 Device front).

The nominal gas flow setting must be higher than the smoke evacuation value; otherwise the efficiency of the evacuation is limited.

The smoke evacuation function is not available during the cardio indication.

### 5.5.1 Replacing the Smoke Evacuation Filter Cassette

### WARNING!


**Visual inspection of the filter cassette**

**Before the operation, perform a visual inspection of the filter cassette. Damaged filter cassette may not be used.**

To insert the filter cassette, align the bottom of the filter cassette on the under-

side of the socket on the front of the insufflator (See Figure A). The smoke evacuation filter cassette is intended for use in multiple surgeries and will remain in the device. The service life of the smoke evacuation filter cassette depends on the quantity of smoke produced during multiple surgeries. Before the smoke evacuation filter cassette is completely clogged, the message "Smoke evacuation will be paused soon - be prepared for filter cassette change" will appear on the device. If the smoke evacuation filter cassette is completely clogged, the smoke evacuation from and to the device turns off and the device depicts the message "Smoke evacuation paused - Change filter cassette." To replace the smoke evacuation cassette, unplug the disposable tube set for the smoke evacuation and press the eject button for the smoke evacuation filter cassette (see Figure D).

If a disposable tube set is not connected to the smoke evacuation, the smoke evacuation filter cassette is to be closed with the protective cap to retain any odor that may emanate from the filter.

### 5.5.2 Replacing the Disposable Smoke Evacuation Tube Set

The disposable tube set used for the smoke evacuation is delivered sterile. It may only be used once for a patient and must be discarded after each surgery. To connect the disposable smoke evacuation tube set to the filter cassette, remove the protective cap from the smoke evacuation filter cassette and plug it in (See Figure B).

To remove the disposable smoke evacuation tube set, unplug it from the smoke evacuation filter cassette (See Figure C).

### 5.6 Displaying/Selecting Insufflation Indication

The display depicts an overview of the insufflation indications. Press the respective function button to select the desired indication (e.g. Standard).

The displayed parameters correspond with the factory settings or the values set in the user menu (see chapter 7 Menu).

## 5.7 Setting the Nominal Pressure - All Indications Except Fieldflooding

To set the nominal pressure, press the (+) or (-) button on the display below where the actual pressure is depicted.

- The nominal pressure can be increased or decreased with running or stopped insufflation.
- The nominal pressure is increased/reduced in steps of 1 mmHg every time the + or - button is pressed. If the + or - button longer is pressed longer than 1.5 seconds, fast scrolling is activated, but only to the safety threshold, which depends on the particular indication (see below).
- When setting the nominal pressure, an indication-dependent safety threshold of 12 mmHg or 15 mmHg is applied. If the pressure is to be increased above these values, the depicted message on the insufflator display must be confirmed accordingly.

### CAUTION!

**Exceeding this safety threshold is to be decided by and the responsibility of the user/operator.**


Indication	Modes	Safety threshold
Standard Indication	Veress	15 mmHg
	Highflow	15 mmHg
Bariatric Indication	Veress	15 mmHg
	Highflow	15 mmHg
Pediatric Indication	Sensitive	12 mmHg
Cardio Indication	Fieldflooding	Not applicable
	Vessel Harvesting	12 mmHg

## 5.8 Setting the Nominal Gas Flow Rate - All Indications

To set the nominal gas flow rate, press the + or - button on the display below where the actual gas flow is depicted.

- The nominal gas flow is increased/reduced in steps of 1 l/min every time the + or - button is pressed. If you are working in the Pediatric indication with a range of 0.1 to 2 l/min, the value is increased or decreased by 0.1 l/min. Hold the - or + button for 1.5 seconds to select one of three flow rates, but only up to the safety threshold, which depends on the respective indication (see below). Please refer to chapter 7 Menu for details on the factory setting. The preset values of the user menu can be customized as well.
- In Pediatric indication, the safety threshold is at 5 l/min, in Vessel Harvesting mode at 6 l/min. If the gas flow is to be increased above these values, the depicted message must be confirmed accordingly.

### CAUTION!

**Exceeding this safety threshold is to be decided by and the responsibility of the user/operator.**

Indication	Modes	Safety threshold
Standard Indication	Veress	Not applicable
	Highflow	Not applicable
Bariatric Indication	Veress	Not applicable
	Highflow	Not applicable
Pediatric Indication	Sensitive	5 l/min
Cardio Indication	Fieldflooding	Not applicable
	Vessel Harvesting	6 l/min

### 5.8.1 Veress Insufflation

**Veress insufflation** is a gentle type of insufflation that prevents exceeding the preset nominal pressure by the actual pressure even in cases of small volumes. To minimize the risks of a faulty incision, the manufacturer recommends using the **Veress insufflation** mode to start a procedure (filling abdomen with CO<sub>2</sub>).

#### CAUTION!

Please note that the automatic venting system is not active during Veress insufflation.

#### Activating Veress insufflation

The activation of Veress Insufflation depends on the chosen indication.

**Bariatric, Standard indication:** With a nominal flow setting of  $\leq 5$  l/min, **Veress insufflation** is displayed after pressing the Start/Stop switch.

**Pediatric indication:** With a nominal flow setting of  $\leq 1$  l/min, **Veress insufflation** is displayed after pressing the Start/Stop switch.

**Cardio indication:** In this indication, Veress insufflation is not available.

When setting higher nominal gas flow values, highflow mode is automatically activated (see chapter 5.8.2 Highflow Insufflation).

### 5.8.2 Highflow Insufflation

**Highflow mode** is the most common insufflation type in Standard and Bariatric indications. Leaks are compensated quickly, the required nominal pressure is reached as fast as possible.

The Highflow mode should be selected if the surgery is not started with a Veress needle, e.g., when using Hasson trocars.

### 5.8.3 Profiles

The Profiles mode offers the possibility to select previously defined user-specific profiles (see chapter 7.1 Profiles on how to create profiles). A profile is selected after selecting the indication on the Standard, Bariatric, Pediatric or Cardio pages.

## 5.9 Function Test before Device Use

Check all the sterile disposable items before removing them from the package to ensure that the packaging is intact and that the expiration date is still valid.

### Preparation

For your own safety and that of your patient, use only original accessories (see chapter 14 Accessories).

#### **WARNING!**

**The function test must be performed prior to each surgical procedure.**

#### **WARNING!**

**For the safety of the patient please fill the tube set with CO<sub>2</sub> gas prior to beginning the insufflation by activating the insufflation for a few seconds and then turning it off again before introducing the insufflation instrument to the cavity and beginning the surgery.**

#### **WARNING!**

**Do not attach any cannula with tube set to patient during the function test.**

### Checking the device

#### **NOTE!**

**This test can be conducted with the tube set and cannula that will be used for the procedure. Consider sterile/non sterile aspects during handling.**

1. The device is switched off, no tube set is connected.
2. Make sure the gas supply is connected and open.
3. Use the ON/OFF switch to turn the device on. The device now conducts an automatic self-test.
4. Select the indication (e.g. Pediatric).
5. Connect a new insufflation tube set to the device.
6. Attach the insufflation tube to the cannula and close the stopcock.
7. If the gas consumption display field depicts a value, press the function button to reset the display to zero **0 L**.
8. Select a nominal pressure of 8 mmHg and a gas flow of 2 l/min.
9. Start insufflation: Press the **START** button.
10. Insufflate approximately 30 s. The display status line depicts **Occlusion** after maximum 4 seconds.
11. Stop insufflation: Press the function button **STOP** (10) Fig. 5-1 Screen displays.

#### **WARNING!**

**If the actual gas consumption is higher than 0.4 l, there is a leak in the system. If this is the case, use steps 12 to 14 outlined below to locate the leak.**

12. Repeat steps 7 to 11 without cannula and with closed tube end. The previously connected cannula has a leak if gas consumption is now below 0.4 l.
13. Repeat steps 8 to 12 with a new tube if another leak becomes apparent. If the gas consumption is now below 0.4 l, the previously used tube set has a leak.
14. If another leak is detected, this leak is then directly within the device. Make sure the device can no longer be operated until an authorized service technician conducts the appropriate checks.

**WARNING!**

**Never work with a leaky insufflation tube, accessory, and/or device. This can lead to an incorrect measurement of the actual pressure values, which can cause an uncontrolled pressure increase in the abdomen.**

EN

Choose desired indication, pressure and flow settings and continue with procedure.

### 5.10 Gas Consumption Display

The gas consumption display indicates the insufflated volume of CO<sub>2</sub> in liters consumed since the last resetting of the display.

#### Resetting the display

The gas consumption display can be reset to **0.0** during running or stopped insufflation.

- ① Press the button to reset the gas consumption display.

### 5.11 Starting/Stopping Insufflation

Depending on the insufflation status, either the **START** or the **STOP** button is displayed.

#### Starting insufflation:

The **START** key is displayed when insufflation is stopped.

Press this key to start insufflation.

The actual gas flow is depicted numerically and with a graphical bar. Here, the actual gas flow rate will range between zero and the set nominal gas flow rate. The actual pressure in the body cavity is depicted numerically and with a graphical bar. In this case, the actual pressure starting at zero approaches the set nominal pressure. Once the actual pressure in the body cavity reaches the nominal pressure, the actual gas flow rate is almost zero, since only the amount of gas required to maintain the pressure must be refilled.

#### Stopping insufflation:

The **STOP** button is displayed while insufflation is running.

- ② Stop insufflation with the **STOP** button.

The actual gas flow rate display will show zero. The actual pressure display continues to show the current pressure in abdomen or cavity.

## 6 Indications of the Flow 50 Insufflator

Indications Overview

	Standard Indication		Bariatric Indication		Pediatric Indication	Cardio Indication	
Mode	Veress	Highflow	Veress	Highflow	Sensitive	Field-flooding	Vessel Harvesting
Gas flow rate settings l/min	1-5	6-40	1-5	6-50	0.1-20	1-10	1-10
Gas flow increment (l/min)	1	1	1	1	0.1 for < 2 l/min 0.5 for 2 to 10 l/min 1 for > 10 l/min	1	1
Gas flow safety threshold (l/min)	Not applicable	Not applicable	Not applicable	Not applicable	5 l/min	Not applicable	6 l/min
Starting pressure	15	15	15	15	3	Not applicable	10
Pressure settings (mmHg)	1-30	1-30	1-30	1-30	1-15	Not applicable	1-20
Pressure change in (mmHg)	1	1	1	1	1	Not applicable	1
Pressure safety threshold (mmHg)	15	15	15	15	12	Not applicable	12

### 6.1 Standard Indication

The **Standard** indication is used to distend the peritoneal cavity during laparoscopic procedures by insufflating CO<sub>2</sub>.

Depending on the surgery, the user can choose one of the following modes. For further details see also chapter 5.8 Setting the Nominal Gas Flow Rate - All Indications.

- Veress
- Highflow

Alternatively, the user can work with a previously created user profile (see chapter 7.1 Profiles).

The **Standard** indication is designed for laparoscopies performed on normal weight and slightly obese (BMI < 30 kg/m<sup>2</sup>) patients over the age of 14.

The maximum nominal pressure can be set to 30 mmHg and the maximum gas flow to 40 l/min.

**Patients**

**Insufflation parameters**

#### WARNING!

##### Idiosyncratic reactions

**Patients with sickle cell anemia or pulmonary insufficiency may have a higher risk of metabolic imbalance related to excessive CO<sub>2</sub> absorption (idiosyncratic reaction).**

**WARNING!****CO<sub>2</sub> absorption**

CO<sub>2</sub> is absorbed during insufflation (intravasation). This means the body absorbs part of the CO<sub>2</sub> gas used for insufflation. CO<sub>2</sub> concentrations in the blood or respiratory system that are too high can lead to death of the patient in extreme cases. To lower this risk, always carefully and closely monitor the patient's vital signs during the entire insufflation process and make sure patient is breathing well. Sufficient respiration can help avoid or limit problems with CO<sub>2</sub>. High pressure or a high gas flow promotes CO<sub>2</sub> absorption. The abdomen is sufficiently distended using a pressure between 10 to 15 mmHg. Pressure values above 15 mmHg are required for only a few cases but do increase the risk of intravasation. Never exceed the maximum intra-abdominal pressure of 30 mmHg.

**WARNING!****Metabolic and cardiac reactions**

Insufflating CO<sub>2</sub> may result in metabolic acidosis. This can lead to cardiac irregularities expressed with the following symptoms:

- Reduced respiration with restricted diaphragm function
- Hypercapnia
- Reduction of venous reflux
- Reduced cardiac output
- Metabolic acidosis

**WARNING!****Hypothermia/monitoring body temperature**

The gas flow can lead to a lowering of the patient's body temperature during insufflation. Hypothermia during insufflation can cause heart and cardiovascular problems. The risk for hypothermia can be significantly reduced with the use of gas that is pre-warmed to body temperature. Therefore, you must always monitor the patient's body temperature during the entire insufflation process. Make especially sure that the following, hypothermia-fostering, surgical conditions are avoided as best as possible:

- High gas flow due to large leaks
- Long surgeries
- Use of irrigation and infusion solutions that are not preheated to body temperature.

**WARNING!****Dehydration**

Insufflation can lead to dehydration of the tissue. This can result in organ tissue damage and cardiovascular reactions of the patient. Long surgeries and large leaks increase the risk of dehydration (especially at the insertion points of the trocars or when changing instruments).

**WARNING!****Embolism/insufflation of internal organs**

Improper placement of the insufflation instrument could cause gas penetrating a vessel or an internal organ, resulting in air or CO<sub>2</sub> embolisms. To reduce the risk, use a low flow rate for the first insufflation and ensure that the insufflation instrument is correctly positioned. Check the position of the insufflation instrument immediately if the actual pressure rapidly reaches the nominal pressure value. CO<sub>2</sub> embolisms can also be caused by a high intra-abdominal pressure. Avoid high-pressure settings and close damaged blood vessels at once.

**WARNING!****Subcutaneous emphysema**

Incorrect placement of a cannula or a trocar into subcutaneous tissue may lead to emphysema. To reduce the risk, use a low gas flow rate for the first insufflation and ensure that the insufflation instrument is correctly positioned. Long surgeries, the use of many access points, duration and size of leaks at these points may also contribute to emphysema. Be sure to close leakages in trocar accesses immediately.

**WARNING!****Additional insufflation sources/automatic venting system**

Make sure the automatic venting system is activated (see chapter 7 Menu). The use of additional insufflation sources increases the intra-abdominal pressure. Continuously monitor intra-abdominal pressure over the course of the entire insufflation if additional sources are used.

**6.2 Bariatric Indication**

The **Bariatric** indication is used to distend the peritoneal cavity during laparoscopic procedures by insufflating CO<sub>2</sub>.

Depending on the surgery, the user can choose one of the following modes. For further details see also chapter 5.8 Setting the Nominal Gas Flow Rate - All Indications.

- Veress
- Highflow

Alternatively, the user can work with a previously created user profile (see chapter 7.1 Profiles).

The **Bariatric** indication is designed for laparoscopies performed on normal weight and overweight patients over the age of 14. **Patients**

The maximum nominal pressure can be set to 30 mmHg and the maximum gas flow to 50 l/min. **Insufflation parameters**

**WARNING!****Idiosyncratic reactions**

**Patients with sickle cell anemia or pulmonary insufficiency may have a higher risk of metabolic imbalance related to excessive CO<sub>2</sub> absorption (idiosyncratic reaction).**

**WARNING!****CO<sub>2</sub> absorption**

CO<sub>2</sub> is absorbed during insufflation (intravasation). This means the body absorbs part of the CO<sub>2</sub> gas used for insufflation. CO<sub>2</sub> concentrations in the blood or respiratory system that are too high can lead to death of the patient in extreme cases. To lower this risk, always carefully and closely monitor the patient's vital signs during the entire insufflation process and make sure patient is breathing well. Sufficient respiration can help avoid or limit problems with CO<sub>2</sub>. High pressure or a high gas flow promotes CO<sub>2</sub> absorption. The abdomen is sufficiently distended using a pressure between 10 to 15 mmHg. Pressure values above 15 mmHg are required for only a few cases but do increase the risk of intravasation. Never exceed the maximum intra-abdominal pressure of 30 mmHg.

**WARNING!****Metabolic and cardiac reactions**

Insufflating CO<sub>2</sub> may result in metabolic acidosis. This can lead to cardiac irregularities expressed with the following symptoms:

- Reduced respiration with restricted diaphragm function
- Hypercapnia
- Reduction of venous reflux
- Reduced cardiac output
- Metabolic acidosis

**WARNING!****Hypothermia/monitoring body temperature**

The gas flow can lead to a lowering of the patient's body temperature during insufflation. Hypothermia during insufflation can cause heart and cardiovascular problems. The risk for hypothermia can be significantly reduced with the use of gas that is pre-warmed to body temperature. Therefore, you must always monitor the patient's body temperature during the entire insufflation process. Make especially sure that the following, hypothermia-fostering, surgical conditions are avoided as best as possible:

- High gas flow due to large leaks
- Long surgeries
- Use of irrigation and infusion solutions that are not preheated to body temperature.

**WARNING!****Dehydration**

Insufflation can lead to dehydration of the tissue. This can result in organ tissue damage and cardiovascular reactions of the patient. Long surgeries and large leaks increase the risk of dehydration (especially at the insertion points of the trocars or when changing instruments).

**WARNING!****Embolism/insufflation of internal organs**

Improper placement of the insufflation instrument could cause gas penetrating a vessel or an internal organ, resulting in air or CO<sub>2</sub> embolisms. To reduce the risk, use a low flow rate for the first insufflation and ensure that the insufflation instrument is correctly positioned. Check the position of the insufflation instrument immediately if the actual pressure rapidly reaches the nominal pressure value. CO<sub>2</sub> embolisms can also be caused by a high intra-abdominal pressure. Avoid high-pressure settings and close damaged blood vessels at once.

**WARNING!****Subcutaneous emphysema**

Incorrect placement of a cannula or a trocar into subcutaneous tissue may lead to emphysema. To reduce the risk, use a low gas flow rate for the first insufflation and ensure that the insufflation instrument is correctly positioned. Long surgeries, the use of many access points, duration and size of leaks at these points may also contribute to emphysema. Be sure to close leakages in trocar accesses immediately.

**WARNING!****Additional insufflation sources/automatic venting system**

Make sure the automatic venting system is activated (see chapter 7 Menu). The use of additional insufflation sources increases the intra-abdominal pressure. Continuously monitor intra-abdominal pressure over the course of the entire insufflation if additional sources are used.

**WARNING!****Altered Respiratory Physiology**

Always monitor the patient's respiratory functions during the entire surgery. The larger body mass supported by the thoracic cage and the larger amount of fat in the abdominal cavity may reduce the elasticity of the thoracic wall. In addition, the increased intra-abdominal pressure secondary to insufflation may alter the normal physiological lung parameters thus resulting in a reduction of the functional lung volume. Shallow, rapid breathing is symptomatic of this condition. Even modest physical stress causes a tremendous increased demand for oxygen, which stands in contrast to the ineffective respiratory musculature that requires more oxygen because it must overcome the reduced elasticity of the thoracic cage. The functional capacity of the lungs is small and even moderate stress can lead to respiratory failure.

**WARNING!****Subcutaneous emphysema**

When puncturing the thicker abdominal wall of morbidly obese patients with the Veress cannula or the trocar, carefully monitor the correct position of the instrument in the abdomen.

**WARNING!****CO<sub>2</sub> supersaturation**

To avoid generating CO<sub>2</sub> supersaturation, an increased level of respiratory activity is required. An overweight patient's oxygen demand and carbon dioxide production are greater and increase faster under physical stress than do those of patients with normal body weight.

**WARNING!****Heart and cardiovascular insufficiency**

Constantly monitor all heart and cardiovascular parameters during surgery since morbidly obese patients have an increased risk of heart and cardiovascular insufficiencies.

### 6.3 Pediatric Indication

The **Pediatric** indication is used to distend the peritoneal cavity during laparoscopic procedures by insufflating CO<sub>2</sub>.

The user can select the following mode:

- Sensitive

Alternatively, the user can work with a previously created user profile (see chapter 7.1 Profiles).

#### Patients

The **Pediatric** indication is designed specifically for use on newborns, infants, and children.

#### Insufflation parameters

The maximum nominal pressure can be set to 15 mmHg and the maximum gas flow rate to 20 l/min.

For use in children, it is recommended to optimize the device control (depending on the nominal flow selection) based on age and weight as follows:

Age group	Weight	Gas flow range
Children younger than 1 year	approximately 1-9 kg	0.1-0.5 l/min
Children from 1 to 3 years	approximately 10-15 kg	0.5-1.0 l/min
Children from 3 to 4 years	approximately 16-19 kg	1.0-2.0 l/min
Children from 4 to 14 years	> 20 kg	> 2.0 l/min

#### WARNING!

##### Recommended work settings

The flow values listed for laparoscopic procedures performed on newborns, infants, and children are only suggested values. The selection of the suitable flow and pressure values is solely the responsibility of the attending physician.

#### WARNING!

##### Gas flow limit

The gas flow may not exceed 14 l/min when performing a laparoscopy on newborns or patients weighing less than 25 kg (approximately 55 US pounds).

#### WARNING!

##### Pneumolabium/pneumoscrotum

Children are at risk of a pneumolabium or pneumoscrotum.

#### NOTE!

##### Low nominal flow setting

The nominal pressure cannot be reached in case of leaks if the nominal flow is set too low.

**WARNING!**

Only specially trained and qualified personnel may use this device on children.

**WARNING!****Increased airway pressure**

When laparoscopic procedures are performed on children, the increased intra-abdominal pressure also increases the risk for higher airway pressures. Always strictly monitor respiration and airway function when performing laparoscopic procedures on children younger than 12 years of age.

**WARNING!****Compression of the vena cava**

When insufflating the abdomen of a child with medical CO<sub>2</sub>, an increased risk of compressing the vena cava exists. This risk can be reduced by monitoring the systolic and diastolic blood pressure during the entire surgery.

**WARNING!****Haemodynamic stability**

A laparoscopy performed on children younger than 12 years of age may result in problems of the haemodynamic system due to the CO<sub>2</sub> content in the blood. It is recommended to increase the breathing rate of the patient and to work with low flow values and pressure values not exceeding 12 mmHg. The patient's circulatory system should be monitored at all times.

**WARNING!****Hypothermia**

The insufflation gas flow usually drops significantly after the target pressure has been reached and it is then only required to maintain the abdominal pressure. However, leaks within the abdomen or the instrument can lead to a constant gas flow of above 1 l/min. When operating on children younger than 12, a gas flow of more than 1 l/min poses an increased risk of hypothermia for the patient. Corresponding measures to prevent hypothermia include the use of blankets or pre-warmed gas. The patient's body temperature has to be monitored at all times during surgery.

**WARNING!**

Laparoscopy with CO<sub>2</sub> should not be performed on children with cardiovascular problems.

**WARNING!**

Because pediatric patients are especially susceptible to hypercapnia, it is recommended to establish an end-tidal CO<sub>2</sub> monitoring routine.

**Patients****6.4 Cardio Indication**

The **Cardiac** indication is suitable for adult patients over the age of 14.

Depending on the surgery, the user can choose one of the following modes. For further details see also chapter 5.8 Setting the Nominal Gas Flow Rate - All Indications.

- Fieldflooding
- Vessel Harvesting

Alternatively, the user can work with a previously created user profile (see chapter 5.8.3 Profiles).

**6.4.1 Fieldflooding**

The **Cardiac** indication is used to replace the ambient air from the surgical site during open and endoscopic-assisted heart surgery. In this indication, no distention is carried out.

**Insufflation parameters**

Pressure is not controlled in Fieldflooding mode. A controlled flow of CO<sub>2</sub> gas is produced. The max. **gas flow rate** is 10 l/min.

**WARNING!**

**A gas flow of 10 l/min should not be used longer than up to 30 seconds.**

**WARNING!**

**When using accessories not specified in the user manual, make sure the CO<sub>2</sub> is not directly applied to soft tissue focusses.**

**WARNING!****Fieldflooding**

**In case of CO<sub>2</sub> flooding of the surgical field (so-called fieldflooding), there is no guarantee of an CO<sub>2</sub> environment.**

**CAUTION!****Flow during fieldflooding**

**For CO<sub>2</sub> flooding of the surgical field (so-called fieldflooding, a continuous flow of 5 l/min for the entire cross-clamp time is recommended.**

**CO<sub>2</sub> diffuser**

When using Cardiac indication to push ambient air from the chest, for example, it is recommended to use a CO<sub>2</sub> diffuser to improve the CO<sub>2</sub> environment.

**6.4.2 Vessel Harvesting**

The **Vessel Harvesting** procedure is used to distend the cavity along the vena saphena magna and radial artery during an endoscopic vessel harvesting procedure. Carrying out the surgery for vessel harvesting requires the use of a special instrument.

**Patients**

The **Vessel Harvesting** procedure is suitable for adult patients over the age of 14.

**Insufflation parameters**

The maximum nominal pressure can be set to 20 mmHg and the maximum gas flow to 10 l/min.

**WARNING!****Idiosyncratic reactions**

**Patients with sickle cell anemia or pulmonary insufficiency may have a higher risk of metabolic imbalance related to excessive CO<sub>2</sub> absorption (idiosyncratic reaction).**

**WARNING!****CO<sub>2</sub> absorption**

**Due to the special surgical procedures - start of the heart bypass operation, and the endoscopic removal of the vessel - special care has to be taken as CO<sub>2</sub> is always absorbed through the tissue of the patient during insufflation (intravasation). This means the body absorbs part of the CO<sub>2</sub> gas used for insufflation. CO<sub>2</sub> concentrations in the blood or respiratory system that are too high can lead to death of the patient in extreme cases. To lower this risk, always carefully and closely monitor the patient's vital signs during the entire insufflation process and make sure patient is breathing well. Sufficient respiration can help avoid or limit problems with CO<sub>2</sub>. High pressure or a high gas flow promotes CO<sub>2</sub> absorption.**

**WARNING!****Metabolic and cardiac reactions**

**Insufflating CO<sub>2</sub> may result in metabolic acidosis. This can lead to cardiac irregularities expressed with the following symptoms:**

- Reduced respiration with restricted diaphragm function
- Hypercapnia
- Reduction of venous reflux
- Reduced cardiac output
- Metabolic acidosis

**WARNING!****Dehydration**

**Insufflation can lead to dehydration of the tissue. This can result in organ tissue damage and cardiovascular reactions of the patient. Long surgeries and large leaks increase the risk of dehydration (especially at the insertion points of the trocars or when changing instruments).**

**WARNING!****Embolism/insufflation of internal organs**

**Improper placement of the insufflation instrument could cause gas penetrating a vessel or an internal organ, resulting in air or CO<sub>2</sub> embolisms. To reduce the risk, use a low flow rate for the first insufflation and ensure that the insufflation instrument is correctly positioned. Check the position of the insufflation instrument immediately if the actual pressure rapidly reaches the nominal pressure value. CO<sub>2</sub> embolisms can also be caused by a high intra-abdominal pressure. Avoid high-pressure settings and close damaged blood vessels at once.**

EN

**WARNING!**

**Only those who are specially trained and qualified for procedures with children or endoscopic vessel harvesting procedures may use this device for these purposes.**

---

**WARNING!**

**Before using the insufflator to endoscopically harvest vessels, please check whether the instrument used is intended for CO<sub>2</sub> insufflation.**

---

**WARNING!**

**Pneumoperitoneum**

**When a vessel is harvested from the leg of a patient with a perforated groin, it is possible for CO<sub>2</sub> to reach the abdomen and cause a pneumoperitoneum. Make sure the abdomen does not fill with CO<sub>2</sub> during surgery.**


---

## 7 Menu

The following page explains access to the user menu and how to utilize this menu.

1. Press the MENU button to access the user menu.
2. Select the desired menu item in the user menu.

**Fig. 7-1 Profiles display**

The "Profiles" menu item is used to create a personal profile for the user. The preferences for the starting pressure, start gas flow, brightness, volume, and timer for the smoke evacuation system can be stored. You can also create profiles for the indications: Standard, Bariatric, Pediatric, Fieldflooding, and Vessel Harvesting. These settings have priority over the default settings for this indication.

Selecting a profile takes place after selecting the indication on the selection pages Standard, Bariatric, Pediatric or Cardio (see chapter 5.8.3 Profiles).

**Starting pressure:** Nominal pressure to be reached at the start of insufflation (setting range is based on the indication).

**Starting gas flow rate:** Nominal gas flow rate to be reached at the start of insufflation (setting range is based on the indication).

**Brightness:** Brightness of the user's screen (20 %, 40 %, 60 %, 80 %, 100 %).

**Volume:** Volume of warning sounds (20 %, 40 %, 60 %, 80 %, 100 %).

**Smoke evac timer:** Time interval after which the smoke evacuation automatically is switched from 12 l/min to 6 l/min (1 to 30 min or indefinitely).

### 7.1.1 Creating Profiles



1. Select the menu item "Profiles" in the user menu.
2. In the Profiles menu, press the (+) button (see ① Fig. 7-1 Profiles display).
3. Enter the profile name using the screen keyboard and continue with the arrow key (->) (see ② Fig. 7-2 Enter profile name).
4. Select the profile settings and save with the save button (see ③ Fig. 7-3 Details display). Use the arrow key (<-) to cancel input at any time.
5. Use the arrow key (<-) to jump to the previous menu.

### 7.1.2 Change Profile

1. Select the menu item "Profiles" in the user menu.
2. If a profile is already selected, the selected profile will have a black background.
3. In the Profiles menu, select the profile to be changed using the Shift key (see ④ Fig. 7-1 Profiles display).
4. Change the profile name, if desired use the on-screen keyboard, and proceed with confirmation using the arrow key (->) (see ② Fig. 7-2 Enter profile name).
5. Change the profile settings, if desired, and save with the save button (see ③ Fig. 7-3 Details display). Use the arrow key (<-) to cancel input at any time.
6. Use the arrow key (<-) to jump to the previous menu.

### 7.1.3 Delete Profile

1. Select the menu item "Profiles" in the user menu.
2. If a profile is already selected, the selected profile will have a black background.
3. In the Profiles menu, select the profile to be deleted using the (-) key (see ⑤ Fig. 7-1 Profiles display).
4. Confirm deletion of the profile with "Delete", or cancel the operation.
5. Use the arrow key (<-) to jump to the previous menu.

**Fig. 7-2 Enter profile name****Fig. 7-3 Details display**

## 7.2 Indication Setting

### 7.2.1 Insufflation

The menu item "Insufflation" is used to adjust the values for the starting pressure, the gas flow rates, the venting system, the maximum pressure, and the occlusion alert. The settings are made for the currently selected indication and the currently selected mode. However, the starting pressure set for a profile is applied when selecting that profile.

Starting pressure: Nominal pressure to be reached at the start of insufflation (setting range is based on the indication).

Gas flow rates: Keeping the (-) or (+) button depressed for 1.5 seconds allows you to select one of 3 levels for the nominal gas flow. The setting ranges for these gas flow rates depend on the indication.

Venting system: The venting valve is opened to relieve excess pressure in the patient. The venting vent can be turned off completely or only during Veress insufflation. The venting valve opens with a time delay after the actual pressure has exceeded the nominal pressure for a specific time. This time, and this pressure difference, can be set here.

Maximum pressure: The maximum pressure adjustable in the OR screen. The setting range for the maximum pressure depends on the indication.

Occlusion alert: The acoustic signal of the occlusion alert can be turned on or off here.

### 7.2.2 Smoke Evacuation

The "Smoke Evacuation" menu item is used to adjust the values for timer and level of the smoke evacuation by means of a foot control switch.

Smoke evac timer: Time interval after which the smoke evacuation automatically is switched from 12 l/min to 6 l/min (1 to 30 min or indefinitely).

Foot control level: Level of the smoke evacuation turned on by pressing the foot pedal (6 l/min or 12 l/min).

## 7.3 System

The "System" menu item is used to adjust the values for volume, brightness, gas supply, and language and as well as the values to reset the settings to the factory defaults.

The settings can be done for all indications and modes (except Profiles).


General volume - volume of beeps (20 %, 40 %, 60 %, 80 %, 100 %), factory setting 80 %.

General brightness - brightness of the display (20 %, 40 %, 60 %, 80 %, 100 %), factory setting 80 %.


Gas supply: House gas supply or bottle gas supply (factory setting is house gas supply). When a bottle is connected, the device automatically switches to bottled gas supply.

Language: Language of displays on the user's screen, factory setting is English.

Factory reset: All settings are reset to the factory settings. Resetting to the factory settings must be confirmed.

**Fig. 7-4 System**

## 7.4 Service

Under the "Service" menu item, the current software and firmware versions as well as the dates of the next maintenance will be shown. Moreover, the access for the authorized service technician is here as well.

Software: The current software version number.

Firmware: The current firmware version number.

Next maintenance: Date on which the next biannual maintenance becomes due.

Technician: password protected access for authorized service technicians. Please consult the service manual for additional information.

## 7.5 Protocol

**Fig. 7-5 Protocol display**

Under the "Protocol" menu item, the last 10 protocols can be displayed, exported to a memory stick, or deleted.

A list of the protocols of the last 10 surgeries is depicted. The protocols are sorted chronologically and the protocol of the most recent surgery is listed at the top. Select a protocol from the list to display the information about the surgery.

The following are displayed: indications used or the selected user profile; the duration of surgery; the insufflation duration; smoke evac duration; minimum, maximum and starting value for pressure and gas flow; the gas consumption; and the gas heating details.

Use the arrow key (<-) to jump back to the protocol list.

Individual protocols can be deleted with the (-) button (see ① Fig. 7-5 Protocol display) or exported to a memory stick with the save button (see ② Fig. 7-5 Protocol display). To export protocols to a memory stick, a memory stick must first be *inserted into the memory stick port* Located on the device rear panel (see chapter 4.4 Rear of the Device).

## 7.6 User Menu (Overview)

The following settings can be made in this menu.

**The following settings can be made under indication settings:**

User Menu	Indications	Modes	Options	Factory Settings
Setting the start pressure	Standard		1-15 mmHg	15 mmHg
	Bariatric		1-15 mmHg	15 mmHg
	Pediatric	Sensitive	1-12 mmHg	3 mmHg
	Cardiac	Vessel Harvesting	1-12 mmHg	10 mmHg
		Fieldflooding	Pressure not adjustable	Not applicable

**The following settings can be made under indication settings:**

User Menu	Indications	Modes	Options	Factory Settings
<b>Gas flow rates</b> (depend on indication)  Keeping the (-) or (+) button depressed for 1.5 seconds allows you to select one of 3 levels for the nominal gas flow.	Standard	Veress	Gas Flow Rate 1: 1-5 l/min	Gas Flow Rate 1: 3 l/min
		Highflow	Gas Flow Rate 2: 6-30 l/min Gas Flow Rate 3: 20-40 l/min	Gas Flow Rate 2: 20 l/min Gas Flow Rate 3: 40 l/min
	Bariatric	Veress	Gas Flow Rate 1: 1-5 l/min	Gas Flow Rate 1: 3 l/min
		Highflow	Gas Flow Rate 2: 6-40 l/min Gas Flow Rate 3: 20-50 l/min	Gas Flow Rate 2: 25 l/min Gas Flow Rate 3: 50 l/min
	Pediatric	Sensitive	Gas Flow Rate 1: 0.1-2 l/min Gas Flow Rate 2: 1-5 l/min Gas Flow Rate 3: 5-20 l/min	Gas Flow Rate 1: 1 l/min Gas Flow Rate 2: 3 l/min Gas Flow Rate 3: 10 l/min
		Vessel Harvesting	Gas Flow Rate 1: 1-10 l/min Gas Flow Rate 2: 1-10 l/min Gas Flow Rate 3: 5-10 l/min	Gas Flow Rate 1: 2 l/min Gas Flow Rate 2: 5 l/min Gas Flow Rate 3: 10 l/min
	Cardiac	Fieldflooding	Not applicable	
<b>Venting valve</b> (depends on indication)	Standard	Veress	- ON - OFF	OFF
		Highflow	- ON - OFF	ON
	Bariatric	Veress	- ON - OFF	OFF
		Highflow	- ON - OFF	ON
	Pediatric	Sensitive	- ON - OFF	ON
		Vessel Harvesting	- ON - OFF	OFF
	Cardiac	Fieldflooding	Not applicable	
<b>Max. pressure</b> (depends on indication)	Standard	Veress	5-30 mmHg	30 mmHg
		Highflow	5-30 mmHg	30 mmHg
	Bariatric	Veress	5-30 mmHg	30 mmHg
		Highflow	5-30 mmHg	30 mmHg
	Pediatric	Sensitive	5-15 mmHg	15 mmHg
		Vessel Harvesting	5-20 mmHg	20 mmHg
	Cardiac	Fieldflooding	Not applicable	

<b>The following settings can be made under indication settings:</b>				
User Menu	Indications	Modes	Options	Factory Settings
<b>Acoustic occlusion alert</b>	Standard	Veress	- ON - OFF	ON
		Highflow	- ON - OFF	ON
	Bariatric	Veress	- ON - OFF	ON
		Highflow	- ON - OFF	ON
	Pediatric	Sensitive	- ON - OFF	ON
	Cardiac	Vessel Harvesting	- ON - OFF	ON
		Fieldflooding	Not applicable	
<b>Timer for smoke evacuation</b>	Standard	Veress	1, 2, 3, 5, 10, 30 min, Off	5 min
		Highflow	1, 2, 3, 5, 10, 30 min, Off	5 min
	Bariatric	Veress	1, 2, 3, 5, 10, 30 min, Off	5 min
		Highflow	1, 2, 3, 5, 10, 30 min, Off	5 min
	Pediatric	Sensitive	1, 2, 3, 5, 10, 30 min, Off	5 min
	Cardiac	Vessel Harvesting	Not applicable	
		Fieldflooding	Not applicable	
<b>Foot control level (for smoke evacuation)</b>	Standard		Evacuation rate 6 and 12 l/min	OFF
<b>Setting the initial evacuation rate for each indication</b>	Bariatric		Evacuation rate 6 and 12 l/min	OFF
	Pediatric	Sensitive	Evacuation rate 6 and 12 l/min	OFF
	Cardiac	Vessel Harvesting	not available	
		Fieldflooding	not available	

**Under System, the following settings can be made:**

User Menu	Indications	Modes	Options	Factory Settings
<b>System settings</b>	All indications	All modes	<b>Brightness</b> 20 %, 40 %, 60 %, 80 %, 100 %	80 %
			<b>Volume</b> 20 %, 40 %, 60 %, 80 %, 100 %	80 %
			<b>Gas Supply</b> House Gas Bottled gas	House Gas
			<b>Language</b> - English - French - German - Spanish - Italian - Chinese - Swedish - Korean - Polish - Russian - Portuguese - Slovak - Czech	English
<b>Reset:</b>	All indications	All modes	This function allows you to reset to factory default all changes done in the user menu.	

**Under Service you will find the following information:**

User Menu	Indications	Modes	Options	Factory Settings
<b>Software</b>	All indications	All modes	Shows the installed software version.	
<b>Firmware</b>	All indications	All modes	Displays the installed firmware version	
<b>Next Maintenance</b>	All indications	All modes	Shows the next maintenance date DD. MM.YYYY	
<b>Technician</b>	All indications	All modes	Access to the service area	

**Automatic venting system****8 Safety Functions****8.1 Automatic Venting System**

The device is equipped with an automatic venting system.

When the insufflator detects an overpressure it automatically activates the venting system. The venting system releases gas from the cavity until the set nominal value has been reached again.

The automatic venting system can be configured individually in the user menu for each indication.

If the venting system is activated, the safety threshold for the release is 3 mmHg above the set pressure with a delay of 3 seconds.

Safety threshold and delay can be set in user menu.

If the venting valve is turned off, the message "**Venting valve deactivated!**" is displayed and you will hear an acoustic signal as soon as the indication is loaded.

Even if the valve is turned off, the venting function runs automatically and without delay, if the abdominal pressure is greater than 45 mmHg.

**Insufflation source**

The manufacturer advises against using additional, non-pressure controlled insufflation sources during minimally invasive surgical procedures.

The use of CO<sub>2</sub>-cooled lasers and argon beamers can lead to values exceeding the recommended and safe pressure rating.

**WARNING!****Additional insufflation sources/automatic venting system**

**Make sure the automatic venting system is activated (see chapter 7 Menu). The use of additional insufflation sources increases the intra-abdominal pressure. Continuously monitor intra-abdominal pressure over the course of the entire insufflation if additional sources are used.**

**WARNING!**

**The venting rate of the automatic venting system is limited. Always monitor the actual pressure when using additional insufflation sources.**

**Overpressure**

If the overpressure cannot be reduced by the automatic venting system within 5 seconds, the display depicts Overpressure. An acoustic signal is emitted.

Once the nominal pressure limit of the used indication has been reached/exceeded the display depicts the message **Maximum pressure reached!** A warning signal is emitted at the same time.

If the actual pressure exceeds 45 mmHg for more than 10 seconds the display depicts the message Excessive Overpressure along with an acoustic signal.

**8.2 Occlusion****"Occlusion" alert**

When tube, Veress needle, or trocar are blocked, the message "**Occlusion!**" is depicted and a warning signal is emitted. The actual pressure display depicts **0**.

The acoustic signal can be turned on/off in the user menu.

**CAUTION!**

**The pressure in the body cavity cannot be measured during an occlusion!**

### 8.3 Contamination

**WARNING!**

**Contamination**

**Do not use device and/or accessories if signs of contamination are detected. Make sure the device or/and accessories can no longer be operated until a qualified service technician conducts the appropriate tests and repairs.**

When fluid has penetrated the device via the insufflation tube connection, the display depicts **Device contaminated!** and gives an acoustic signal. It is possible to conclude any currently ongoing surgery with this device.

Error message "Device contaminated!"

If insufflation is stopped subsequently, it can be restarted within 20 minutes and the surgery can be continued. A change of indication is no longer possible. If insufflation is not restarted during 20 minutes or if the device is turned off and on again, further use of the device is not possible. This is to prevent cross-contamination.

When you turn on an already contaminated device, a corresponding message appears on the display. The device can no longer be used. The contaminated device has to be clearly marked as contaminated and sealed in two separate protective layers of safety foil. Make sure the device can not be used until a qualified service technician conducts the appropriate tests and repairs.

**Contamination**

### 8.4 Leakage Detection

If the device is not able to achieve the set pressure a leakage notification is displayed.

Leakages can occur at port sites, port sealings or due to excessive use of a suction device.

### 8.5 Gas Heating - Overtemperature

Use the gas heating and the heating tube to insufflate lukewarm gas (37 °C).

**CAUTION!**

**Do not subject the heating tube to direct heat (e.g., endoscope connected with light source) or high room temperatures.**

The status bar depicts the following in case of an overtemperature: Gas heating - overtemperature! A warning beep can be heard. Gas heating and insufflation are switched off.

**WARNING!**

**If the temperature sensor measures a gas temperature > 42 °C, unplug the heating tube plug from the device. Hot gas in the abdomen can lead to serious injuries.**

**Manufacturer's specifications****9 Care and Maintenance**

Special care is necessary when servicing, maintaining, and storing the device and its accessories to maintain the functionality of the device and its accessories.

**9.1 Cleaning the Device**

The manufacturer recommends that the device is cleaned after every procedure as follows:

1. Use the **ON/OFF** switch to turn the device off.
2. Remove the power cable.
3. Wipe the surface of the device with a soft cloth moistened with the surface disinfectant based on alcohol or aldehyde (for example Meliseptol® rapid). The concentration of the used disinfectant depends on the information provided by the manufacturer of the disinfectant. Make sure moisture does not enter the device.

**NOTE!**

**Do not sterilize the device.**

**Manufacturer's specifications****9.2 Maintenance Intervals**

The manufacturer stipulates that qualified personnel or hospital technicians must regularly test the device to assess its functionality and technical safety. This check must be carried out once a year. The tests are described in chapter 10 Annual Inspection.

Regular inspections will assist in early detection of possible malfunctions. This helps preserve the device and increases its safety and service life.

**NOTE!**

**Service or maintenance work may not be carried out during surgery.**

**Two-year maintenance interval****9.3 Maintenance Carried out by Authorized Service Technician**

An authorized service technician must inspect and service the device at appropriate intervals to ensure the safety and functionality of the device. The device must be serviced at least once every two years depending on frequency and duration of use. If the service interval is not maintained, the manufacturer does not assume any liability for the operational safety of the device.

The device displays the next maintenance due date during the automatic self-test.

A sticker located on the rear panel of the device lists the next service or maintenance date.

Authorized service technicians are trained and certified only by the manufacturer.

**Authorized trained personnel**

All of the service tasks, such as changes, modifications, repairs, calibrations, etc. may be carried out only by the manufacturer or manufacturer-approved trained and skilled technicians. The technical documentation (service instructions, values for calibration, etc.) may be provided by the manufacturer to authorized service technicians.

**Unauthorized personnel**

The manufacturer is not liable for the operational safety of the device if unauthorized persons conduct this maintenance or any other service tasks.

**Liability**

Unauthorized opening of the device and repairs performed by unauthorized personnel or third parties and/or changes or modifications release the manufacturer of any liability concerning the operational safety of the device.

Receiving technical documentation from the manufacturer does not authorize individuals to perform repairs, adjustments, or alterations on the device or accessories/peripherals.

#### Technical documents

Ask the service technician for a certificate after he or she has inspected the unit or performed any service tasks. This certificate lists the type and scope of the service as well as the date and name of the servicing company together with the signature of the service technician.

#### Certification

### 9.4 Replacing the Fuse

#### WARNING!

##### Replace the fuse

**Replace the fuse only with a fuse provided by the manufacturer (see chapter 14 Accessories).**

The fuse may be defective and is in need of replacement if:

- displays and LEDs do not light up,
- the device does not function.

Check to make sure

- the main power supply cable is properly connected to the power supply input and to a safety socket,
- the fuse of the building's electrical supply is functioning.


#### WARNING!

**Unplug the power cable from the device before checking the fuse.**

The device does **not** have to be opened to replace the fuse.

1. Switch device off.
2. Disconnect device from power supply. To do so, pull the mains plug from the wall socket.
3. Remove power connection cable from mains socket.
4. The fuse holder is located right next to the mains socket.
5. Remove fuse holder as depicted in Fig. 9-1 Opening the fuse holder.
6. **A** Undo the latch of the fuse holder with your fingers.
7. **B** Remove the fuse holder.
8. **C** Check fuse.
9. Insert a new fuse. Use only the specified type of fuse (see chapter 13 Technical Data).
10. Insert the fuse holder until it can be heard snapping into place.
11. Reconnect device with the mains power. To do so, plug the mains plug into the wall socket.

**Fig. 9-1 Opening the fuse holder**

## **9.5 Care and Maintenance of Reusable Tube Sets (NOT APPROVED FOR SALE IN THE UNITED STATES)**

The following instructions apply when using reusable tube sets.

### **General rules**

The reusable tube systems must be cleaned, disinfected and sterilized before each use. This applies in particular also to the first use after delivery, because all reusable tube systems are supplied non-sterile (cleaning and disinfecting after removing transport protection packing, sterilization after removing packaging). Effective cleaning and disinfection is an essential prerequisite for effective sterilization. As part of your responsibility for the sterility of products during use, please make sure that

- only cleaning/disinfection and sterilization methods sufficiently validated for the device and specifically designed for the product are used.
- the used devices (RDG, sterilizer) are regularly maintained and checked.
- adherence to the validated parameters during each cycle is ensured.

Please also comply with the legal or statutory requirements applicable to your country as well as the hygiene rules of the physician's office or hospital. This applies in particular to the differing requirements regarding effective prion inactivation.

### **WARNING!**

**The reusable tube system consists of silicone, polysulfone (PSU). For cleaning and disinfection, use only pH-neutral or slightly alkaline cleaners (e.g., neodic MediClean 2.0 %), disinfectants (e.g., Lysetol V 8%), drying, and clear rinsing substances approved for the utilized materials.**

**If using unsuitable substances (e.g., the neodic clear rinsing agent MediKlar), the tube system and especially the PSU connection points may become damaged.**

### **Recycling limitations**

The manufacturer has verified the reusable tube system for a limited number of recycling cycles. Always comply with the instructions on the label.

Never exceed the number of uses indicated by the manufacturer.

**WARNING!**

**Use the tear-off tabs attached to the tube to keep track of the number of sterilization cycles. The heating tube may not be sterilized after the last tab has been removed.**

**Please check the reusable tube for signs of damage after sterilization and before use. Never use a tube showing signs of damage, especially brittleness and perforation.**

### 9.5.1 Cleaning and Disinfecting

A mechanical process (cleaning and disinfecting equipment) should be used if possible. A manual process - even when using an ultrasonic bath - should be used only if a mechanical method is not available or not possible due to the significantly reduced efficiency and reproducibility. Pretreatment is required in both cases.

Directly after use (within maximum 2 h) large impurities must be removed from the reusable tube system:

**Basics****Pretreatment****Process**

- **Rinse the products at least 1 min under running water (temperature <35 °C / 95 °F). Rinse all lumina of the products three times using a disposable syringe (minimum volume 100 ml).**
- Place the products in a freshly prepared cleaning solution<sup>1</sup>. Observe relevant requirements as listed in Manual cleaning and disinfecting concerning the selection, concentration, temperature, and exposure time. Use only freshly prepared solutions. Rinse all lumina of the products three times using a disposable syringe (minimum volume 100 ml).
- Manually remove all externally visible particles with a clean soft cloth and soft brush, which you are using only for this purpose.
- Rinse the products at least 1 min under running water (temperature <35 °C / 95 °F). Rinse all lumina of the products three times using a disposable syringe (minimum volume 100 ml).

<sup>1</sup> If your are using a combination cleaner and disinfectant for this step - e.g., for reasons of occupational safety and health - please keep in mind that this product should be free of aldehyde (otherwise blood particles may become permanently deposited), should have a certified efficacy (e.g., VAH/DGHM or FDA/EPA approval/clearance/registration or CE marking), is suitable for the disinfection of products, and is compatible with the products being cleaned and disinfected (see "material resistance"). Please keep in mind that the disinfectant used for the pre-treatment serves only personal protection and cannot replace the subsequent disinfection carried out after the cleaning step.

When selecting a cleaning and disinfecting equipment, make sure

**Mechanical cleaning/disinfecting (cleaning and disinfecting equipment)**

- this device has a certified efficacy (e.g., DGHM or FDA approval/clearance/registration or CE marking according to DIN EN ISO 15883).
- to use a certified program for thermal disinfection (AO value > 3000 or - for older devices - at least 5 min at 90 °C/194 °F), if possible (with chemical disinfection risk of disinfectant residues on the products).
- the used program for the products is suitable and contains sufficient rinsing cycles.
- rinsing is done only with water (e.g., purified water/highly purified water) that is sterile or low-germ (max. 10 germs/ml) and low in endotoxins (max. 0.25 endotoxin units/ml).
- the air used for drying is filtered (oil-free, low in particles and germs).
- the cleaning and disinfecting equipment will be regularly maintained and checked.

When choosing the cleaning agent system, make sure

- it is generally suitable for the cleaning of instruments made of metal and plastics.
- an additional, suitable disinfectant with certified efficacy (e.g., VAH/DGHM or FDA/EPA approval/clearance/registration or CE marking) is used and that this disinfectant is compatible with the utilized cleaning agent - unless a thermal disinfection process is applied.
- the chemicals used are compatible with the products (see "material resistance").

It is imperative to comply with the concentrations, temperatures, and exposure times, as well as the post-rinsing requirements specified by the manufacturer of the cleaning agent and/or disinfectant.

#### **Process**

1. Place the products into the cleaning and disinfecting equipment. Make sure the products do not touch one another. Close all lumina of the products using the existing Luer lock adapter at the rinsing port of the cleaning and disinfecting equipment.
2. Start the program.
3. Remove the products from the cleaning and disinfecting equipment after the end of the program.
4. Check and pack the products as soon as possible after removal (see "Inspection," "Maintenance" and "Packaging," possibly after additional drying in a clean location).

Documentation of the fundamental suitability of the reusable heating tube system for an effective mechanical cleaning and disinfection was prepared by an independent accredited and recognized (§15 (5) MPG) testing laboratory using the cleaning and disinfecting equipment type G 7836 CD (thermal disinfection, Miele & Cie. GmbH & Co., Gütersloh) and the pre-cleaning/cleaning agent neodisher MediZym (Dr. Weigert GmbH & Co. KG, Hamburg). This was done by taking into account the procedure described above.

#### **Manual cleaning and disinfecting**

When choosing the cleaning and disinfecting agent to be used, make sure

- they are generally suitable for the cleaning or disinfection of tube systems made of metal and plastics.
- a disinfectant with certified efficacy (e.g., VAH/DGHM or FDA/EPA approval/clearance/registration or CE marking) is used and that this disinfectant is compatible with the utilized cleaning agent.
- the chemicals used are compatible with the products (see Material resistance).

Combined cleaning/disinfecting agents should not be used if possible. Combined cleaning/disinfecting agents can be used only in cases of very low contamination (no visible impurities).

It is imperative to comply with the concentrations, temperatures, and exposure times, as well as the post-rinsing requirements specified by the manufacturer of the cleaning agent and disinfectant. Use only freshly prepared solution, only water (e.g., purified water/highly purified water) that is sterile or low-germ (max. 10 germs/ml) and low in endotoxins (max. 0.25 endotoxin units/ml) and only a soft, clean and lint-free cloth and/or filtered air for drying.

#### **Cleaning procedure**

1. Insert the products into the cleaning solution for the prescribed contact time and make sure the products are sufficiently submerged. Make sure the products do not touch one another.
2. Rinse all lumina of the products at least five times at the beginning or at the end of the contact time using a disposable syringe (minimum volume 100 ml). Clean the exterior with a clean, soft cloth and a soft brush that you use just for this purpose.
3. Remove the products from the cleaning solution and rinse at least three

- times thoroughly (for at least 1 min) with water. Rinse all lumina of the products at least five times using a disposable syringe (minimum volume 100 ml).
4. Check the products (see Inspection and Maintenance).

#### Disinfecting procedure

1. Insert the cleaned and checked products into the disinfecting solution for the prescribed contact time and make sure the products are sufficiently submerged. Make sure the products do not touch one another. Rinse all lumina of the products at least five times at the beginning or at the end of the contact time using a disposable syringe (minimum volume 100 ml).
2. Remove the products from the disinfecting solution and rinse at least five times thoroughly (at least 1 min) with water. Rinse all lumina of the products at least five times using a disposable syringe (minimum volume 100 ml).
3. Dry the products by blowing off/blowing out with filtered compressed air.
4. Package the products as soon as possible after removal (see Packaging, possibly after additional drying in a clean location).

Documentation of the fundamental suitability of the reusable heating tube system for an effective manual cleaning and disinfection was prepared by an independent accredited and recognized (§15 (5) MPG) testing laboratory using the pre-cleaning/cleaning agent neodisher Medizym and the disinfecting agent neodisher Septo 3000 (Dr. Weigert GmbH & Co. KG, Hamburg). This was done by taking into account the procedure described above.

Check the reusable tube system after cleaning or cleaning/disinfection for corrosion, damaged surfaces, chipping, soiling and discoloration and discard damaged products. (See label of the product for the numerical restriction of reuse). In case of damage (scratches, flaking) on the chrome-plated Luer-lock adapter, the products must be replaced immediately (no reuse permitted). Products still contaminated or dirty must be once again cleaned and disinfected.

Instrument oils or grease may not be used.

#### Inspection

#### Maintenance

Please pack the tube systems into sterile disposable packaging (single or double packaging) that comply with the following requirements (material/process):

#### Packaging

- DIN EN ISO/ANSI AAMI ISO 11607
- Suitable for steam sterilization (temperature resistance up to min. 142 °C (288 °F) with sufficient steam permeability)
- Sufficient protection of the products or sterilization packaging against mechanical damage

#### 9.5.2 Sterilization of Reusable Tube Sets (NOT APPROVED FOR SALE IN THE U.S.)

Sterilization must be carried out with only the following sterilization method listed below; any other sterilization procedures are not permitted.

- Fractionated vacuum method<sup>2,3</sup>(with sufficient product drying<sup>4</sup>)
- Steam sterilizer according to DIN EN 13060/DIN EN 285 or ANSI AAMI ST79
- validated according to DIN EN ISO 17665 (valid IQ/OQ (picking) and product-specific performance assessment (PQ))
- maximum sterilization temperature 138 °C (280 °F; plus tolerance according to DIN EN ISO 17665)
- Sterilization time (exposure time at sterilization temperature) at least 5 minutes at 132 °C (270 °F) / 134 °C (273 °F)

#### Steam sterilization

<sup>2</sup>at least three vacuum stages

<sup>3</sup>Use of the less effective gravitation procedure is allowed only in case of non-availability of the fractionated vacuum procedure; requires significantly longer sterilization times and must be validated specifically for the products, devices, and parameters in sole responsibility of the user.

<sup>4</sup>The actual required drying time depends on parameters, which are the sole responsibility of the user (loading configuration and density, sterilizer condition, etc.) and must therefore be determined by the user. However, drying times should never be less than 10 minutes.

Flash sterilization procedures are not permitted. Furthermore, do not use hot air sterilization, radiation sterilization, formaldehyde or ethylene oxide sterilization, or plasma sterilization.

*Documentation of the general suitability of the reusable heating tube system for effective sterilization was provided by an independent accredited and recognized (§ 15 (5) MPG) testing laboratory using the autoclave Systec HX-320.*

This was done by taking into account the procedure described above.

#### **Storage**

After sterilization, the products in the sterilization packaging must be stored dry and dust-free.

#### **Material resistance**

When selecting the cleaning and disinfecting agents, please ensure that the following ingredients are not included:

- organic, mineral and oxidizing acids (minimum admissible pH 6.5)
- strong alkalies (maximum admissible pH 9.5, enzymatic/neutral or slightly alkaline cleaner recommended)
- organic solvents (alcohols, ethers, ketones, petroleum)
- oxidizing agents (such as hydrogen peroxide)
- halogens (chlorine, iodine, bromine)
- aromatic/halogenated hydrocarbons

Never clean the reusable tube with metal brushes or steel wool. Rinse aids or acid neutralizers are not permitted. All products may not be exposed to temperatures in excess of 142 °C (288 °F)! Failure to comply with these instructions voids any liability.

## 10 Annual Inspection

Each test conducted must be documented with date and signature on the test log. **Measured values and tolerances**

The following measuring tools and resources were used by the manufacturer to determine the listed measurements and tolerances:

- (1) Manometer Range 0-100 mmHg, error class 1.6
- (2) Syringe 60 ml
- (3) PG012 (Disposable tube set)
- (4) Silicone tube 4 mm x 40 cm
- (5) T adapter 6/4/6
- (6) Veress cannula Length 100 mm  
opening diameter 1.4 mm,  
inner cannula diameter 1.6 mm

An authorized service technician must check the device if the specified parameters and tolerances are exceeded.

### 10.1 Electrical Safety Test

1. Perform a visual inspection. Make sure that
  - the fuse corresponds with the specifications indicated by the manufacturer,
  - labels and stickers on device are legible,
  - the mechanical condition of the device allows for its safe use,
  - the device is clean to ensure proper and safe functionality.
2. Carry out the measurements for the ground leakage current, short-circuit current/housing leakage current, and the protective conductor resistance as per IEC 62353 in the current version or according to the applicable national standard.

### 10.2 Basic Function Test (in Standard Indication)

1. Remove insufflation tube from device.
2. Use the On/Off switch to turn the device on. The device now conducts a device check. A acoustic signal can be heard. Set the Standard operating mode.
3. A nominal pressure of 15 mmHg and a nominal flow of 3 l/min are set at the factory.
4. The following values are displayed:

Nominal pressure 15\* (mmHg)

Nominal gas flow 3\* (l/min)

Actual pressure 0 (mmHg)

Gas consumption 0.0 (l)

\*These values correspond with the factory setting. If values in the user menu were changed, these changed values are displayed. Insert an PG012 (Disposable tube set).

5. Start insufflation: Press the **START** button. The following values are displayed:

Actual pressure 0 (mmHg)

**Veress Insufflation** is displayed. Streaming gas can be heard at the insufflation tube.

6. Select a nominal gas flow of 20 l/min. The following value is displayed:

Nominal gas flow 20 l/min

Actual pressure 0 (mmHg)

**Insufflation Active** is displayed. Streaming gas can be heard at the insufflation tube.


7. Stop insufflation: Press the **STOP** button.  
The following values are displayed:  
Actual pressure 0.0 (mmHg)  
Gas consumption >0.0 (l)  
8. Press the gas volume button.  
Gas consumption 0.0 (l)

The basic function check of the device is complete.

### 10.3 Pressure Sensor Test (in Standard Indication and Veress Mode)

**Fig. 10-1 Tools for pressure sensor test**

- ① Manometer
- ② Syringe
- ③ PG012 (disposable tube set)
- ④ Silicone tube
- ⑤ T adapter

1. Set the Standard indication.
2. Select a nominal gas flow rate of 1.0 l/min.  
Do not press the **START/STOP** button.

#### CAUTION!

**Never use the syringe to extract gas from the device.**

3. Connect a manometer and an air-filled syringe to the insufflation tube connection.
4. Use the syringe to generate a pressure of at least 10 mmHg, which is indicated on the manometer.  
Actual pressure display:  $10 \pm 2$  (mmHg)
5. Use the syringe to generate a pressure of at least 20 mmHg, which is indicated on the manometer.  
Actual pressure display:  $20 \pm 2$  (mmHg)
6. Use the syringe to generate a pressure of at least 30 mmHg, which is indicated on the manometer.  
Actual pressure display:  $30 \pm 2$  (mmHg)

### 10.4 Pressure Monitoring Test (in Standard Indication and Veress Mode)

1. Set the Standard indication.
2. Select a nominal pressure of 15 mmHg and a nominal gas flow of 3 l/min.
3. Activate the **Venting Valve** in the user menu.
4. Use the syringe to generate a pressure of at least 27 mmHg and maximum 30 mmHg, indicated on the manometer. Start insufflation: Press the **START** button. An acoustic signal is emitted (for 5 s) and the display depicts **Overpressure** if the pressure exceeds 19 mmHg.
5. Reduce the pressure. The warning ends when the pressure falls below 19 mmHg (nominal pressure plus 4 mmHg). Stop insufflation: Press the **STOP** button.
6. Select a nominal pressure of 29 mmHg.
7. Use the syringe to generate a pressure of at least 30 mmHg, which is indicated on the manometer. Start insufflation: Press the **START** button. An acoustic signal is emitted with a pressure of more than 30 mmHg and the display de-

picts **Maximum Insufflation Pressure Reached**.


8. Reduce the pressure. The warning ends when the pressure falls below 31 mmHg. Stop insufflation: Press the **STOP** button.

## 10.5 Venting Valve Test

See also 10.3 Pressure Sensor Test (in Standard Indication and Veress Mode).

1. In the user menu under Mode Settings/Insufflation, adjust the venting system so that at an overpressure of 3 mmHg the venting will open after three seconds (see chapter 7.2.1 Insufflation).
2. Select a nominal pressure of 15 mm Hg and a nominal gas flow rate of 10 l/min.
3. Use the syringe to generate a pressure of at least 18 mmHg, which is indicated on the manometer. Start insufflation. The venting valve is activated with a pressure of more than 18 mmHg (for 3 s).

## 10.6 Device Maximum Pressure Test

1. Select Standard indication.
2. Select the maximum nominal gas flow.
3. Connect a manometer (1) and an open Veress cannula (3) to the insufflation tube connection (2).
4. Select a gas flow of > 5 l/min.
5. Start insufflation:  
Press the **START** button. The manometer registers a pulsing pressure increase. When the pressure stabilizes, the manometer registers a pressure between 55 and 65 mmHg.
6. Stop insufflation:  
Press the **STOP** button.

### **10.7 Gas Flow Rate Test**

- Connect a tube set.
- Select a nominal gas flow rate of 15 l/min.
- Start insufflation:  
    Press the **START** button.
- Press the gas volume button (0.0 l must be displayed).  
    Now start measuring for one minute.
- Stop the insufflation after one minute: Press the **STOP** button.  
    The gas consumption should be at least 11 l.

Each successfully conducted test must be documented with the test log.

## 11 Electromagnetic Compatibility

Medical devices are subject to special safety and protective measures concerning electromagnetic compatibility (hereafter abbreviated as EMC).

### Precautionary measures

This device is to be used only for the purposes described in the manual and has to be installed, set up, and operated in compliance with the EMC notes and instructions.

#### CAUTION!

##### Accessories

**To ensure compliance with the requirements of IEC 60601-1-2 in the current version, the device Flow 50 must be used only with the accessories listed in Chapter 14.**

To ensure the basic safety and essential functionality in relation to electromagnetic interference over the life of the device, the device must be restarted after 24 hours so that a diagnostic self-test can be performed. The maintenance intervals indicated in Chapter 9.2 Maintenance Intervals must also be observed.

#### CAUTION!

##### Continuous operation

**After 24 hours of continuous operation, a device self-test must be carried out.**

- Switch device off and on again.

### 11.1 Impact of Mobile and Portable HF Communication Devices

The emission of high frequency energy by mobile communication devices may impact the function of the electrical medical device. Operating such devices (e.g., cell phones, GSM phones) in the proximity of the electrical medical device is prohibited.

### 11.2 Electrical Connections

The following are ESD precautionary measures:

### ESD (Electrostatic Discharge) precautionary measures

- Apply potential equalization (PE), if available on your equipment, to all devices to be connected.
- Use only the listed equipment and accessories.

Hospital employees should be informed about and trained in ESD precautionary measures.

### 11.3 Guidance and Manufacturer's Declaration - Electromagnetic Immunity

The Flow 50 is intended for use in the electromagnetic environment as described below. The user/operator of the insufflator should make sure the device is operated within such an environment.

IMMUNITY tests	IEC 60601 Test level	Compliance level	Electromagnetic environment - guidance
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	3 V	<p>Portable and mobile RF communications equipment should be used no closer to any part of the insufflator, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended safety distance:</p> <p><b>Recommended separation distance</b></p> $d = 1.2\sqrt{P} \text{ for } 150 \text{ KHz to } 80 \text{ MHz}$ $d = 1.2\sqrt{P} \text{ for } 80 \text{ MHz to } 800 \text{ MHz}$ $d = 2.3\sqrt{P} \text{ for } 800 \text{ MHz to } 2.7 \text{ GHz}$ <p>where <math>P</math> is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and <math>d</math> is the recommended separation distance in meters (m).</p>
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.7 GHz  80% AM by 1 KHz	3 V/m	<p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey<sup>a</sup> should be less than the compliance level in each frequency range.<sup>b</sup></p> <p>Interference may occur in the vicinity of equipment marked with the following symbol:</p> 

Note 1: At 80 MHz and 800 MHz, the higher frequency range applies.

Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

<sup>a</sup> Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the insufflator is used exceeds the applicable RF compliance level above, the insufflator should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the insufflator.

<sup>b</sup> Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

**IMMUNITY TEST****Radiated RF IEC 61000-4-3**

<b>Test frequency (MHz)</b>	<b>Band (MHz)</b>	<b>Service</b>	<b>Modulation</b>	<b>Maximum power (W)</b>	<b>Distance (m)</b>	<b>IMMUNITY TEST LEVEL (V/m)</b>
385	380-390	TETRA 400	Pulse modulation 18 Hz	1.8	0.3	27
450	430-470	GMRS 460, FRS 460	FM $\pm 5$ kHz deviation	2	0.3	28
710						
745	704-787	LTE Band 13, 17	Pulse modulation 217 Hz	0.2	0.3	9
780						
810						
870	800-960	GSM 800/900, TETRA 800, IDEN 820, CDMA 850, LTE Band 5	Pulse modulation 18 Hz	2	0.3	28
930						
1720						
1845	1700-1990	GSM 1800; CDMA 1900; GSM 1900; DECT; LTE Band 1,3, 4, 25; UMTS	Pulse modulation 217 Hz	2	0.3	28
1970						
2450	2400-2570	Bluetooth, WLAN, 802.11 b/g/n RFID 2450, LTE Band 7	Pulse modulation 217 Hz	2	0.3	28
5240						
5500	5100 5800	WLAN 802.11 a/n	Pulse modulation 217 Hz	0.2	0.3	9
5785						

**11.4 Guidance and Manufacturer's Declaration - Electromagnetic Emissions**

The Flow 50 is intended for use in an environment as described below. The user/operator of the insufflator should make sure the device is operated within such an environment.

<b>Emissions test</b>	<b>Compliance</b>	<b>Electromagnetic environment - guidance</b>
RF emissions CISPR 11	Group 1	The Flow 50 uses RF energy solely for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The Flow 50 is suitable for use in all establishments, other than domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/flicker emissions IEC 61000-3-3	Complies	

## 11.5 Recommended Separation Distances between Portable and Mobile RF Communications Equipment and the Flow 50

The Flow 50 is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The user/operator of the insufflator can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the insufflator as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter (W)	Separation distance according to frequency of transmitter (m)		
	<b>150 kHz to 80 MHz outside ISM bands</b> $d = 1.2\sqrt{P}$	<b>80 MHz to 800 MHz</b> $d = 1.2\sqrt{P}$	<b>800 MHz to 2.7 GHz</b> $d = 2.3\sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance  $d$  in metres (m) can be determined using the equation applicable to the frequency of the transmitter, where  $P$  is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

Note 1: At 80 MHz and 800 MHz, the higher frequency range applies.

Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

### WARNING!

Portable HF communication devices can have an effect on the performance characteristics of the device. Therefore such devices must be kept a minimum distance of 30 cm (independent of any calculation) from the device Flow50 and its accessories and cables

## 12 Error and Warning Messages

Error and warning messages	Cause	Troubleshooting
<b>Electronic error! Restart the device. If the error occurs again, call service.</b>	Malfunctions in the electronic system.	Restart the device. If the error occurs again, contact technical service.
<b>Calibration error! The device must be re-calibrated. Call service.</b>	The device is not calibrated properly.	The device must be re-calibrated. Contact technical service.
<b>Device temperature error! Do not use device. Turn device off!</b>	Temperature in device is too low or too high.	Serious fault. Do not use device. Contact technical service.  Make sure the device has been acclimatized to its environment before turning on.
<b>Device contaminated! Fluid has penetrated the device! Call service!</b>	The device is contaminated with fluid.	The device must be checked by an authorized service technician. It is to be clearly marked and wrapped twice in safety foil before sending to the service department for repairs.
<b>Device service overdue! Flow rate and pressure might not be accurate anymore.</b>	Last maintenance is more than 24 months ago	The device must be serviced by an authorized service technician.
<b>Device service due!</b>	Maintenance is due in 4 weeks	Continued use of the device is possible. However, the device must be serviced by an authorized service technician soon.
<b>Vacuum pump error! Restart the device. If the error occurs again, call service! Insufflation is still possible but smoke evacuation function will not work.</b>	Suction pump defective	Contact the technical service to check the smoke evacuation. The device may still be used for the insufflation. The smoke evacuation function is deactivated.
<b>Excessive overpressure!</b>	The actual pressure > 45 mmHg	Determine the cause for exceeding the nominal pressure. If the venting valve is disabled, the overpressure must be reduced manually by opening the trocar.
<b>Maximum insufflation pressure reached!</b>	Nominal pressure limit 20/30 mmHg has been reached/exceeded.	Reduce the abdominal pressure by reducing the set pressure or opening a stopcock.
<b>Overpressure!</b>	Abdominal pressure exceeds set pressure and can not be reduced by the automatic venting system within 5 seconds.	Determine the cause for exceeding the set pressure. If the venting valve is disabled, the overpressure must be reduced manually by opening the trocar.
<b>Gas heating - Over temperature!</b>	The temperature of the gas exceeds 42 °C / 108 °F. The device terminates gas heating and insufflation.	There are three options to continue insufflation: <ul style="list-style-type: none"><li>• Wait until the temperature has decreased below 43 °C, insufflation and heating will resume automatically.</li><li>• Disconnect and reconnect tube set if temperature is below 43 °C insufflation with heating can be resumed (manual start).</li><li>• Exchange tube set with a new one.</li></ul> If this error message is displayed again, have the gas heating checked by a qualified service technician.
<b>Safety threshold</b>	Nominal pressure > 15 mmHg (Standard, Bariatric).  Nominal pressure > 12 mmHg (Pediatric, Vessel Harvesting).  Actual flow > 5 l/min (Pediatric).  Actual flow > 6 l/min (Vessel Harvesting).	Confirm the query to further increase the pressure. Exceeding this safety limit is to be decided by and the responsibility of the user/operator.
<b>Change gas bottle!</b>	Low insufflation gas supply (bottle or central gas supply).	Check gas supply.
<b>Gas level low. Prepare to change gas bottle.</b>	Gas supply < 15 bar/217.5 psi.	Prepare for changing the gas bottle.
<b>Check gas supply!</b>	Low insufflation gas supply (bottle or central gas supply) during insufflation.	Insufflation starts automatically after restoring the gas supply.

Error and warning messages	Cause	Troubleshooting
<b>Heating tube set not valid</b>	The tube set is not valid.	Insert a valid tube set.
<b>Gas heating defective!</b>	Gas heater malfunction.	Insufflation can be continued in this case but the heating function is not available.  Restart the device if you want to use the heating function.  If the gas heating is still defective, have an authorized service technician check or repair the device.
<b>Venting valve deactivated</b>	The venting system has been deactivated	The user can activate the venting valve in the user menu.
<b>Check for leakage!</b>	The device cannot provide the predetermined pressure.	Check for leaks at the connection points and at the connection seals, or where a suction instrument has been used excessively.
<b>Smoke evacuation paused</b>	Smoke evacuation tube is blocked.	Make sure that the valve on the trocar is open and that there are no kinks in the smoke evacuation tube.
<b>Suction line not open or clogged</b>		
<b>Occlusion!</b>	Tube, Veress cannula or trocar is blocked.	Make sure that the Veress cannula or the trocar is positioned correctly and that there are no kinks in the insufflation tube.

### 13 Technical Data

Type or model designation:	PG150 (Flow 50)	
Software version:	See user menu/ Service / Software	
Mains voltage range:	AC 100-240 V ~	
Supply frequency range:	50/60 Hz	
Fuse designation:	2x T 4 A H, 250 V, UL recognized	
Power consumption:	Current (A)	Power consumption (VA)
Normal operation 100V AC/60 Hz	1.05	105
Normal operation 240V AC/50 Hz	0.50	120
Maximum Power consumption	-	120
Protection class (I, II, III):	I	
Application part type (B, BF, CF):	CF	
Protection type (IP code):	IP21	
Classification (I, IIa, IIb, III):	IIa	
Operating conditions:	15-37 °C/59 - 98.6°F 20 - 90 % rel. humidity, non condensing 700 to 1060 hPa air pressure 3000 m/9843 ft maximum altitude above sea level for device use	
Storage and transportation conditions:	-20 to +60 °C/-4 to +140 °F 10-95% rel. humidity, non condensing 500 to 1060 hPa air pressure	
Maximum device noise emission	50 dB	
Heat radiation	≤ 43 °C	
Inlet pressure range	Maximum inlet pressure	80 bar/1160.3 psi
	Minimum inlet pressure for central gas	2.7 bar/39.2 psi
Maximum output pressure	75 mmHg (1 mmHg = 1.33 mbar = 133 Pa)	
Maximum gas flow	50 l/min	
Max. suction capacity	12 l/min	
Adjustable pressure range	1-30 mmHg	
Pressure display resolution	1 mmHg	
Dimensions	Width x Height x Depth	
Weight (kg)	330 mm x 200 mm x 400 mm	
Weight (kg)	12 kg	
Interfaces/ports/connections	Memory stick interface (USB Typ A) Service interface (Mini USB) Mains power connection (IEC-60320-1 C14)	
The device complies with the following standards	IEC 60601-1:1988 + A1:1991 + A2:1995 (Edition 2.0) IEC 60601-1:2005 + CORR. 1:2006 + CORR. 2:2007 + A1:2012 (Ed. 3.1) EN 60601-1:2006/A1:2013 (Edition 3.1) EN 60601-1-2:2007/AC:2010 (Edition 3) IEC 60601-1-2:2014 (Edition 4) IEC 60601-1-6:2010 (Edition 3) + A1:2013 HF emissions according to CISPR 11 Class B	

## 14 Accessories

### **WARNING!**

**Original accessories**

**For your own safety and that of your patient, use only original accessories.**

### Article

#### **Tube sets**

PG012	Disposable tube set incl. CO <sub>2</sub> filter, sterile
PG014	Reusable tube set for insufflation
PG096SU	Disposable heating tube set for insufflation, sterile
PG097	Reusable heating tube set for insufflation
PG106	Smoke Evacuation Filter Cassette
PG107SU	Disposable tube set for Smoke Evacuation, sterile
PG019	Disposable CO <sub>2</sub> filter, sterile

#### **Hoses for gas supply**

PG056	Hose for CO <sub>2</sub> gas bottle supply DIN:US, 1.5 m
PG057	Hose for CO <sub>2</sub> gas bottle supply ISO:US, 1.5 m
PG058	Hose for CO <sub>2</sub> gas bottle supply PIN Index:US, 1.5 m
PG061	Hose for CO <sub>2</sub> central gas supply DIN:US, 3 m
PG091	Hose for CO <sub>2</sub> central gas supply DIN:US, 5 m
PG092	Hose for CO <sub>2</sub> central gas supply NF:US, 5 m
PG093	Hose for CO <sub>2</sub> central gas supply UNI:US, 5 m
PG094	Hose for CO <sub>2</sub> central gas supply AGA:US, 5 m
PG065	Hose for CO <sub>2</sub> central gas supply DISS:DISS, 3.0 m

#### **Miscellaneous**

PG066	DISS adapter for central gas supply
PG067	NIST adapter for central gas supply
PG068	Switching valve
PG150200	Sinter filter
PG090	Gas connector CO <sub>2</sub> angled
PG014200	ISO adapter
PG124	One-pedal footswitch
TA022006	Fuse 5x20mm 4,0 A H T 250 V
TA014996	Instructions for use

## 15 Glossary

Term	Meaning
ME Device	Medical electrical device used for therapy, monitoring or diagnosis of patients, equipped with no more than one connection to a supply network and which necessarily comes into physical or electrical contact with the patient or which transfers energy to or from the patient or which records or captures such energy transfer to or from the patient.
ME system	Medical electrical system consisting of a combination of devices, of which at least one is classified as a medical electrical device and specified by the manufacturer as such and which are joined together by a functional connection or by using a power strip.
EMC	The electromagnetic compatibility describes the capability of an electrical device to function satisfactorily in an electromagnetic environment without affecting the environment in any manner by electromagnetic interference, which can be problematic for other nearby devices or systems.

EN

## **16 Appendix**

### **16.1 Test Log**

<b>Date</b>	<b>Result</b>	<b>Comment</b>	<b>Signature</b>

**Index****"**

"Occlusion" alert 50

**A**

Authorized service technician 4  
Authorized trained personnel 52  
Automatic venting system 50

**C**

Care and maintenance 4  
Central gas supply 22  
Certification 53  
Connecting the gas heating 25  
Contamination 4, 51  
Contraindications 5

**D**

Delivery inspection 16  
Disconnect device from power supply 23  
Disposal 4

**E**

Error message "Contamination" 51  
ESD (Electrostatic Discharge) precautionary measures 63  
Exclusion of liability 4

**F**

Federal Law 4

**G**

Gas heating error 26  
Gas heating ON/OFF 26  
Gas supply with gas bottle 22  
Gas temperature exceeds 42°C 26  
Grounding contact 18

**I**

Inspection 57  
Insufflation source 50

**L**

Liability 52

**M**

Mains connection 18  
Maintenance 57  
Manual cleaning and disinfecting 56  
Manufacturer's specifications 52  
Material resistance 58  
Measured values and tolerances 59

**O**

Only for U.S. operators 18

**P**

Packaging 57  
Patients 40  
Potential equalization 18  
Precautionary measures 63  
Preparation 31

**R**

Resetting the display 32  
Returning the device 16

**S**

Setting up 17  
Subject to technical changes 3

**T**

Technical documents 53  
Two-year maintenance interval 52

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