**en Instructions for use/Technical description**

2D camera heads PV481, PV482, PV485

**USA Note for U.S. users**

This Instructions for Use is NOT intended for United States users. Please discard. The Instructions for Use for United States users can be obtained by visiting our website at [www.aesculapusaifus.com](http://www.aesculapusaifus.com). If you wish to obtain a paper copy of the Instructions for Use, you may request one by contacting your local Aesculap representative or Aesculap's customer service at 1-800-282-9000. A paper copy will be provided to you upon request at no additional cost.

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

## 2D camera heads PV481, PV482, PV485

### Legend

- 1 Safety lock
- 2 Endocoupler latch
- 3 Focus ring
- 4 Zoom ring
- 5 Camera head buttons
- 6 Lock

### Contents

1.	About this document.....	2
1.1	Scope .....	2
1.2	Safety messages.....	3
2.	Clinical use.....	3
2.1	Product description.....	3
2.1.1	Operating principle .....	3
2.1.2	Scope of Supply .....	4
2.1.3	Components required for use .....	4
2.2	Areas of use and limitations of use .....	4
2.2.1	Intended use.....	4
2.2.2	Endoscopic indications .....	4
2.2.3	Contraindications.....	4
2.3	Safety information.....	4
2.3.1	Clinical user .....	4
2.3.2	Product.....	5
2.3.3	Sterility.....	6
2.4	Preparation.....	6
2.4.1	Combination with medical electrical equipment .....	6
2.4.2	Inspections.....	6
2.5	First use .....	7
2.6	Application.....	7
2.6.1	General notes on use.....	7
2.6.2	System set-up .....	7
2.6.3	Function checks.....	7
2.6.4	Safe operation .....	8
2.7	Troubleshooting .....	9
2.7.1	Repairs .....	9
3.	Reprocessing procedure .....	10
3.1	General safety instructions.....	10
3.2	General information .....	10
3.3	Dismantling prior to carrying out the reprocessing procedure .....	10
3.4	Preparations at the place of use.....	10
3.5	Preparation before cleaning .....	10
3.6	Cleaning/disinfection.....	10
3.6.1	Product-specific safety instructions for the reprocessing procedure .....	11
3.6.2	Validated cleaning and disinfection procedure.....	11
3.7	Manual cleaning and disinfecting .....	12
3.7.1	Post use and pre-cleaning activities .....	12
3.7.2	Manual cleaning with immersion disinfection .....	12
3.8	Mechanical cleaning/disinfection.....	14
3.8.1	Mechanical alkaline cleaning and thermal disinfection.....	14

3.9	Mechanical cleaning/disinfecting with manual pre-cleaning.....	15
3.9.1	Post use and pre-cleaning activities.....	15
3.9.2	Manual pre-cleaning .....	15
3.9.3	Mechanical alkaline cleaning and thermal disinfecting .....	16
3.10	STERRAD Sterilization.....	16
3.10.1	Performing sterilization .....	16
3.11	Inspection .....	17
3.11.1	Visual inspection.....	17
3.11.2	Functional test .....	17
3.12	Storage .....	17
3.12.1	Ambient conditions.....	17
4.	Repair, maintenance and service.....	17
4.1	Repair.....	17
4.2	Maintenance .....	17
4.3	Technical service.....	18
4.4	Accessories / spare parts .....	18
5.	Disposal .....	18
6.	Technical data.....	18
6.1	Classification acc. to Regulation (EU) 2017/745 .....	18
6.2	Performance data, information about standards.....	18
7.	Symbols on product and packaging .....	19

### 1. About this document

#### Note

*General risk factors associated with surgical procedures are not described in these instructions for use.*

#### 1.1 Scope

These instructions for use apply for the following products:

Art. no.	Designation
PV481	Full HD CMOS camera head with pendulum coupler
PV482	Full HD CMOS camera head with zoom coupler
PV485	Full HD 3CMOS camera head with zoom coupler
JF441R	Tray for 2D camera heads

These instructions for use are an integral component of the device and contain all the information required by users and operators for safe and proper use.

#### Target group

These instructions for use are intended for physicians, medical assistants, medical technicians and employees of sterile services who are entrusted with the installation, operation, maintenance and reprocessing of the device.

#### Using and storing this document

These instructions for use must be stored in a defined location so that it may be accessed at all times by the target group.

In the event of the sale of this device or its relocation, this document must be handed over to the new owner.

## Supplementary documents

For the use and assignment of the camera head buttons and the related options and possible functions, observe the instructions for use of the camera control unit concerned.

The instructions for use of the camera head used and instructions for use of all other devices used must be observed to ensure the safe use of the device.

- For article-specific instructions for use as well as information on material compatibility see B. Braun eIFU at [eifu.bbraun.com](http://eifu.bbraun.com)

## 1.2 Safety messages

Safety messages make clear the dangers to patient, user and/or product that could arise during the use of the product. Safety messages are labeled as follows:

**DANGER**  
Indicates a possible threat of danger. If not avoided, death or serious injury may result.

**WARNING**  
Indicates a possible threat of danger. If not avoided, minor or moderate injury may result.

**CAUTION**  
Indicates a possible threat of material damage. If not avoided, the product may be damaged.

## 2. Clinical use

### 2.1 Product description

#### 2.1.1 Operating principle

The camera heads within the scope of these instructions for use can be used with the PV480 and PV630 camera control units and may only be connected to this unit.

The combination of camera control unit and camera head result in a camera based on CMOS technology.

In combination with a suitable monitor, the camera provides two-dimensional images.

The camera head features four configurable camera head buttons that can be used to control image display, record images and video, or to navigate and change the settings in the configuration menu of the camera control unit.

The camera heads are intended for use with the applicable sterile cover. (B. Braun JG904).

#### Camera head buttons

The camera head buttons have different functions depending on the operating mode.

The assignment of the buttons is always synchronous with the assignment of the buttons on the front of the camera control unit.

Different functions can be assigned for a short press of the button (<2 sec.) and a long press of the button (≥2 sec.).

The procedure for changing the assigned functions is described in the instructions for use of the respective camera control unit.

### Button assignment in live mode

The button assignment in live mode depends on whether a standard profile or user profile is set.

If a user profile is set, the buttons can be assigned individually, however the assignment for entering the camera control unit menu cannot be changed. In standard profile, the buttons are assigned as follows:

	Short press of button	Long press of button
	not working	Open menu
	Zoom	White balance
	Light source on	Light source on/off
	Single frame image	Video recording start/stop

### Button assignment when on-screen display is active

When on-screen display is active, the buttons are used to navigate the menu until the menu is left.

	Short press of button	Long press of button
	up / control +	Leave the menu
	right / next level across / save control setting	
	down / control - / on	
	left / back a level / save control setting	

### Button assignment when on-screen keyboard is active

When the on-screen keyboard is active, the buttons are used to navigate the keyboard until the keyboard is left.

	Short press of button	Long press of button
	up	not working
	right	Confirm/enter
	down	not working
	left	not working

## 2.1.2 Scope of Supply

Art. no.	Designation
PV481	Full HD CMOS camera head with pendulum coupler – or –
PV482	Full HD CMOS camera head with zoom coupler – or –
PV485	Full HD 3CMOS camera head with zoom coupler
TA014624	Instructions for use

## 2.1.3 Components required for use

### ⚠ WARNING

If the recommended components are not used, the following errors/complications may occur, among others:

- Incorrect color display
- No/incorrect display of the camera menu
- Restricted operating function from the camera head
- Increased fogging on the optics

The camera head (PV481/PV482/PV485) is used for visualization during endoscopic surgical procedures.

The camera head is intended for use with the following components:

- Camera control unit PV480 or PV630
- Disposable sterile cover JG904
- LED light source OP950

These components are optimally aligned with one another and therefore offer the best possible quality and unlimited functionality.

## 2.2 Areas of use and limitations of use

### 2.2.1 Intended use

#### Camera head (PV481/PV482/PV485)

The Full HD CMOS camera control unit (PV480) serves the 2D HD visualization of the interior of the body during minimally invasive surgical procedures and endoscopic examinations.

The application in combination with a camera head serves 2D visualization of the intracorporeal surgical area during endoscopic diagnostic and surgical procedures.

The Full HD CMOS camera control unit is only used for visualizing and not for diagnosis.

#### Tray for 2D camera heads

Aesculap sterilizer baskets incl. accessories made of silicone/plastic or with silicone/plastic components are used by OP and AEMP personnel and are dimensionally stable containers intended for repeated use. They are used for holding goods to be sterilized, for transport and storage in a sterile barrier system (e.g. sterilization containers, etc.).

### 2.2.2 Endoscopic indications

#### Note

The manufacturer is not responsible for any use of the product against the specified indications and/or the described applications.

#### Camera head (PV481/PV482/PV485)

- Minimally invasive procedures
- General surgery
- General endoscopy

#### Tray for 2D camera heads

For indications, see Intended use.

### 2.2.3 Contraindications

#### Camera head (PV481/PV482/PV485)

The use of the Full HD CMOS camera components and its accessory is contraindicated if endoscopic procedures are contraindicated for any reason. As is appropriate with any surgical procedure, consideration must be given to patient size and workspace volume when using the Full HD CMOS camera components.

Depending on the patient's disease, there may exist contraindications that rely on the general condition of the patient or the specific disease pattern. The decision to perform an endoscopic procedure rests with the responsible surgeon and should be made on the basis of an individual risk-benefit analysis.

#### Tray for 2D camera heads

No known contraindications.

## 2.3 Safety information

### 2.3.1 Clinical user

#### General safety information

To prevent damage caused by improper setup or operation, and to not compromise the manufacturer warranty and liability:

- ▶ Use the product only according to these instructions for use.
- ▶ Follow the safety and maintenance instructions.
- ▶ Ensure that the product and its accessories are operated and used only by persons with the requisite training, knowledge and experience.
- ▶ Store any new or unused products in a dry, clean, and safe place.
- ▶ Prior to use, check that the product is in good working order.
- ▶ Keep the instructions for use accessible for the user.

#### Note

The user is obligated to report all severe events in connection with the product to the manufacturer and the responsible authorities of the state in which the user is located.

#### Notes on surgical procedures

It is the user's responsibility to ensure that the surgical procedure is performed correctly.

Appropriate clinical training as well as a theoretical and practical proficiency of all the required operating techniques, including the use of this product, are prerequisites for the successful use of this product.

The user is required to obtain information from the manufacturer if there is an unclear preoperative situation regarding the use of the product.

## 2.3.2 Product

### Product-specific safety information

Risk to user and patient from non-observance of instructions, warnings and precautions!

- ▶ Use the product only according to these instructions for use.
- ▶ Check the product is in good working order before use.
- ▶ Carry out function checks before each use.
- ▶ Do not use the product if irregularities are observed during inspections and checks.

Risk to patient due to incorrect application!

- ▶ Use the product only after instruction by the manufacturer or by an authorized person.
- ▶ The instructions for use enclosed with the individual components and all products used (e.g. high frequency surgery) must be observed.
- ▶ Endoscopic procedures may only be carried out by specialists who have the necessary medical training, knowledge and experience.

Risk to patient and user due to premature wear!

- ▶ Handle and maintain the product properly.
- ▶ Use the product only in accordance with the intended purpose.

Possible restricted function if third party devices are used!

- ▶ Use the product with the recommended components and accessories.
- ▶ The full functionality of the product can only be guaranteed if the recommended accessory components are used.

Risk of injury due to electric shock!

- ▶ When installing a medical electrical system, there is a risk of fire, short-circuit or electric shock. Installation must be performed by qualified personnel only.
- ▶ If a number of electrical devices are being used in combination, observe Annex I of IEC 60601-1:2006. Non-medical devices that comply with the applicable IEC safety standards can be connected only via a medical isolating transformer. Do not connect any additional non-medical devices to a medical electrical system.
- ▶ Signal lines from devices with a functional connection, which are connected to different branches of the mains power supply, must be galvanically isolated at both ends.
- ▶ Connect the devices only to a power supply with a protective earth conductor.
- ▶ After installing a medical electrical system, it must be inspected according to IEC 62353.

Risk to the patient if the product is not used gently/a damaged product is used!

- ▶ Handle the product with appropriate care.
- ▶ Do not use the product if it has been subjected to strong mechanical stress or if it has been dropped, and send the product to the manufacturer or an authorized repair center for inspection.

Risk to patient if the device power fails!

- ▶ Use the device only with an uninterruptible power supply.
- ▶ To guarantee a continuous power supply a medical emergency supply (USP) is suggested.

Device failure due to incorrect storage and usage conditions!

- ▶ Store and operate the product only within the specified environmental conditions.

Restricted function if third party devices are used!

- ▶ Use the product only with the recommended components and accessories.
- ▶ Full functionality is only guaranteed if the recommended components and accessories are used.

#### Note

*Accessories and/or peripheral devices which are connected to the interfaces of the medical device must be verified to comply with the relevant specifications (e.g. IEC 60601-1).*

#### Note

*The medical electrical device may not be modified in any way.*

#### Note

*To guarantee optimal function of the product, use in controlled environmental conditions is recommended (e.g. air-conditioned operating room).*

#### Note

*Always handle the medical device very carefully as it contains sensitive optic, mechanic and electronic components. Do not strike or drop the camera head.*

#### Note

*Ensure that all devices operated in the vicinity meet their relevant EMC requirements.*

#### Note

*The camera components are used to visualize the inside of the body during minimally invasive procedures. Do not use the components for diagnostic purposes. This applies in particular with the use of image optimization algorithms.*

#### Note

*Before putting into operation, check the compatibility of all components using the accessories list.*

#### Note

*All accessories and spare parts must only be procured from the manufacturer.*

- ▶ Clean (manually or mechanically) any new product that is supplied from the factory after it has been removed from the transport packaging and before first sterilization.
- ▶ Prior to use, check that the product is in good working order.
- ▶ Observe "Notes on the electromagnetic compatibility (EMC) for the PV480 2D camera platform" TA022467, see B. Braun eIFU at eifu.bbraun.com
- ▶ Only combine Aesculap products with each other.
- ▶ Always adhere to applicable standards.

#### Ambient conditions

The following environmental conditions apply to the use of the product:

Temperature	10 °C to 35 °C (PV481) 10 °C to 37 °C (PV482/PV485)
Relative humidity	30 % to 90 %
Atmospheric pressure	700 hPa to 1 060 hPa

### 2.3.3 Sterility

The product is delivered in an unsterile condition.

- Clean the new product after removing its transport packaging and prior to its initial sterilization.

## 2.4 Preparation

If the following instructions are not observed, the manufacturer assumes no responsibility for possible consequences.

Before installation and use, ensure that the:

- electrical installations comply with the relevant technical regulations,
- relevant regulations concerning fire and explosion protection are observed.

#### *Note*

*For the safety of patients and users it is essential that the mains power cord and, especially, the protective ground connection are intact. In many cases defective or missing protective ground connections are not recognized immediately.*

### ⚠ WARNING

Risk to patient from misdirected leakage current due to incorrect or defective earthing!

- Do not touch the product and patient at the same time.

### ⚠ WARNING

Risk to patient due to failure of a device!

- Keep an operationally ready replacement device on hand and if necessary change to conventional surgical methods.

### ⚠ WARNING

Danger of explosion from a device that is improperly set up!

- Make certain that the power plug is connected to the power supply outside of any areas where there is a risk of explosion.
- Do not use the product in potentially explosive areas or close to highly flammable or explosive gases (e.g. oxygen, anesthetic gases).

### ⚠ WARNING

Risk to patient and user due to condensation and short-circuit!

- Before putting into operation make sure that all components used have enough time to adjust to the changed environmental conditions.

### ⚠ WARNING

Risk of injury to the eyes!

- Do not insert the light guide into the light source during the check of the fibers.

### ⚠ WARNING

Risk of infection and injury!

- Do not use products with damaged fiber optics, damaged glass surfaces or stubborn deposits which cannot be removed by cleaning.

### ⚠ WARNING

Risk to patient due to use of damaged products!

- Do not use products if they have sharp edges or other hazardous surface damage.

### ⚠ CAUTION

Risk to persons and risk of damage to equipment due to improper cable routing!

- Lay all cables and leads so that they do not present a tripping hazard.
- Do not place objects on the cables.

### ⚠ CAUTION

Risk to patient due to poor/lack of vision!

- Adjust the monitors and display elements so that they are highly visible to the user.

### ⚠ CAUTION

Interference of high frequency energy with the device!

- Do not use any mobile or portable equipment that emit high-frequency energy (e.g. mobile phones, cell phones, GSM phones) in the vicinity of the product.

#### *Note*

*The product including the camera cable is fragile when bent, kinked, twisted, pulled, or pressed.*

#### *Note*

*Handle the camera head and camera cable with care. Do not exert any mechanical forces and avoid mechanical impacts.*

#### *Note*

*Connect the equipotential bonding terminals of all equipment used with the equipotential bonding strip see IEC 60601-1-1 / EN 60601-1-1 or in accordance with national standards.*

#### *Note*

*Ensure that the corresponding interconnection conditions, standards and respective national deviations are observed.*

### 2.4.1 Combination with medical electrical equipment

The product can be combined with components from other manufacturers provided that all components are compliant with the medical electrical equipment safety requirements according to IEC 60601-1.

It is the operator's responsibility to check and make sure that the system is and remains fully operational.

When using devices from different manufacturers and when operating an endoscope and/or endoscopic accessories together with medical electrical equipment, it must be ensured that the applied part is properly isolated: Type CF, defibrillation-proof.

### 2.4.2 Inspections

Carry out the following inspection steps before reprocessing and immediately before use of the device:

#### Inspection of the glass surfaces

- Before each use, inspect the camera head for damage such as rough surfaces, sharp edges or protrusions to avoid injuries to the patient.
- Perform a visual inspection of the glass surfaces. The surfaces must be clean and smooth.

If there are impairments/damage, see Troubleshooting.

## 2.5 First use

### ⚠ WARNING

Risk of injury and/or product malfunction due to incorrect operation of the medical electrical system!

- ▶ Follow the instructions for use of all products required.

## 2.6 Application

### 2.6.1 General notes on use

The camera cable is fragile when bent, kinked, twisted, pulled, or pressed. Observe the maximum bend radius of 7 cm and handle the cable with care. Do not make any changes to the cable and do not load with sharp-edged objects.

### 2.6.2 System set-up

### ⚠ WARNING

Risk of infection due to unsterile parts!

- ▶ Ensure that unsterile components do not enter the sterile area.
- ▶ Process products and accessories that are delivered unsterile before use and only use them with sterile accessories.

### ⚠ WARNING

Risk of infection due to unsterile, soiled or contaminated products!

- ▶ Only use the camera head in the sterile area if it is in sterile condition. To do so, use the camera head with the correctly applied sterile cover or after the properly performed sterilization process.

### ⚠ WARNING

Risk to patient due to addition of leakage currents!

- ▶ If the product is used with medical electrical equipment and/or power-driven endoscope accessories, the leakage currents can multiply.

### ⚠ WARNING

Risk to patient and user from burns, sparks or explosion!

- ▶ When using a high frequency surgical device during an endoscopic procedure, observe the safety instructions of the corresponding instructions for use.

### ⚠ WARNING

Risk of injury due to improper use!

- ▶ Observe the intended use of the product.
- ▶ Do not use the endoscope as a lever.

### ⚠ CAUTION

Interference of electromagnetic emissions with the image quality (e.g. minor banding, minor color changes on the monitor image)!

- ▶ Check the image quality if used in combination with additional peripheral devices (e.g. monitor, video equipment).

### Note

In combination with the camera control unit, PV480 or PV630, the endoscope is classed as a type CF defibrillation-proof applied part.

### Note

Before application of endoscopic high frequency surgery (HF surgery) prepare the patient accordingly.

### Note

*Take measures to remove or avoid formation of combustible gases (e.g. gas-tro-intestinal tract/colonoscopy, urinary bladder/transurethral resection).*

### Connecting the accessories

### ⚠ DANGER

Risk of injury due to unapproved configuration using additional components!

- ▶ For all components used, ensure that their classification matches that of the applied part (e.g. type CF defibrillation-proof) of the respective device.

Combinations of accessories that are not mentioned in the present instructions for use may only be employed if they are specifically intended for the respective application, and if they do not compromise the performance and safety characteristics of the products.

Any equipment connected at the interfaces must also demonstrably meet the respective IEC standards (e.g. IEC 60950 for data processing equipment, and IEC/DIN EN 60601-1 medical electrical equipment).

All configurations must meet the basic standard IEC/DIN EN 60601-1. The person connecting the devices is responsible for the configuration and must ensure that the basic standard IEC/DIN EN 60601-1 or the corresponding national standards are met.

- ▶ Please contact your B. Braun/Aesculap partner or Aesculap Technical Service, (address: see Technical service) with any inquiries in this respect.

### 2.6.3 Function checks

### ⚠ WARNING

Risk to patient due to incorrectly displayed image!

- ▶ Remove any contaminations on the optical surfaces (camera head) before and during the procedure. Clean the endoscope tip with a cloth soaked in alcohol (70 % ethanol) or with a neutral cleaning agent.
- ▶ Prior to application and after changing the settings, check the correct display of the live image. If necessary, perform white balance.

### ⚠ WARNING

Risk to patient and user from functional failure or electric shock!

- ▶ Use only dry components (e.g. plug connector for camera control unit, optical lens).

### Note

*Only perform the surgical procedure if all of the components being used are in a flawless condition.*

### Note

*Before each use, after a power failure or any interruption, the functionality of all connected devices must be checked and all connections must be checked to see they are correct.*

### Note

*The endoscope image must be focused, bright and clean at an appropriate working distance.*

### Note

*Do not continue to use damaged products.*

## 2.6.4 Safe operation

### **⚠ WARNING**

Risk of burns due to development of heat at the tip of the endoscope!

- ▶ Do not put down the endoscope on the patient during use.
- ▶ In the patient's abdominal cavity, always leave a sufficient distance between the optical lens and the tissue surfaces and mucous membranes of the patient.
- ▶ Use automatic light intensity control or set the light source so that a bright, well-lit image is visible with the lowest possible light intensity.
- ▶ Switch off the light source if the illumination is no longer needed or if the endoscope is out of the patient for a longer period of time.

### Excessive temperatures in combination with light sources

Light sources, especially high-power light sources, emit large amounts of light energy and thermal energy. The light guide connector and the distal end of the endoscope may therefore become extremely hot.

### **⚠ WARNING**

Risk of injury due to excessive temperature!

- ▶ Do not touch the light guide connector or distal end of the endoscope during use or immediately after use.

Risks from use of light sources:

- Irreversible tissue damage or unwanted coagulation to the patient or user
- Burns or thermal damage to surgical equipment (e.g. surgical drapes, plastic materials, etc.)
- If the light source fails during use, this may endanger the patient. Therefore keep an operationally ready replacement light source on hand.

### Safety precautions

- ▶ Do not illuminate the inside of the patient for longer than necessary with the light source.
- ▶ Use the automatic light intensity control or set the light source so that a bright, well-lit image is visible with the lowest possible light intensity.
- ▶ Do not allow the distal end of the endoscope or light guide connector to come into contact with patient tissue or with combustible or heat-sensitive materials.
- ▶ Do not touch the distal end of the endoscope.
- ▶ Remove contaminations on the distal end surface or light emission surface.

### **⚠ WARNING**

Risk of injury and/or malfunction!

- ▶ Always carry out a function check prior to using the product.

### **⚠ WARNING**

Risk of injury from burns and unwanted deep penetration and risk of damage to product!

- ▶ Only switch on high frequency current when the corresponding applied part (electrode) can be seen through the endoscope and there is no contact between the two.

### **⚠ WARNING**

Risk of infection due to contamination of the sterile area!

- ▶ Fix the camera cable (unsterile) with sterile cover (sterile) sufficiently nearby in the operation area.

### **⚠ WARNING**

Risk of injury due to misdirected leakage currents!

- ▶ Before thoracic or cardiothoracic surgery, deactivate implanted defibrillators (ICD).
- ▶ Before any kind of defibrillation, remove the used camera head from the patient.

### **⚠ WARNING**

Risk of infection/injury to patients and users!

- ▶ Only use the sterile cover once.
- ▶ Do not clean the sterile cover with ultrasound.
- ▶ Do not reprocess the sterile cover.

### **⚠ WARNING**

Risk to patient from gas embolism!

- ▶ Avoid overinsufflation (e.g. with air or inert gas) before high frequency surgery.

### **⚠ WARNING**

Risk to patient from damage to the optical lens!

- ▶ Withdraw the endoscope slowly from the used trocar.
- ▶ If a trocar with a manual opening valve is being used, only withdraw the endoscope from the trocar when the valve is open.

### **⚠ WARNING**

Risk of infection due to unsterile cover or medical device!

- ▶ Replace the sterile cover or the camera head if it falls on the ground or if it is touched by unsterile objects/persons.

### **⚠ CAUTION**

Malfunction in connection with magnetic resonance!

- ▶ Do not use the product in a magnetic resonance environment.

### Note

The current assignment of the buttons is displayed on the screen as "Camera head information" after opening the camera menu.

### Connection to camera control unit

- ▶ Insert the connection plug fully into the connection socket of the camera control unit until it engages in place.

### Connecting the endoscope

- ▶ Press both endocoupler latches 2 at the same time in order to spread apart the holding clamps.
- ▶ Insert the endoscope and release the endocoupler latches 2.
- ▶ In order to avoid inadvertent opening of the endocoupler: Rotate safety lock 1 to the locked position.

### Connecting the light cable

- ▶ Connect the light cable to the endoscope.
- ▶ Switch on the light source at a low light intensity setting.

## Adjusting the image focus and enlargement

- Adjust the intensity of the light source until sufficient illumination is achieved.
- In order to sharpen the image, turn the focus ring **3** on the camera head.
- In order to set the desired image size, turn the zoom ring **4** on the camera head (only PV482 and PV485).

## Pendulum function (only PV481)

This camera head is self-aligning. The image always remains upright as the camera head remains aligned perpendicularly to the endoscope axis due to the force of gravity.

- In order to inhibit rotatability around the endoscope axis: Actuate lock **6**.

## 2.7 Troubleshooting

Problem	Possible cause	Remedy
Image cloudy	Glass surfaces soiled	Clean the glass surfaces according to the Cleaning and Disinfection section.
	Stubborn residue on the glass surfaces	Remove residues according to the Cleaning and Disinfection section; check water quality.
Image too dark, too little illumination	Glass surfaces soiled	Clean the glass surfaces according to the Cleaning and Disinfection section.
	Light guide dirty, defective	Check light guide (e. g. illuminate onto a white surface).
	Excessive chloride concentration	Check water quality.
	Heavy metal ions and/or silicates, increased content of iron, copper, manganese in water	Check water quality; if applicable, only use deionized (fully desalinated) water.
	Too high concentration of mineral substances (e.g. calcium) or organic substances	
	Contaminated cleaning/disinfecting solutions, too frequent use	The cleaning and disinfecting solutions should be replaced regularly.
	Extraneous rust, (e.g. resulting from previously damaged or non-corrosion-resistant instruments being reprocessed at the same time)	Check supply systems; in case of joint reprocessing, check material compatibility and for existing damage, and avoid mutual contact.
	Contact corrosion	Avoid contact with other components.

### 2.7.1 Repairs

For repairs, contact the manufacturer or authorized repair center. Authorized repair centers can be inquired about from the manufacturer.

For a fast processing of your service requests, send in the product indicating:

- Article number (REF)
- Serial number (SN)
- Detailed description of defects

#### Note

If a component of the visualization system is replaced (e. g. service, upgrade), the system start-up process must be performed again. It is also recommended to perform servicing after any maintenance or replacement of a device.

### 3. Reprocessing procedure

#### 3.1 General safety instructions

**Note**

Adhere to national statutory regulations, national and international standards and directives, and local, clinical hygiene instructions for sterile processing.

**Note**

For patients with Creutzfeldt-Jakob disease (CJD), suspected CJD or possible variants of CJD, observe the relevant national regulations concerning the reprocessing of products.

**Note**

The camera head is not autoclavable. Therefore, the camera head and the corresponding sterile cover should not be used in patients who are suspected of having Creutzfeldt-Jacob disease or who have already been diagnosed with Creutzfeldt-Jacob disease.

**Note**

Mechanical reprocessing should be favored over manual cleaning as it gives better and more reliable results.

**Note**

Successful processing of this medical device can only be ensured if the processing method is first validated. The operator/sterile processing technician is responsible for this.

**Note**

For up-to-date information about reprocessing and material compatibility, see B. Braun eifu at [eifu.bbraun.com](http://eifu.bbraun.com)

The validated steam sterilization procedure was carried out in the Aesculap sterile container system.

#### 3.2 General information

Dried or affixed surgical residues can make cleaning more difficult or ineffective and lead to corrosion. Therefore the time interval between application and reprocessing should not exceed 6 h; also, neither fixating pre-cleaning temperatures >45 °C nor fixating disinfecting agents (active ingredient: aldehydes/alcohols) should be used.

Excessive neutralizing agents or basic cleaners may result in a chemical attack and/or fading and the laser marking becoming unreadable either visually or by machine in case of non stainless steel.

In the case of stainless steel, residues containing chlorine or chlorides (e. g., in surgical residues, medicines, saline solutions, and in the service water used for cleaning, disinfection, and sterilization) will cause corrosion damage (pitting, stress corrosion) and result in damage to the products. These must be removed by rinsing thoroughly with demineralized water and then drying.

Perform additional drying, if necessary.

Only process chemicals that have been tested and approved (e.g. VAH or FDA approval or CE mark) and which are compatible with the product's materials according to the chemical manufacturers' recommendations may be used for reprocessing the product. All the chemical manufacturer's application specifications must be strictly observed. Failure to do so can result in the following problems:

- Optical material changes (e.g. fading or discoloration) in titanium or aluminum. For aluminum, the application/process solution only needs to be pH >8 to cause visible surface changes.
- Material damage (e.g. corrosion, cracks, fracturing, premature aging, or swelling)
- ▶ Do not use metal cleaning brushes or other abrasives that would damage the product surface and could cause corrosion.
- ▶ Further detailed advice on hygienically safe and material-/value-preserving reprocessing can be found at [www.a-k-i.org](http://www.a-k-i.org), link to "AKI-Brochures", "Red brochure".

#### 3.3 Dismantling prior to carrying out the reprocessing procedure

- ▶ Unplug the camera plug from the camera control unit.
- ▶ Remove the sterile cover from the camera head and dispose of it according to regulations. The sterile cover is only intended for single use.
- ▶ Remove the endoscope, light guide and all detachable parts.

#### 3.4 Preparations at the place of use

- ▶ Remove any visible surgical residues as much as possible with a damp, lint-free cloth. Use a soft cloth for this to avoid scratching the optical lens.
- ▶ Pre-clean all used parts of the fully disassembled product at the operating table at the conclusion of the surgery by wiping with a lint-free wipe wetted with an enzymatic cleaning solution until visibly clean.
- ▶ Place the dry product in a sealed waste container and forward it on for cleaning and disinfection within 6 hours.

#### 3.5 Preparation before cleaning

- ▶ Non-fixating/NaCl-free pre-cleaning is mandatory immediately after use.

#### 3.6 Cleaning/disinfection

The effectiveness of the following procedures:

- Manual cleaning and manual disinfection
- Automated cleaning and thermal disinfection
- Sterilization

As described in this document, has been fully validated.

It is the operator's responsibility to introduce, document, implement, and maintain a validated reprocessing procedure. Make sure that the equipment used for reprocessing is properly maintained.

The reprocessing procedure described in this document consists of the following steps:

- Pre-cleaning immediately after use
- Cleaning and disinfection (manual or automated) and
- Sterilization

### 3.6.1 Product-specific safety instructions for the reprocessing procedure

Risk of electric shock and fire hazard!

- ▶ Unplug the device before cleaning.
- ▶ Do not use flammable or explosive cleaning or disinfecting solutions.

Risk of infection for patient and/or user!

- ▶ Ensure that there is no residue from cleaning and disinfecting agents on the product.
- ▶ Avoid insufficient or incorrect cleaning and disinfection of the product and accessory.

Damage to the product due to inappropriate cleaning/disinfecting agents and/or excessive temperatures!

- ▶ Use cleaning and disinfecting agents according to the manufacturer's instructions. The cleaning and disinfecting agent must be approved for plastics and high-grade steel, and it must not attack softeners (e.g. silicone).
- ▶ Observe specifications regarding concentration, temperature and exposure time.
- ▶ Do not exceed the maximum permitted cleaning temperature of 60 °C.

Damage to the product due to improper reprocessing!

- ▶ Do not, under any circumstances, clean or disinfect the product in an ultrasonic cleaning bath.
- ▶ Use cleaning and disinfecting agents which are suitable and approved for the present product.
- ▶ Observe the manufacturer's cleaning and disinfecting instructions regarding concentration, temperature and exposure time.
- ▶ Only reconnect the system to power when all cleaned parts are completely dry.
- ▶ Never autoclave the product.

Risk of defects, consequential damage or a shortened product service life!

- ▶ Follow and observe the manufacturer's requirements for reprocessing.

#### Note

*There may be significant discoloration of the camera head when using Cidex OPA. However, this discoloration has no influence on the functionality or safety of the product.*

#### Note

*The silicone elements of the JF441R reprocessing basket can get brittle on the surface over the course of its service life.*

*If signs of brittleness appear, replace with the replacement set JF441500 the silicone elements.*

### 3.6.2 Validated cleaning and disinfection procedure

Validated procedure	Specific requirements	Reference
Manual cleaning with immersion disinfection ■ Camera head PV481 / PV482 / PV485	Drying phase: Use a lint-free cloth or medical compressed air ■	Chapter Manual cleaning and disinfecting and subsection: ■ Chapter Manual cleaning with immersion disinfection
Mechanical alkaline cleaning and thermal disinfection ■ Tray for 2D camera heads JF441R	None	Chapter Mechanical cleaning/disinfection and subsection: ■ Chapter Mechanical alkaline cleaning and thermal disinfection
Manual pre-cleaning with brush and subsequent mechanical alkaline cleaning and thermal disinfection ■ Camera head PV481 / PV482 / PV485	Place the product on tray JF441R (avoid rinsing blind spots). ■	Chapter Mechanical cleaning/disinfecting with manual pre-cleaning and subsection: ■ Chapter Manual pre-cleaning ■ Chapter Mechanical alkaline cleaning and thermal disinfecting

### 3.7 Manual cleaning and disinfecting

#### 3.7.1 Post use and pre-cleaning activities

- Pre-clean the device immediately after use.
- Remove and discard the sterile drape, if present.
- Decouple the endoscope from the optocoupler.
- Uncouple the camera head from the controller.
- Pre-clean all used parts of the fully disassembled product at the operating table at the conclusion of the surgery by wiping with a lint-free wipe wetted (i.e. wetted but not dripping; squeeze out excess solution) with an enzymatic cleaning solution until visibly clean. Before wetting the lint-free wipe with the enzymatic clean solution, the solution should be prepared according to the manufacturer's instruction.
- Arrange for reprocessing, ensuring that all of the device components are reprocessed within 6 hours.

#### 3.7.2 Manual cleaning with immersion disinfection

##### Note

*On the basis of material tests carried out, full functionality can be ensured up to 300 reprocessing cycles.*

*Careful visual and functional testing prior to next use is the best way to identify a malfunctioning product, see Inspection.*

Phase	Step	T [°C/°F]	t [min]	Water quality	Chemicals
I	Cleaning	as per manufacturer's instruction	2-5	Tap water	Enzymatic cleaning agent (Cidezyme/Enzol)
II	Rinse 2x	<45/113	2x ≥1	Tap water	-
III	Disinfection <sup>1)</sup>	as per manufacturer's instruction	12	as per manufacturer's instruction	0.55 % Orthophtaldehyde solution (Cidex OPA)
IV	Rinse 2x <sup>1)</sup>	<45/113	2x ≥1	Tap water	-
V	Final Rinse	<45/113	≥1	FD-W	-
VI	Drying	-	-	-	-

FD-W Fully desalinated water (demineralized, low microbiological contamination, max. 10 germs/ml and low in endotoxins, max. 0.25 Endotoxin units/ml)

<sup>1)</sup> Phase does not apply to US market

## **⚠ CAUTION**

**Damage to the product due to improper handling!**

**Scratch-sensitive surface.**

- ▶ Handle the device with appropriate care.
- ▶ Do not use metal brushes, metal objects or abrasive cleaners.

### **Phase I: Cleaning**

- ▶ Prepare the cleaning solution according to the manufacturer's instructions.
- ▶ Fully immerse all parts of the fully disassembled device that were used in the cleaning solution bath for two (2) to five (5) minutes. All accessible surfaces must remain immersed in the cleaning solution bath throughout the entire cleaning time.
- ▶ After soaking in the cleaning solution, clean all external surfaces - while the items are fully immersed in the cleaning solution – using a clean, lint-free, soft, non-sterile cloth or a soft bristled brush until all visible soiling is removed:
  - Brush for at least one (1) minute or until no more residues can be removed.
  - During cleaning, move non fixed components 3 times in each direction and as far as possible.
- ▶ Thoroughly rinse all surfaces of products with hidden crevices, lumens with working channels, or complex geometry five (5) times. Use a disposable 50 ml syringe.

### **Phase II: Rinse**

- ▶ Fully immerse all parts in a tap water bath (< 45 °C/113 °F) and thoroughly rinse all accessible surfaces twice for a minimum of one minute each.
- ▶ Move non fixed components 3 times in each direction and as far as possible.
- ▶ Thoroughly rinse all hidden crevices, lumens or complex geometry at least three times:
  - Use a disposable 50ml syringe.
  - Use fresh water for each rinse.
  - Allow water to drip off for a sufficient length of time.

### **Phase III: Disinfection**

- ▶ Immerse all parts completely in disinfectant solution for at least 12 minutes. All accessible surfaces must remain immersed in the disinfectant solution bath throughout the entire disinfecting time.
- ▶ Remove all adherent air bubbles from the component surfaces
- ▶ Thoroughly rinse all surfaces of products with hidden crevices, lumens with working channels, or complex geometry five (5) times. Use a disposable 50 ml syringe.
- ▶ Move non fixed components 3 times in each direction and as far as possible.

### **Phase IV: Rinse**

- ▶ Fully immerse all parts in a tap water bath (<45 °C/113 °F) and thoroughly rinse all accessible surfaces twice for a minimum of one (1) minute each.
- ▶ Move non fixed components 3 times in each direction and as far as possible.
- ▶ Thoroughly rinse all hidden crevices, lumens or complex geometry at least three times.
  - Use a disposable 50ml syringe.
  - Use fresh water for each rinse.
  - Allow water to drip off for a sufficient length of time.

### **Phase V: Final rinse**

- ▶ Fully immerse all parts in a fully desalinated water bath (<45 °C/113 °F) and thoroughly rinse all accessible surfaces for a minimum of one minute.
- ▶ Move non fixed components 3 times in each direction and as far as possible.
- ▶ Thoroughly rinse all hidden crevices, lumens or complex geometry at least three times:
  - Use a disposable 50 ml syringe.
  - Use fresh water for each rinse.
  - Allow water to drip off for a sufficient length of time.

### **Phase VI: Drying**

- ▶ Thoroughly dry all parts with a clean, lint-free wipe or a lint-free surgical towel.
- ▶ Dry all accessible surfaces, especially the channels with medical-quality filtered compressed air ( $p_{max} = 0,5$  bar).
- ▶ Visually inspect the device in a well-lit area; they should be completely clean and dry. If necessary please use a magnifying glass and repeat the manual cleaning performance.

This concludes the manual cleaning and disinfection process.

### 3.8 Mechanical cleaning/disinfection

*Note*

The cleaning and disinfection device must be tested and approved effectiveness (e.g. FDA approval or CE mark according to DIN EN ISO 15883).

*Note*

The cleaning and disinfection device used for processing must be serviced and checked at regular intervals.

Machine type: single-chamber cleaning/disinfection device without ultrasound

#### 3.8.1 Mechanical alkaline cleaning and thermal disinfection

Phase	Step	T [°C/°F]	t [min]	Water quality	Chemicals/Note
I	Prerinse	<25/77	3	Tap water	-
II	Cleaning	55/131	10	Deionized water	Alkaline Detergent Neodisher Mediclean forte 0.5 %
III	Intermediate rinse	>10/50	1	Deionized water	-
IV	Thermal disinfection	90/194	5	Deionized water	-
V	Drying	-	-	-	According to the program of the washer/disinfector

- Check visible surfaces for residues after mechanical cleaning/disinfecting.

### 3.9 Mechanical cleaning/disinfecting with manual pre-cleaning

#### Note

*On the basis of material tests carried out, full functionality can be ensured up to 240 reprocessing cycles.*

*Careful visual and functional testing prior to next use is the best way to identify a malfunctioning product, see Inspection.*

#### 3.9.1 Post use and pre-cleaning activities

- If applicable, flush not observably surfaces primarily with deionized water e.g. with a single use syringe.
- Remove as far as possible visible residues with a moist lint free wipe. Use an enzymatic cleaning solution to moisten the cloth that was prepared in accordance with the instructions of the manufacturer.
- Place the device in a dry and closed disposable container and transport it within 6 hours to cleaning and disinfection.

#### 3.9.2 Manual pre-cleaning

Phase	Step	T [°C/°F]	t [min]	Water quality	Chemicals
I	Cleaning	<45/113	10-30	Tap water	Enzymatic cleaning agent (Cidezyme/Enzol)
II	Rinse 2x	<45/113	2x ≥1	Tap water	-

#### Phase I: Cleaning

- Fully immerse all parts of the fully disassembled device that were used in a cleaning solution bath (<40°C) for ten (10) to thirty (30) minutes.
- Immerse all accessible surfaces in the cleaning solution bath throughout the entire cleaning time.
- While the components soak in the solution, use a soft wipe or soft brush to remove all visible residues from all exterior surfaces of the components.
- Move non fixed components 5 times in each direction and as far as possible.

#### Phase 2: Rinse

- Fully immerse all parts in a tap water bath (<45 °C/113 °F) and thoroughly rinse all accessible surfaces twice for a minimum of one (1) minute each:
  - Use fresh water for each rinse.
  - Move non fixed components 5 times in each direction and as far as possible.
  - Allow water to drip off for a sufficient length of time.

### 3.9.3 Mechanical alkaline cleaning and thermal disinfecting

Phase	Step	T [°C/F]	t [min]	Water quality	Chemicals/Note
I	Prerinse	<25/77	3	Tap water	-
II	Cleaning	55/131	10	Deionized water	Alkaline Detergent Neodisher Mediclean forte 0.5 %
III	Intermediate rinse	>10/50	1	Deionized water	-
IV	Thermal disinfection	90/194	5	Deionized water	-
V	Drying	-	-	-	According to the program of the washer/disinfector

- ▶ Use a washer/disinfector that has been validated as effective and meets the performance requirements of ISO 15883-1 (or the respective country specific version thereof)
- ▶ For reprocessing, use the camera tray, which is used without a lid. This tray is only suitable for the camera and not for instruments.
- ▶ Place the device in the provided reusable camera tray as shown in the Fig. below and than lay the cable in the brackets. Avoid areas that cannot be reached during rinsing, see Fig. 1.

Fig. 1

- ▶ Place the tray in the automatic washer/disinfector in a suitable carrier that ensures both the protection of the instrument and access to surfaces for cleaning. No lid will be used for the tray.
- ▶ Start the cleaning cycle as described (see table) in accordance with the manufacturer's instructions and instruction for use for the washer/disinfector.
- ▶ Remove the tray and the devices from automatic washer.
- ▶ Visually inspect the device in a well-lit area; they should be completely clean, dry and undamaged.
- ▶ Use a magnifying glass. Repeat the cleaning process if necessary. Immediately set aside any damaged components.

### 3.10 STERRAD Sterilization

The STERRAD® sterilization systems, manufactured by Advanced Sterilization Products (ASP), use low-temperature, hydrogen peroxide gas plasma technology for terminal sterilization of properly cleaned, rinsed, and dried reusable medical devices.

- ▶ Refer to the STERRAD® Sterilization Systems User Guide for detailed instructions for use of any STERRAD® unit and ASP's STERRAD Sterility Guide (SSG) at [www.sterradsterilityguide.com](http://www.sterradsterilityguide.com) or contact ASP customer service.

#### 3.10.1 Performing sterilization

##### Note

*On the basis of material tests carried out, full functionality can be ensured up to 240 reprocessing cycles.*

*Careful visual and functional testing prior to next use is the best way to identify a malfunctioning product, see Inspection.*

##### Note

*STERRAD® sterilization may result in cosmetic device changes that do not necessarily affect the device function.*

#### ⚠ CAUTION

**Loads containing moisture may cause a cycle cancellation!**

- ▶ Thoroughly dry the product before loading into the STERRAD® sterilizer.
  - ▶ Place all parts in the processing basket JF441R, see Mechanical alkaline cleaning and thermal disinfection.
  - ▶ Place a STERRAD® indicator strip in the basket.
  - ▶ Assemble sterilization container JM441 according to manufacturer's instruction.
  - ▶ Place the basket in the sterilization container JM441.
  - ▶ Load the container in the sterilizer:
    - Position the container such that the plasma can fully surround it.
    - Load only one container per cycle in the chamber.
    - Leave the other shelf empty.
  - ▶ Start the sterilization cycle in accordance with the manufacturer's instructions and instructions for use for the steam sterilizer.
  - ▶ STERRAD® sterilization was validated for the following cycle:
    - STERRAD® 100S short cycle
  - ▶ Remove the sterilized products from the sterilizer.
  - ▶ Ensure that the products remain sterile after reprocessing.
- This concludes sterilization.

### 3.11 Inspection

- Allow the product to cool down to room temperature.
- Dry the product if it is wet or damp.

#### 3.11.1 Visual inspection

- Make certain that all soiling has been removed. In particular, pay attention to mating surfaces, hinges, shafts, recessed areas, drill grooves and the sides of the teeth on rasps.
- If the product is dirty: repeat the cleaning and disinfection process.
- Check the product for damage, e.g. insulation or corroded, loose, bent, broken, cracked, worn or severely scratched and fractured components.
- Check the product for missing or faded labels.
- Check the products with long, slim shapes (in particular rotating instruments) for deformities.
- Check the surfaces for rough spots.
- Check the product for burrs that could damage tissue or surgical gloves.
- Check the product for loose or missing parts.
- Immediately put aside damaged or inoperative products and send them to Aesculap Technical Service, see Technical service.

#### 3.11.2 Functional test

- Check that the product functions correctly.
- Check that all moving parts are working properly (e.g. hinges, locks/latches, sliding parts etc.).
- Check the product for any atypical running noise, overheating or excessive vibration.
- Check for compatibility with associated products.
- Immediately put aside inoperative products and send them to Aesculap Technical Service, see Technical service.

### 3.12 Storage

#### ⚠ CAUTION

##### Damage to product due to improper storage!

- Store the product protected from dust in a dry, well ventilated, and temperature-controlled room.
- Store the product protected from direct sunlight, high temperatures, high air humidity or radiation.
- Protect the product from direct UV light, radioactivity or strong electromagnetic radiation.
- Store the product separately or use containers in which it can be fixed in place.
- Always transport the product with care.

#### ⚠ CAUTION

##### Damage to product due to improper handling!

- Do not drop the camera head and handle it with care.
- When transporting the camera head over thresholds or uneven ground, store the camera head securely.

#### ⚠ CAUTION

##### Damage to the optical lens due to inappropriate transport packaging!

- Use the tray only for the camera head and not for instruments.
- Use the tray for the camera head only within the clinic.
- Only place a disinfected camera head in the tray.
- Store sterile products in germ-proof packaging, protected from dust, in a dry, dark, temperature-controlled area.

#### 3.12.1 Ambient conditions

The following environmental conditions apply to the transport and storage of the product:

Temperature	-10 °C to 50 °C
Relative humidity	10 % to 90 %
Atmospheric pressure	500 hPa to 1 060 hPa

## 4. Repair, maintenance and service

### 4.1 Repair

Send damaged products to the manufacturer or authorized repair center. Authorized repair centers can be inquired about from the manufacturer.

#### ⚠ WARNING

##### Risk of infection due to soiled or contaminated products!

- Clean, disinfect and sterilize the product and/or any accessories thoroughly before sending.
- Prior to sending, remove the sterile cover from the camera head.
- Select suitable and safe packaging (ideally the original packaging).
- Package the product in such a way that the packaging does not become contaminated.

#### Note

If for urgent reasons it is not possible to thoroughly clean/disinfect and sterilize the product, reprocess the product as far as possible and label it accordingly.

The specialized repair company may refuse to repair soiled or contaminated products for safety reasons.

The manufacturer reserves the right to return contaminated products to the sender.

### 4.2 Maintenance

This medical product contains no components or parts that must be replaced in the form of service at regular intervals specified by the manufacturer.

A repeat inspection of the medical product must be performed every 12 months. This must also be carried out after any repair and after it was dropped, damaged or misused.

The repeat inspection may only be carried out by persons authorized to do so by the manufacturer, using the service manual.

- Comply with applicable national and international standards.

For services to that end, please contact your national B. Braun/Aesculap agency, see Technical service.

## 4.3 Technical service

### ⚠ DANGER

Danger to life of patients and users if the product malfunctions and/or protective measures fail or are not used!

- ▶ Do not perform any servicing or maintenance work under any circumstances while the product is being used on a patient.

### ⚠ CAUTION

Modifications carried out on medical technical equipment may result in loss of guarantee/warranty rights and forfeiture of applicable licenses.

- ▶ Do not modify the product.
- ▶ For service and repairs, please contact your national B. Braun/Aesculap agency.

### ⚠ CAUTION

Malfunction of the product due to damage in transport!

- ▶ Select suitable and safe packaging (ideally the original packaging).
- ▶ Keep the original packaging for possible returns in case of service.
- ▶ Package the product in such a way that the packaging does not become contaminated.

#### Service addresses

Aesculap Technischer Service

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78532 Tuttlingen / Germany

Phone: +49 7461 95-1601

Fax: +49 7461 16-2887

E-Mail: ats@aesculap.de

Other service addresses can be obtained from the address indicated above.

## 4.4 Accessories / spare parts

Art. no.	Designation
JG904	Disposable sterile cover, pack of 25
JF441R	Tray for 2D camera heads
JF441500	Spare parts set silicone for JF441R

## 5. Disposal

### ⚠ WARNING

Risk of infection due to contaminated products!

- ▶ Adhere to national regulations when disposing of or recycling the product, its components and its packaging.

#### Note

The user institution is obliged to reprocess the product before its disposal, see Reprocessing procedure.

The recycling pass can be downloaded from the Extranet as a PDF document under the respective article number. (The recycling pass includes disassembling instructions for the product, as well as information for proper disposal of components harmful to the environment.)

Products carrying this symbol are subject to separate collection of electrical and electronic devices. Within the European Union, disposal is taken care of by the manufacturer as a free-of-charge service.

- ▶ Detailed information concerning the disposal of the product is available through your national B. Braun/Aesculap agency, see Technical service.

## 6. Technical data

### 6.1 Classification acc. to Regulation (EU) 2017/745

Art. no.	Designation	Class
PV481	Full HD CMOS camera head with pendulum coupler	I
PV482	Full HD CMOS camera head with zoom coupler	I
PV485	Full HD 3CMOS camera head with zoom coupler	I

### 6.2 Performance data, information about standards

Protection rating	IPX7
Sensor format	Native Full HD 1/3"
Scanning system	Progressive scan
Frequency	50 Hz or 60 Hz
Protection class (according to IEC/DIN EN 60601-1)	I
Weight (excl. cable)	<ul style="list-style-type: none"> <li>■ PV481: 180 g</li> <li>■ PV482: 294 g</li> <li>■ PV485: 380 g</li> </ul>
Length of camera cable	3,5 m
Applied part	Type CF, defibrillation-proof (in combination with camera control unit PV480 or PV630)
Dimensions (L x W x H)	<ul style="list-style-type: none"> <li>■ PV481: 54 mm x 52 mm x 113,9 mm</li> <li>■ PV482: 149 mm x 52,5 mm x 50 mm</li> <li>■ PV485: 149 mm x 52,5 mm x 50 mm</li> </ul>
EMC	IEC/EN/DIN 60601-1-2
Conforming to standard	IEC/EN/DIN 60601-1
CISPR 11	In combination with PV480: Class B In combination with PV630: Class A

## 7. Symbols on product and packaging

	Indicates a hazard. If not avoided, the hazard can result in death or serious injury.
	Indicates a potential hazard. If not avoided, this hazard may result in injury and/or damage to the product.
	Follow the instructions for use
	Caution (IEC 60601-1 3rd edition) / Attention, please take note of the accompanying document (IEC 60601-1 2nd edition)
	Defibrillation-proof type CF rated part according to IEC 60601-1
	Permissible storage temperature
	Permissible relative air humidity during storage
	Permissible atmospheric pressure during storage
	Not permitted for use in a magnetic resonance environment
	Caution, fragile
	Article number
	Serial number
	Manufacturer
	Date of manufacture
<b>Rx only</b>	Federal Law (US) restricts this device to sale to or on behalf of a physician.
	Marking of electric and electronic devices according to directive 2012/19/EU (WEEE)
	Medical device

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Aesculap AG | Am Aesculap-Platz | 78532 Tuttlingen | Germany