

en	Instructions for use/Technical description Charnley Retractor
USA	Note for U.S. users This Instructions for Use is NOT intended for United States users. Please discard. The Instructions for Use for United States users can be obtained by visiting our website at www.aesculapusa.com and clicking the "Products" menu. If you wish to obtain a paper copy of the Instructions for Use, you may request one by contacting your local Aesculap representative or Aesculap's customer service at 1-800-282-9000. A paper copy will be provided to you upon request at no additional cost.
de	Gebrauchsanweisung/Technische Beschreibung Charnley Retraktor
fr	Mode d'emploi/Description technique Rétracteur Charnley
es	Instrucciones de manejo/Descripción técnica Retractor Charnley
it	Istruzioni per l'uso/Descrizione tecnica Retratore Charnley
pt	Instruções de utilização/Descrição técnica Retrator Charnley
nl	Gebruiksaanwijzing/Technische beschrijving Charnley-retractor
sv	Bruksanvisning/Teknisk beskrivning Charnley Retraktor
ru	Инструкция по применению/Техническое описание Ретрактор Чарнли
cs	Návod k použití/Technický popis Retraktor Charnley
pl	Instrukcja użytkowania/Opis techniczny Retraktor Charnley
sk	Návod na použitie/Technický opis Charnleyho retractor
tr	Kullanım Kılavuzu/Teknik açıklama Charnley rekraktörü

Aesculap®

Charnley Retractor

Legend

- 1 Retractor arm for fixed valve
- 2 Retractor arm for moving valve
- 3 Moving valve
- 4 Fixed valve
- 5 Locking bolt with spring plate
- 6 Weight with chain

1. About this document

Note
General risk factors associated with surgical procedures are not described in these instructions for use.

1.1 Scope

These instructions for use apply for the following products:

Art. no.	Designation
MB191R	Charnley Basic Retractor Kit
MG697R	Frame Charnley 300 mm x 240 mm
MG081R	Weight Charnley 1 740 g with chain
MB156R, MB157R, MB177R, MB178R, MB192R, MB193R, MG565R, MG566R, MG567R, MG568R, MG569R, MG570R, MG698R, MG699R, MG700R, MG701R, MG702R, MG703R	Sheet Charnley
► For item-specific instructions for use and information on material compatibility, see also the Aesculap Extranet at https://extranet.bbraun.com	

1.2 Safety messages

Safety messages make clear the dangers to patient, user and/or product that could arise during the use of the product. Safety messages are labeled as follows:

⚠ WARNING
Indicates a potential danger. If these are not prevented, minor or minor injury or damage to the product may result.

⚠ CAUTION
Indicates a possible threat of material damage. If not avoided, the product may be damaged.

2. Clinical use

2.1 Areas of use and limitations of use

2.1.1 Intended use

The Charnley Retractor is a retractor system used for the retraction of soft parts during hip replacement (arthroplasty of the hip). It can be used in both standard and minimal or less invasive surgery. The aim is to provide an unobstructed view of the hip joint.

2.1.2 Indications

Note
The manufacturer is not responsible for any use of the product against the specified indications and/or the described applications.

Application for the retraction of soft tissue in arthroplasty of the hip for both standard and minimal or less invasive surgery.

2.1.3 Absolute contraindications

No known absolute contraindications.

2.1.4 Relative contraindications

The following conditions, individual or combined, can lead to delayed healing or compromise the success of the operation:

- Medical or surgical conditions (e.g. comorbidities) which could hinder the success of the operation.

In the presence of relative contraindications, the user decides individually regarding the use of the product.

2.2 Safety information

2.2.1 Clinical user

General safety instructions

To prevent damage due to improper delivery and application and to avoid jeopardizing warranty and liability:

- Use product only in accordance with these instructions for use.
- Observe safety information and maintenance instructions.
- Only combine Aesculap products with each other.
- Only persons who have the necessary training, knowledge and experience may operate and apply the product and accessories.
- Store brand new or unused product in a dry, clean and protected place.
- Check the functionality and proper condition before using the product.
- Keep instructions for use accessible to the user.
- Always adhere to applicable standards.

Note
The user agrees to report all serious incidents occurring in connection with the product to the manufacturer and to the competent authority of the country in which the user is registered.

Notes on surgical interventions

The health care professional shall make the decision on the specific applicability based on the warranted properties and technical data.

- The user shall be responsible for the proper execution of the surgical intervention.
- The successful use of the product requires appropriate clinical training in addition to theoretical and practical mastery of all required surgical techniques, including the use of this product.
- The user is required to obtain information from the manufacturer if there is an ambiguous preoperative situation regarding the use of the product.

2.2.2 Sterility

The product is delivered in an unsterile condition.

- Clean the new product after removing its transport packaging and prior to its initial sterilization.

2.3 Application

⚠ WARNING
Risk of injury and/or malfunction!
► Prior to each use, inspect the product for loose, bent, broken, cracked, worn, or fractured components.
► Always carry out a function test prior to each use of the product.

3. Validated reprocessing procedure

3.1 General safety information

Note
Adhere to national statutory regulations, national and international standards and directives, and local, clinical hygiene instructions for sterile processing.

Note
For patients with Creutzfeldt-Jakob disease (CJD), suspected CJD or possible variants of CJD, observe the relevant national regulations concerning the reprocessing of products.

Note
Mechanical reprocessing should be favored over manual cleaning as it gives better and more reliable results.

Note
Successful processing of this medical device can only be ensured if the processing method is first validated. The operator/sterile processing technician is responsible for this.

Note
If there is no final sterilization, then a virucidal disinfectant must be used.

Note
For up-to-date information about reprocessing and material compatibility, see also the Aesculap Extranet at <https://extranet.bbraun.com>
The validated steam sterilization procedure was carried out in the Aesculap sterile container system.

3.2 General information

Dried or affixed surgical residues can make cleaning more difficult or ineffective and lead to corrosion. Therefore the time interval between application and processing should not exceed 6 h; also, neither fixating pre-cleaning temperatures >45 °C nor fixating disinfecting agents (active ingredient: aldehydes/alcohols) should be used. Excessive measures of neutralizing agents or basic cleaners may result in a chemical attack and/or to fading and the laser marking becoming unreadable visually or by machine for stainless steel.

Residues containing chlorine or chlorides e.g. in surgical residues, medicines, saline solutions and in the service water used for cleaning, disinfection and sterilization will cause corrosion damage (pitting, stress corrosion) and result in the destruction of stainless steel products. These must be removed by rinsing thoroughly with demineralized water and then drying.

Additional drying, if necessary.
Only process chemicals that have been tested and approved (e.g. VAH or FDA approval or CE mark) and which are compatible with the product's materials according to the chemical manufacturers' recommendations may be used for processing the product. All the chemical manufacturer's application specifications must be strictly observed. Failure to do so can result in the following problems:

- Optical changes of materials, e.g. fading or discoloration of titanium or aluminum. For aluminum, the application/process solution only needs to be of pH >8 to cause visible surface changes.
- Material damage such as corrosion, cracks, fracturing, premature aging or swelling.
- Do not use metal cleaning brushes or other abrasives that would damage the product surfaces and could cause corrosion.
- Further detailed advice on hygienically safe and material-/value-preserving reprocessing can be found at www.ak-lorg, link to Publications, Red Brochure – Proper maintenance of instruments.

3.3 Restriction of re-usability

- There is no maximum number of applications and processing cycles for the product.
- The service life of the product is limited by damage, normal wear, type and duration of the application, handling, storage and transportation of the product.
- Careful visual and functional testing prior to next use is the best way to identify a malfunctioning product.

3.4 Preparations at the place of use

- If applicable, rinse non-visible surfaces preferably with deionized water, with a disposable syringe for example.
- Remove any visible surgical residues to the extent possible with a damp, lint-free cloth.
- Transport the dry product in a sealed waste container for cleaning and disinfection within 6 hours.

3.5 Preparing for cleaning

- Disassemble the product prior to cleaning, see Disassembly.

3.6 Disassembly

- If applicable: Remove weight with chain 6 from moving valve 3.
- Release moving valve 3 from retractor arm 2.
- Push the locking bolt with spring plate 5 using thumb in the direction of the arrow and slide the fixed valve 4 from retractor arm 1 using other thumb, see Fig. A.

3.7 Cleaning/Disinfection

3.7.1 Product-specific safety information on the reprocessing method

Damage to or destruction of the product due to inappropriate cleaning/disinfecting agents and/or excessive temperatures!

- Use cleaning agents and disinfectants approved for stainless steel according to the manufacturer's instructions.
- Observe specifications regarding concentration, temperature and exposure time.
- Do not exceed a disinfection temperature of 96 °C.
- Use suitable cleaning/disinfecting agents if the product is disposed of in a wet condition. To prevent foaming and degradation of the efficacy of the process chemicals: prior to mechanical cleaning and disinfection, rinse the product thoroughly with running water.

3.7.2 Validated cleaning and disinfection procedure

Validated procedure	Specific requirements	Reference
Manual cleaning with immersion disinfection	<ul style="list-style-type: none">■ Suitable cleaning brush■ Drying phase: Use a lint-free cloth or medical compressed air	Chapter Manual cleaning/disinfection and subsection: <ul style="list-style-type: none">■ Chapter Manual cleaning with immersion disinfection
Mechanical alkaline cleaning and thermal disinfection	<ul style="list-style-type: none">■ Place the product on a tray that is suitable for cleaning (avoid rinsing blind spots).	Chapter Mechanical cleaning/disinfection and subsection: <ul style="list-style-type: none">■ Chapter Mechanical alkaline cleaning and thermal disinfection

3.8 Manual cleaning/disinfection

- Prior to manual disinfecting, allow water to drip off for a sufficient length of time to prevent dilution of the disinfecting solution.
- After manual cleaning/disinfection, check visible surfaces visually for residues.
- Repeat the cleaning/disinfection process if necessary.

3.8.1 Manual cleaning with immersion disinfection

Phase	Step	T [°C/°F]	t [min]	Conc. [%]	Water quality	Chemical
I	Disinfecting clean- ing	RT (cold)	>15	2	D–W	Aldehyde-free, phenol-free, and QUAT-free concentrate, pH ~ 9*
II	Intermediate rinse	RT (cold)	1	–	D–W	–
III	Disinfection	RT (cold)	5	2	D–W	Aldehyde-free, phenol-free, and QUAT-free concentrate, pH ~ 9*
IV	Final rinse	RT (cold)	1	–	FD–W	–
V	Drying	RT	–	–	–	–

D–W: Drinking water
FD–W: Fully desalinated water (demineralized, low microbiological contamination: drinking water quality at least)
RT: Room temperature
*Recommended: BBraun Stabimed fresh

► Note the information on appropriate cleaning brushes and disposable syringes, see Validated cleaning and disinfection procedure.

Phase I
► Fully immerse the product in the cleaning/disinfectant for at least 15 min. Ensure that all accessible surfaces are moistened.
► Clean the product with a suitable cleaning brush in the solution until all discernible residues have been removed from the surface.
► If applicable, brush through non-visible surfaces with an appropriate cleaning brush for at least 1 min.
► Mobilize non-rigid components, such as set screws, links, etc. during cleaning.
► Thoroughly rinse through these components with the cleaning disinfectant solution (at least five times), using a disposable syringe.

Phase II
► Rinse/flush the product thoroughly (all accessible surfaces) under running water.
► Mobilize non-rigid components, such as set screws, joints, etc. during rinsing.
► Drain any remaining water fully.

Phase III
► Fully immerse the product in the disinfectant solution.
► Mobilize non-rigid components, such as set screws, joints, etc. during rinsing.
► Rinse lumens at least 5 times at the beginning of the exposure time using an appropriate disposable syringe. Ensure that all accessible surfaces are moistened.

Phase IV
► Rinse/flush the product thoroughly (all accessible surfaces).
► Mobilize non-rigid components, such as set screws, joints, etc. during final rinse.
► Rinse lumens with an appropriate disposable syringe at least five times.
► Drain any remaining water fully.

Phase V
► Dry the product in the drying phase with suitable equipment (e.g. cloth, compressed air), see Validated cleaning and disinfection procedure.

3.9 Mechanical cleaning/disinfection

Note
The cleaning and disinfection device must be of tested and approved effectiveness (e.g. FDA approval or CE mark according to DIN EN ISO 15883).

Note
The cleaning and disinfection device used for processing must be serviced and checked at regular intervals.

3.9.1 Mechanical alkaline cleaning and thermal disinfection

Machine type: single-chamber cleaning/disinfection device without ultrasound

Phase	Step	T [°C/°F]	t [min]	Water quality	Chemical/Note
I	Pre-rinse	<25/77	3	D–W	–
II	Cleaning	55/131	10	FD–W	<div><div>■</div>Concentrate, alkaline:<div><div>–</div>pH = 13<div>–</div><5 % anionic surfactant</div><div><div>■</div>0.5 % working solution<div>–</div>pH = 11*</div></div>
III	Intermediate rinse	>10/50	1	FD–W	–
IV	Thermal disinfecting	90/194	5	FD–W	–
V	Drying	–	–	–	According to the program for cleaning and disinfection device

D–W: Drinking water
FD–W: Fully desalinated water (demineralized, low microbiological contamination: drinking water quality at least)
*Recommended: BBraun Helimatic Cleaner alkaline

► Check visible surfaces for residues after mechanical cleaning/disinfecting.

3.10 Inspection

► Allow the product to cool down to room temperature.
► Dry the product if it is wet or damp.

3.10.1 Visual inspection

► Ensure that all soiling has been removed. In particular, pay attention to mating surfaces, hinges, shafts, recessed areas, drill grooves and the sides of the teeth on rasps.
► If the product is dirty: repeat the cleaning and disinfection process.
► Check the product for damage, e.g. insulation or corroded, loose, bent, broken, cracked, worn or severely scratched and fractured components.
► Check the product for missing or faded labels.
► Check the surfaces for rough spots.
► Check the product for burrs that could damage tissue or surgical gloves.
► Check the product for loose or missing parts.
► Immediately put aside damaged or inoperative products and send them to Aesculap Technical Service, see Technical service.

3.10.2 Functional test

► Assemble disassembled products, see Assembly.
► Check that the product functions correctly.
► Check that all moving parts are working property (e.g. hinges, locks/latches, sliding parts etc.).
► Check rotating products (e.g. reusable drills and cutters) for bends and deformities. To do this, roll the product on an even surface.
► Check for compatibility with associated products.
► Immediately put aside inoperative products and send them to Aesculap Technical Service, see Technical service.

3.11 Assembly

► Push fixed valve 4 onto retractor arm 1 until locking bolt snaps into spring plate 5.
► Stick moving valve 3 on bolt of retractor arm 2.
► If necessary: Hang weight with chain 6 on moving valve 3.

3.12 Packaging

► Place the product in its holder or on a suitable tray. Ensure that sharp edges are covered.
► Package trays appropriately for the sterilization process (e.g. in Aesculap sterile containers).
► Ensure that the packaging provides sufficient protection against contamination of the product during storage.

3.13 Steam sterilization

► Check to ensure that the sterilizing agent will come into contact with all external and internal surfaces (e.g., by opening any valves and faucets).
► Validated sterilization process

–

Steam sterilization in fractionated vacuum process

–

Steam sterilizer in accordance with DIN EN 285 and validated in accordance with DIN EN ISO 17665

–

Sterilization in fractionated vacuum process at 134 °C, holding time 5 min.

► If several devices are sterilized at the same time in the same steam sterilizer: Ensure that the maximum permitted load according to the manufacturers' specifications is not exceeded.

3.14 Storage

► Store sterile products in germ-proof packaging, protected from dust, in a dry, dark, temperature-controlled area.

4. Technical service

 **CAUTION**
Modifications carried out on medical technical equipment may result in loss of guarantees/warranty rights and forfeiture of applicable licenses.
► Do not modify the product.
► For service and repairs, please contact your national B. Braun/Aesculap agency.

Service addresses
Aesculap Technischer Service
Am Aesculap-Platz
78532 Tuttlingen / Germany
Phone: +49 7461 95-1601
Fax: +49 7461 16-2887
E-Mail: ats@aesculap.de
Other service addresses can be obtained from the address indicated above.

5. Disposal

 **WARNING**
Risk of infection due to contaminated products!
► Adhere to national regulations when disposing of or recycling the product, its components and its packaging.

Note
The user institution is obliged to reprocess the product before its disposal, see Validated reprocessing procedure.

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