# Aesculap®

## Aesculap Endoscopic Technology

### en Instructions for use/Technical description

LED light source OP950

### USA Note for U.S. users

This Instructions for Use is NOT intended for United States users. Please discard. The Instructions for Use for United States users can be obtained by visiting our website at [www.aesculapusa.com](http://www.aesculapusa.com). If you wish to obtain a paper copy of the Instructions for Use, you may request one by contacting your local Aesculap representative or Aesculap's customer service at 1-800-282-9000. A paper copy will be provided to you upon request at no additional cost.

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


### Contents

1. About this document.....	3	2.8 Application .....	6
1.1 Scope .....	3	2.8.1 Turning the device on and off .....	6
1.2 Safety messages.....	3	2.8.2 Brightness control.....	6
2. Clinical use.....	3	2.8.3 Activation/standby mod .....	6
2.1 Product description.....	3	2.9 Ending the operation of the device.....	6
2.1.1 Scope of supply .....	3	2.10 Troubleshooting list .....	7
2.1.2 Components required for operation.....	3	2.10.1 Error with error code .....	7
2.1.3 Operating principle .....	3	2.10.2 Malfunctions.....	7
2.2 Areas of use and limitations of use .....	4	2.10.3 Replacing fuses.....	7
2.2.1 Intended use.....	4	3. Validated reprocessing procedure .....	8
2.2.2 Indications .....	4	3.1 General safety instructions .....	8
2.2.3 Contraindications.....	4	3.2 General information .....	8
2.3 Safety information.....	4	3.3 Reusable products.....	8
2.3.1 Clinical user .....	4	3.4 Cleaning/Disinfection .....	8
2.3.2 Product.....	4	3.4.1 Product-specific safety information on the reprocessing method.....	8
2.3.3 Sterility.....	4	3.5 Wipe disinfection .....	8
2.4 Preparation.....	4	3.6 Inspection .....	9
2.4.1 Set-up environment/set-up location .....	5	3.6.1 Visual inspection.....	9
2.4.2 Stacking of units .....	5	3.6.2 Functional test .....	9
2.5 System set-up .....	5	3.7 Storage .....	9
2.5.1 Connecting devices for communication.....	5	3.7.1 Ambient conditions.....	9
2.5.2 Connecting the optical cable .....	5	4. Maintenance and service .....	9
2.5.3 Separating the LED light source from the optical cable.....	5	4.1 Maintenance .....	9
2.5.4 Connecting the sterile adapter (optional) .....	5	4.2 Technical service.....	9
2.5.5 Removing the sterile adapter from the LED light source .....	5	4.3 Accessories / spare parts .....	9
2.5.6 Connecting the power supply.....	5	4.3.1 Accessories for sterile adapter .....	9
2.5.7 Connecting the accessories.....	6	4.3.2 Accessories for power cord .....	9
2.5.8 Putting out of operation.....	6	4.3.3 Spare part for fuse .....	10
2.6 Functional test.....	6	4.3.4 Communications accessories .....	10
2.7 Optical cable test.....	6	4.3.5 Light lead accessories .....	10
		4.3.6 Equipotential bonding cable accessories .....	10
		5. Disposal .....	10
		6. Technical data .....	10
		6.1 Classification acc. to Regulation (EU) 2017/745 .....	10
		6.2 Performance data, information about standards .....	10
		7. Symbols on product and packaging .....	11

# 1. About this document

## Note

*General risk factors associated with surgical procedures are not described in these instructions for use.*

## 1.1 Scope

These instructions for use apply for the following products:

Art. no.	Designation
OP950	LED light source

- For article-specific instructions for use as well as information on material compatibility and lifetime see B. Braun eifu at [eifu.bbraun.com](http://eifu.bbraun.com)

## 1.2 Safety messages

Safety messages make clear the dangers to patient, user and/or product that could arise during the use of the product. Safety messages are labeled as follows:

### ⚠ DANGER

Indicates a possible threat of danger. If not avoided, death or serious injury may result.

### ⚠ WARNING

Indicates a possible threat of danger. If not avoided, minor or moderate injury may result.

### ⚠ CAUTION

Indicates a possible threat of material damage. If not avoided, the product may be damaged.

# 2. Clinical use

## 2.1 Product description

### 2.1.1 Scope of supply

Designation	Art. no.
LED light source	OP950
Instructions for use	TA014627

### 2.1.2 Components required for operation

Designation	Art. no.
LED light source	OP950
Optical cable with a diameter of 4.8 mm of type Aesculap/Storz, Olympus or Wolf	-
Power cord	-

### 2.1.3 Operating principle

The LED light source OP950 is very similar to a point light source.

The light is propagated through the optical cable and directed to the field of view (operative field).

### Brightness control

The brightness is adjusted gradually (5 % increments) via the brightness control 4.

If a camera system with auto light quantity control is connected via the communication interface cable OP942, the light intensity can be controlled via the camera.

When the brightness control 4 is adjusted during active auto light quantity control, the system is immediately switched to manual mode and the light intensity is displayed in %.

### Activation/standby mode

On pressing the brightness control with push function 4, the active light source is put into standby mode. On pressing the brightness control with push function 4 again, the light source goes back to the previously set value.

If a camera system with auto light quantity control is connected via the communication interface cable OP942, the LED light source can also be activated or deactivated via the camera.

### Multi optical cable socket

Various optical fiber cables from Aesculap/Storz, Olympus and Wolf with an active diameter of 3.5 mm or 4.8 mm can be plugged into the multi-optical cable socket 7.


### Light source test unit with sterile adapter

The optical cable test unit can be used in the sterile field with a sterile adapter OP941.

The optical cable test unit is used for testing the light transmission of the 4.8 mm optical fiber cable connected in each case.

The testing of 3.5 mm optical fiber cables causes false measurement results.

### Modes shown on the display

Display	Explanation
STANDBY	Standby Device is switched on, LED is off. Circular green light on the brightness control is deactivated.
100 %	Light intensity % (0 % to 100 %)
AUTO	Automatic light quantity control is active Is only displayed if a corresponding camera system is automatically controlling the light quantity.
	Optical cable is not plugged in or not plugged in properly
	Optical cable is being tested
	Result of optical cable test: Optical cable OK.
	Result of optical cable test: Optical cable should be replaced.

## 2.2 Areas of use and limitations of use

### 2.2.1 Intended use

The LED light source OP950 is intended for use in endoscopy systems for human medicine. In combination with a light cable, endoscope, camera and monitor, its purpose is to illuminate the inside of the human body.

### 2.2.2 Indications

There are no medical indications associated with the use of the LED light source OP950 itself if separated from the endoscopic system.

The indications for an endoscopic application depend on the condition of the patient and the individual risk/benefit analysis by the surgeon.

### 2.2.3 Contraindications

There are no medical contraindications associated with use of the LED light source OP950 itself if separated from the endoscopic system.

Depending on the intended use of the applied part, contraindications may exist in respect of the patient's general condition or the specific presenting disease.

## 2.3 Safety information

### 2.3.1 Clinical user

#### General safety information

To prevent damage caused by improper setup or operation, and to not compromise the manufacturer warranty and liability:

- ▶ Use the product only according to these instructions for use.
- ▶ Follow the safety and maintenance instructions.
- ▶ Ensure that the product and its accessories are operated and used only by persons with the requisite training, knowledge and experience.
- ▶ Store any new or unused products in a dry, clean, and safe place.
- ▶ Prior to use, check that the product is in good working order.
- ▶ Keep the instructions for use accessible for the user.

#### Note

*The user is obligated to report all severe events in connection with the product to the manufacturer and the responsible authorities of the state in which the user is located.*

#### Notes on surgical procedures

It is the user's responsibility to ensure that the surgical procedure is performed correctly.

Appropriate clinical training as well as a theoretical and practical proficiency of all the required operating techniques, including the use of this product, are prerequisites for the successful use of this product.

The user is required to obtain information from the manufacturer if there is an unclear preoperative situation regarding the use of the product.

### 2.3.2 Product

#### Product-specific safety information

Risk of death by electric shock!

- ▶ Do not open the product.
- ▶ Connect the product only to a grounded power supply.

Risk of injury from current leakage when touching the patient!

- ▶ Never touch the patient and the device housing or accessible contacts of electric devices at the same time.
- ▶ Never touch the patient and non-medical electric devices at the same time.

Possibly hazardous optical radiation emitted from this product!

- ▶ Do not stare at operating lamp. May be harmful for the eye.

Device failure due to incorrect storage and usage conditions!

- ▶ Store and operate the product only within the specified environmental conditions.
- ▶ Observe "Notes on electromagnetic compatibility (EMC) for LED light source OP950" TA022463, see B. Braun eIFU at [eifu.bbraun.com](http://eifu.bbraun.com)
- ▶ Only combine Aesculap products with each other (except for optical fiber cable).
- ▶ Always adhere to applicable standards.
- ▶ In the sterile field only use sterile processed sterile adapter OP941 (optional).
- ▶ Only allow products and accessories to be operated and used by individuals who have been trained in the use of endoscopic procedures.
- ▶ Ensure that the electrical facilities in the room in which the device is being used conform to IEC norms.
- ▶ Unplug the device by pulling on the plug, and never on the power cord.
- ▶ Do not use LED light sources in areas where there is a risk of explosion.
- ▶ Make certain that the ambient temperature does not exceed 40 °C.
- ▶ Make certain that the device's ventilation ducts are not covered.

#### Ambient conditions

The following environmental conditions apply to the use of the product:

Temperature	10 °C to 40 °C
Relative humidity	30 % to 90 %
Atmospheric pressure	700 hPa to 1 060 hPa

### 2.3.3 Sterility

The product is supplied non-sterile and must be used under non-sterile conditions.

- ▶ Inspect the new product after removing its transport packaging and prior to first use to ensure it is in good working order.

## 2.4 Preparation

Non-compliance with these rules will result in complete exclusion of liability on the part of Aesculap.

- ▶ When setting up and operating the product, adhere to:
  - national regulations for installation and operation,
  - national regulations on fire and explosion protection.

#### Note

*For the safety of patients and users it is essential that the mains power cord and, especially, the protective ground connection are intact. In many cases defective or missing protective ground connections are not recognized immediately.*

- ▶ Connect the device via the potential equalization terminal at the rear panel of the device to the potential equalization system of the room used for medical purposes (see IEC 60601-1/EN 60601-1 or the relevant national standards).
- ▶ Ensure that all devices operated in the vicinity meet their relevant EMC requirements.

#### Note

*The potential equalization lead can be ordered from the manufacturer as art. no. GK535 (4 m long) or TA008205 (0.8 m long).*

#### 2.4.1 Set-up environment/set-up location

##### **⚠ DANGER**

**Fire and explosion hazard!**

- ▶ Use the product outside of explosion hazard zones (e.g. areas with highly pure oxygen or anesthetic gases).
- ▶ Make certain that the power plug is connected to the power supply outside of any areas where there is a risk of explosion.
- ▶ Have ready replacement products/light source.

##### **⚠ WARNING**

**Risk of injury and/or product malfunction due to incorrect operation of the electromedical system!**

- ▶ Adhere to the instructions for use of any medical device.
- ▶ Do not use the product in a magnetic resonance environment.

The device is authorized for operation in the operating room.

##### **Note**

*After being installed and put into operation, the device must not be transported or moved to a different set-up location.*

##### **Note**

*If the equipment cart is set up outside of the sterile area, measure the cable length to give enough room for movement so that the sterility is not compromised.*

- ▶ Make certain that the ventilation slots in the bottom of the housing and in the rear panel of the device are not covered, e.g. by an OR cloth.
- ▶ Make certain that the following minimum distances from other devices are maintained:
  - on the sides 50 mm
  - from the rear 60 mm
- ▶ Make certain that the control elements, mains power switch and power socket are freely accessible for the user.
- ▶ Make certain the system is set up on a sufficiently stable support (e.g. a table, ceiling support, equipment cart, etc.).
- ▶ Follow the instructions for use of the support.

#### 2.4.2 Stacking of units

- ▶ Do not exceed the maximum stack height of 450 mm.
- ▶ Position units so that they are stable.
- ▶ Stack Aesculap devices evenly on top of each other.
- ▶ Never stack in an offset arrangement.

### 2.5 System set-up

##### **⚠ WARNING**

**Risk of burns from excessively high operating temperatures!**

**The optical cable socket and optical cable connections get hot.**

- ▶ Do not situate optical cables near any flammable objects (e.g. drape).
- ▶ Select the lowest lighting level necessary to provide satisfactory illumination of the target area.

##### **⚠ WARNING**

**Danger of blinding from light rays!**

- ▶ Never look into the open end of an optical cable or endoscope that is connected.

##### **⚠ CAUTION**

**Risk of injury if the light source fails!**

- ▶ Have a ready-to-operate spare light source ready or use a light source that has a spare bulb.

#### 2.5.1 Connecting devices for communication

The communication interfaces 8 connect the light source to a superordinate device. Through these, the light source receives control signals from the connected device.

The connection of a superordinate device to the communication interfaces leads to the construction of a system.

##### **⚠ WARNING**

**Risk of injury due to unapproved configuration using additional components!**

- ▶ Make sure that accessory equipment and/or peripheral devices that are connected to the device interfaces verifiably satisfy the IEC 60601-1.
- ▶ Make sure that all configurations comply with the system standard IEC 60601-1 chapter 16.
- ▶ Do not bend the communications interface cable, to prevent the cable from breaking and the connection from failing.

#### 2.5.2 Connecting the optical cable

##### **⚠ WARNING**

**Risk of injury if the optical cable slips out of place!**

- ▶ Ensure the optical cable is sufficiently fixed in the direct operating field.
- ▶ Secure the strain relief cord grip on the optical cable.
- ▶ Plug the optical cable into the optical cable socket 7 until it audibly snaps into place.
- ▶ Do not crush or crimp the optical cable or wind it too tightly (minimum diameter 150 mm).

#### 2.5.3 Separating the LED light source from the optical cable

- ▶ To unlock the optical cable, turn the activator 6 in the marked direction.

#### 2.5.4 Connecting the sterile adapter (optional)

- ▶ Before the optical cable test, screw sterile adapter OP941 onto the optical cable test unit 5.
- ▶ If necessary, screw on the adapter for connecting the optical cable.

#### 2.5.5 Removing the sterile adapter from the LED light source

- ▶ After the optical cable test, unscrew sterile adapter OP941 from the optical cable test unit 5.

#### 2.5.6 Connecting the power supply

##### **⚠ DANGER**

**Risk of death by electric shock!**

- ▶ Connect the product only to a grounded power supply.

##### **Note**

*The connection of electrical devices to a multiple socket creates a system and can lead to a reduced level of safety.*

- ▶ Plug in the power cord into the power cord connector 9 on the rear of the LED light source.
- ▶ Plug in the mains plug at a socket of the building mains.
- ▶ Turn on the LED light source with the power ON 2 switch.  
The power ON 3 indicator lamp shows that the device is activated.
- ▶ When operating multiple devices concomitantly, connect the devices to the potential equalization terminal.

## 2.5.7 Connecting the accessories

### **⚠ DANGER**

Risk of injury due to unapproved configuration using additional components!

- For all applied components, ensure that their classification matches that of the application component (e.g. Type BF or Type CF) of the respective device.

Combinations of accessories that are not mentioned in the present instructions for use may only be employed if they are specifically intended for the respective application, and if they do not compromise the performance and safety characteristics of the products.

Any equipment connected at the interfaces must demonstrably meet the respective IEC standards (e.g. IEC 60950 for data processing equipment and IEC/DIN EN 60601-1 for medical electrical equipment).

All configurations must fulfill the basic standard IEC/DIN EN 60601-1. Any individual connecting devices with one another is responsible for such configuration and must ensure compliance with basic standard IEC/DIN EN 60601-1 or applicable national standards.

- In case of questions, contact your B. Braun/Aesculap partner or Aesculap Technical Service, address see Technical service.

## 2.5.8 Putting out of operation

### *Note*

The safe and all-pole disconnection of the product from the main power supply is only guaranteed when the power cord is unplugged.

- Switching off the device: Activate the power OFF switch 1.
- Unplug the power cord from the power socket.

The operation of the device is safely terminated.

## 2.6 Functional test

- Prior to each use, perform a functionality test for the LED light source.
- Make certain that the self-test runs in accordance with legal regulations.

If the self-test is successful, the word "STANDBY" appears after the start logo. After activating the brightness control, the request appears to plug in the optical fiber cable.

In case of an error, an error code appears in the display, see Error with error code.

If the device is functioning properly, the LED light source is activated when the optical fiber cable is plugged in and the brightness control is activated. After the device is switched on, the power ON 2 indicator lamp lights up on the front of the device.

## 2.7 Optical cable test

If necessary, carry out an optical cable test:

- Turn on the LED light source with the power ON 3 switch.
- Plug the optical cable into the optical cable socket 7 until it audibly snaps into place.
- Screw sterile adapter OP941 onto the optical cable test unit 5.

### *Note*

The sterile adapter OP941 consists of three parts. Of these, two bushes serve the connection of the different optical fiber connections.

- According to the optical fiber connection, unscrew or screw off the bushes of the sterile adapter OP941. Keep the bushes.
- Screw the distal end of the optical cable to be tested onto the sterile adapter.
- Activate the brightness control 4 to switch on the light.
- Press the optical fiber unit 5 against the light source.

The results are displayed in the display.

Display	Explanation
	Optical cable can still be used. (100 % - 41 %)
	Optical cable should be replaced. (40 % - 0 %)

### *Note*

The optical cable test takes place regardless of the selected light strength (e.g. 80 %).

## 2.8 Application

### **⚠ WARNING**

Risk of injury and/or malfunction!

- Always carry out a function check prior to each use of the product.

### **⚠ WARNING**

Risk of injury due to application of the product outside the field of view!

- Apply the product only under visual control.

### 2.8.1 Turning the device on and off

- Turn on the LED light source with the power ON 3 switch.  
The power ON indicator 2 light is on.
- Turn off the LED light source with the power OFF 1 switch.

### 2.8.2 Brightness control

- To adjust the brightness in increments (5 % increments), turn the brightness control 4.  
Turning it in a clockwise direction increases the light intensity.

### *Note*

The brightness can also be automatically controlled using a corresponding camera system.

### 2.8.3 Activation/standby mod

- To activate the LED light source, press the brightness control 4.
- To put the active LED light source (optical fiber cable is plugged in, light is emitted) in standby mode, press the brightness control 4.
- To take it out of Standby mode and put it back to the previously set value, press the brightness control 4 again.

### *Note*

The LED light source can also be activated/deactivated using a corresponding camera system.

## 2.9 Ending the operation of the device

- To end the operation of the device, press the power OFF switch 1.
- Remove the optical cable from the sterile adapter 5.
- To completely disconnect the device from the power supply, remove the plug from the power cord connector 10 on the rear of the device.

## 2.10 Troubleshooting list

### 2.10.1 Error with error code

Error code in display	Cause	Remedy
Error Code 1	LED light intensity cannot be adjusted or LED is defective.	Replace the LED light source as soon as possible with a replacement.
Error Code 2	Failure of housing fan	Replace the LED light source as soon as possible with a replacement.
Error Code 3	Housing temperature is more than 60 °C. The light output is limited to 70 % light output. The device will not switch off, but is not guaranteed to remain fully functional.	Replace the LED light source as soon as possible with a replacement.
Error Code 4	The device is fitted with two internal fans for cooling the LEDs. The first fan is defective.	Replace the LED light source as soon as possible with a replacement.
Error Code 5	The device is fitted with two internal fans for cooling the LEDs. The second fan is defective.	Replace the LED light source as soon as possible with a replacement.

### 2.10.2 Malfunctions

Malfunction	Detection	Cause	Remedy
LED light source and fan do not work	Device not powered	Plug not inserted in the socket	Insert plug in socket
	Indicators not illuminated	Fuses blown	Replace fuses, see Replacing fuses
	Indicator light not illuminated even after change of fuse	Malfunction in power supply unit	Have LED light source repaired by manufacturer
Fan does not work	No air comes out of the air outlet	Malfunction in fan	Have LED light source repaired by manufacturer
Brightness significantly lower	The optical cable cannot be removed	Optical cable not completely inserted	Use the correct optical cable connection
Light source cannot be controlled	Light source does not respond to signals from the camera	Interface cable not correctly plugged in or defective	Check connection or replace cable

### 2.10.3 Replacing fuses

**Risk of death by electric shock!**

- Unplug the device before changing the fuses.

Ask your B. Braun/Aesculap agency about the prescribed fuse set.

- Use a small screwdriver to release the clip on the fuse holder 9.
- Remove fuse holder 9.
- Replace both fuses.
- Reinsert fuse holder 9 so that it audibly snaps into place.

#### Note

If the fuses burn out frequently, the device is faulty and should be repaired, see Technical service.

### 3. Validated reprocessing procedure

#### 3.1 General safety instructions

**Note**

Adhere to national statutory regulations, national and international standards and directives, and local, clinical hygiene instructions for sterile processing.

**Note**

Successful processing of this medical device can only be ensured if the processing method is first validated. The operator/sterile processing technician is responsible for this.

**Note**

For up-to-date information about reprocessing and material compatibility, see also the Aesculap Extranet at <https://extranet.bbraun.com>

#### 3.2 General information

**Note**

The housing of the LED light source is only suitable for wiping with a moist cloth (wipe disinfection). Do not immerse into liquid.

#### 3.3 Reusable products

Influences of the reprocessing which lead to damage to the product are not known.

A careful visual and functional inspection before the next use is the best opportunity to recognize a product that is no longer functional, see Inspection.

#### 3.4 Cleaning/Disinfection

##### 3.4.1 Product-specific safety information on the reprocessing method

Risk of electric shock and fire hazard!

- Unplug the device before cleaning.
- Do not use flammable or explosive cleaning or disinfecting solutions.
- Ensure that no fluids will penetrate the product.

Risk of electric shock and damage to the device from moisture!

- Carry out the wipe disinfection with damp cloths only. Never use dripping wet cloths.
- Remove excess cleaning solution from the device immediately.
- Only reconnect the LED light source to power when all cleaned parts are completely dry.

Damage to the product due to improper reprocessing!

- Use only wipe disinfection for cleaning and disinfecting the LED light source.
- Do not, under any circumstances, clean or disinfect the product in an ultrasonic cleaning bath.
- Use cleaning and disinfecting agents which are suitable and approved for the present product.
- Observe the manufacturer's cleaning and disinfecting instructions regarding concentration, temperature and exposure time.
- Do not sterilize the product under any circumstances.

#### 3.5 Wipe disinfection

Phase	Step	T [°C/°F]	t [min]	Conc. [%]	Water quality	Chemical
I	Cleaning	RT	1	-	-	Alcohol(s), quaternary compound(s) *
II	Wipe disinfection	RT	≥1	-	-	Alcohol(s), quaternary compound(s) *

RT: Room temperature

\* Recommended: Meliseptol® wipes sensitive (B. Braun)

**Phase I**

- Remove any visible residues with a disposable disinfectant wipe.

**Phase II**

- Wipe all surfaces of the optically clean product with a fresh, disposable disinfectant wipe.
- Observe the application time (1 min minimum).

### 3.6 Inspection

- Allow the product to cool down to room temperature.
- Dry the product if it is wet or damp.

#### 3.6.1 Visual inspection

- Make certain that all soiling has been removed. In particular, pay attention to mating surfaces, hinges, shafts, recessed areas, drill grooves and the sides of the teeth on rasps.
- If the product is dirty: repeat the cleaning and disinfection process.
- Check the product for damage, e.g. insulation or corroded, loose, bent, broken, cracked, worn or severely scratched and fractured components.
- Check the product for missing or faded labels.
- Check the products with long, slim shapes (in particular rotating instruments) for deformities.
- Check the surfaces for rough spots.
- Check the product for burrs that could damage tissue or surgical gloves.
- Check the product for loose or missing parts.
- Immediately put aside damaged or inoperative products and send them to Aesculap Technical Service, see Technical service.

#### 3.6.2 Functional test

- Check that the product functions correctly.
- Check that all moving parts are working properly (e.g. hinges, locks/latches, sliding parts etc.).
- Check the product for any atypical running noise, overheating or excessive vibration.
- Check for compatibility with associated products.
- Immediately put aside inoperative products and send them to Aesculap Technical Service, see Technical service.

### 3.7 Storage

- Store sterile products in germ-proof packaging, protected from dust, in a dry, dark, temperature-controlled area.

#### 3.7.1 Ambient conditions

The following environmental conditions apply to the transport and storage of the product:

Temperature	-10 °C to 50 °C
Relative humidity	10 % to 90 %
Atmospheric pressure	500 hPa to 1 060 hPa

## 4. Maintenance and service

### 4.1 Maintenance

This medical product contains no components or parts that must be replaced in the form of service at regular intervals specified by the manufacturer.

A repeat inspection of the medical product must be performed every 12 months. This must also be carried out after any repair and after it was dropped, damaged or misused.

The repeat inspection may only be carried out by persons authorized to do so by the manufacturer, using the service manual.

- Comply with applicable national and international standards. For services to that end, please contact your national B. Braun/Aesculap agency, see Technical service

### 4.2 Technical service

#### DANGER

Danger to life of patients and users if the product malfunctions and/or protective measures fail or are not used!

- Do not perform any servicing or maintenance work under any circumstances while the product is being used on a patient.

Modifications carried out on medical technical equipment may result in loss of guarantee/warranty rights and forfeiture of applicable licenses.

- Do not modify the product.
- For service and repairs, please contact your national B. Braun/Aesculap agency.

#### Service addresses

Aesculap Technischer Service

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78532 Tuttlingen / Germany

Phone: +49 7461 95-1601

Fax: +49 7461 16-2887

E-Mail: ats@aeculap.de

Other service addresses can be obtained from the address indicated above.

### 4.3 Accessories / spare parts

All accessories and spare parts must only be procured from the manufacturer.

#### 4.3.1 Accessories for sterile adapter

Art. no.	Designation
OP941	Sterile adapter

#### 4.3.2 Accessories for power cord

Art. no.	Approvals	Color	Length
TE676	IEC power cable	gray	1.0 m
TE730	Europe	black	5.0 m
TE734	Great Britain	black	5.0 m
TE735	USA, Canada, Japan	gray	3.5 m
TE736	IEC power cable	gray	2.5 m
TE780	Europe	black	1.5 m

### 4.3.3 Spare part for fuse

Art. no.	Designation
TA022249	Fuse set

### 4.3.4 Communications accessories

Art. no.	Designation
OP942	Communications interface cable, 0.75 m

### 4.3.5 Light lead accessories

Art. no.	Designation
EV2-000022	Light lead CF, Ø 4.8 mm, 5.0 m
OP906	Light lead, Ø 4.8 mm, 1.8 m
OP913	Light lead, Ø 4.8 mm, 2.5 m
OP914	Light cable, Ø 4.8 mm, 3.5 m
OP923	FullHD light lead, Ø 4.8 mm, 2.5 m

### 4.3.6 Equipotential bonding cable accessories

Art. no.	Designation
TA008205	Equipotential bonding cable, 0.8 m
GK535	Equipotential bonding cable, 4.0 m

## 5. Disposal

### ⚠ WARNING

Risk of infection due to contaminated products!

- Adhere to national regulations when disposing of or recycling the product, its components and its packaging.

### Note

The user institution is obliged to reprocess the product before its disposal, see Validated reprocessing procedure.

The recycling pass can be downloaded from the Extranet as a PDF document under the respective article number. (The recycling pass includes disassembling instructions for the product, as well as information for proper disposal of components harmful to the environment.)

Products carrying this symbol are subject to separate collection of electrical and electronic devices. Within the European Union, disposal is taken care of by the manufacturer as a free-of-charge service.

- Detailed information concerning the disposal of the product is available through your national B. Braun/Aesculap agency, see Technical service.

## 6. Technical data

### 6.1 Classification acc. to Regulation (EU) 2017/745

Art. no.	Designation	Class
OP950	LED light source	I
OP941	Sterile adapter	I

### 6.2 Performance data, information about standards

Protection class (acc. to IEC/EN/DIN 60601-1)	I
Housing protection rating according to IEC/EN/DIN 60529	IP21
Applied part	Type CF, defibrillator-safe
Mains voltage ranges	100–240 V~
Current consumption (ready for operation)	0.07 A (at 100 V~) 0.09 A (at 240 V~)
Current consumption (maximum load)	2.0–1.0 A
Frequency	50/60 Hz
Operating mode	Continuous operation
Device fuse (acc. to IEC 60127-1)	T 2.0 AH, 250 V Structural design: 5 x 20 mm
Weight	8.5 kg
Dimensions (W x H x D)	330 mm x 146 mm x 362 mm
Conforming to standard	IEC/EN/DIN 60601-1
EMC	IEC/EN/DIN 60601-1-2
CISPR11	Class B

## 7. Symbols on product and packaging

	General warning sign	
	Meaning dependent on the standard that applies in each case: 2nd edition IEC 60601-1: Please take note of the accompanying document! 3rd edition IEC 60601-1: Caution	
	Follow the instructions for use	
	"OFF" (power supply)	
	"ON" (power supply)	
	Defibrillation-protected applied part of type CF	
	Optical cable connection	
	Equipotentialization connector	
	Unlock	
	Hazardous electrical voltage	
	Fuse	
	Alternating current	
	Manufacturer	
	Labeling of electrical and electronic devices according to directive 2012/19/EU (WEEE), see Disposal	
	Temperature limits during transport and storage	
	Air humidity limits during transport and storage	
	Atmospheric pressure limits during transport and storage	
	0123 0482	CE marking according to regulation (EU) 2017/745
		Possibly hazardous optical radiation emitted from this product! Do not stare at operating lamp. May be harmful for the eye. TPE031-121-01
		Medical Device

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Aesculap AG | Am Aesculap-Platz | 78532 Tuttlingen | Germany